EN
3
Colon Decompression Sets Marcon / Brodmerkel

Instructions for Use

DA
4
Marcon / Brodmerkel kolondekompressionssæt

Brugsanvisning

NL
5
Kolondecompressieset Marcon / Brodmerkel

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
Bruksanvisning

INTENDED USE

This device is used for treatment of acute non-toxic megacolon, pseudo-obstruction (Ogilvie's syndrome) and colonic strictures.

NOTES

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

Set includes a drainage catheter and a .035 inch wire guide.

Do not use this device for any purpose other than stated intended use.

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Store in a dry location, away from temperature extremes.

Use of this device restricted to a trained healthcare professional.

CONTRAINDICATIONS

Those for colonoscopy include, but are not limited to: fulminant colitis, possible existing perforation, acute severe diverticulitis, presence of barium or inadequate bowel preparation.

Those for colon decompression include, but are not limited to: ischemic bowel, recent rectal surgery, organic colon obstruction, recent surgical anastomosis, acute myocardial infarction.

POTENTIAL COMPLICATIONS

Those associated with colonoscopy include, but are not limited to: perforation, hemorrhage, fever, infection, allergic reaction to medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

Those associated with colon decompression include, but are not limited to: clogging of tube, dislodgment of tube as colonoscope is withdrawn, perforation.

PRECAUTIONS

Careful history and physical examination must be performed prior to tube insertion in order to determine underlying cause of pseudo-obstruction.

Colon should not be prepped unless there is inordinate amount of stool in rectum.

Refer to package label for minimum channel size required for this device.

INSTRUCTIONS FOR USE

Illustrations

1. Advance colonoscope, insufflating minimal air, until it has reached cecum or as far in colon as possible. **Note:** Suction must be performed intermittently throughout procedure to help decompress colon and avoid further colonic distention.

2. Fluoroscopically observe location of colonoscope tip.
3. Preload the decompression tube with wire guide until wire exits tip of tube.
4. Advance preloaded tube through colonoscope to desired location.
5. Remove drainage connector and colonoscope. (See fig. 1)
6. Remove wire guide and attach drainage connector.
7. Secure external portion of decompression tube to patient.

Upon completion of procedure, dispose of device(s) per institutional guidelines for biohazardous medical waste.

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© COPYRIGHT COOK 2013 ©
COOK IRELAND LTD.