COOK®

MEDICAL

CE
0088

EN
2

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Bruksanvisning

F M - 2 7 9 5 - R E V 1

Evolution® Shortie Mechanical Dilator Sheath Set

Federal law restricts this device to sale by or on the order of a physician.

DEVICE DESCRIPTION

The Evolution® Shortie Mechanical Dilator Sheath Set is comprised of an inner polymer sheath connected to a handle capable of mechanically rotating the sheath and an outer telescoping polymer sheath. The inner sheath has a stainless steel tip at its distal end.

INTENDED USE

Evolution® Shortie Mechanical Dilator Sheath Set is intended for use in patients requiring the percutaneous dilation of tissue surrounding cardiac leads, indwelling catheters and foreign objects.

Use With Other Devices

The Evolution® Shortie Mechanical Dilator Sheath Set may be used in conjunction with the following catheter/lead extraction devices from Cook:

Liberator® Locking Stylet

Byrd Telescoping Stainless Steel Dilator Sheath Set

Byrd Workstation™

Needle's Eye Snare®

Be sure to closely follow the "Suggested Instructions for Use" for each device used.

CONTRAINDICATIONS

None known

WARNINGS

Evolution® Shortie Mechanical Dilator Sheath Set should only be used to minimally enter the vessel. Do not attempt to negotiate the Evolution® Shortie Mechanical Dilator Sheath Set past a bend in the vessel as vessel wall or cardiac lead insulation damage may occur. When using dilator sheaths or sheath sets, do not insert more than one sheath set into a vein at a time. Severe vessel damage, including venous wall laceration requiring surgical repair, may occur.

When using a Locking Stylet:

Do not abandon a catheter/lead in a patient with a Locking Stylet still in place inside the catheter/lead. Severe vessel or endocardial wall damage may result from the stiffened catheter/lead or from fracture or migration of the abandoned stylet wire.

Do not apply weighted traction to an inserted Locking Stylet as myocardial avulsion, hypotension, or venous wall tearing may result.

Be aware that a lead that has a J-shape retention wire that occupies its inner lumen (rather than being outside the coil) may not be compatible with the Locking Stylet. Insertion of the Locking Stylet into such a lead may result in protrusion and possible migration of the J-shape retention wire.

Weigh the relative risks and benefits of intravascular catheter/lead removal procedures in cases when:

The item to be removed is of a dangerous shape or configuration,
The likelihood of catheter/lead disintegration resulting in fragment embolism is high or vegetations are attached directly to the catheter/lead body.

Catheter/lead removal devices should be used only at institutions with thoracic surgical capabilities.

Catheter/lead removal devices should be used only by physicians knowledgeable in the techniques and devices for catheter/lead removal.

PRECAUTIONS

When activating the Evolution® Shortie Mechanical Dilator Sheath Set be careful not to cut or ensnare lead suture-ties that are attached to the locking stylet targeted lead or to companion leads.

Position the distal tip of the Evolution® Shortie Mechanical Dilator Sheath Set so that it does not engage companion leads during activation.

Prior to the procedure, consider the size of the catheter/lead in relation to the size of the LEAD EXTRACTION™ devices to determine possible incompatibility.

If selectively removing catheters/leads with the intent to leave one or more chronic catheters/leads implanted intact, the non-targeted catheters/leads must be subsequently tested to ensure that they were not damaged or dislodged during the extraction procedure.

As a result of the risk of complications, doctors highly experienced in this procedure have advised the following precautions:

PREPARATIONS

Obtain a thorough patient history, including patient blood type. Appropriate blood products should be rapidly available.

Ascertain the manufacturer, model number, and implant date of the catheter/lead to be removed.

Perform radiographic/echocardiographic evaluation of catheter/lead condition, type, and position. Use a procedure room that has

high quality fluoroscopy, pacing equipment, defibrillator, thoracotomy tray, and pericardiocentesis tray.

Echocardiography and cardiothoracic surgery should be rapidly available.

Prep/drape the patient's chest for possible thoracotomy; prep/drape the patient's groin for femoral approach or possible femoral approach.

Establish back up pacing as necessary.

Have available an extensive collection of sheaths, Locking Stylets, stylets to unscrew active fixation leads, snares, and accessory equipment.

PROCEDURE

Use fluoroscopic monitoring during ALL catheter/lead and sheath manipulations. Monitor ECG and arterial blood pressure continuously throughout the procedure and during recovery.

If using sheaths or sheath sets, including the Evolution® Shortie Mechanical Dilator Sheath Set, the following precautions should be followed:

Prior to using sheaths including the Evolution® Shortie Mechanical Dilator Sheath Set it is essential to carefully inspect the extravascular catheter/lead tract to ensure removal of all suture sleeves, sutures, and tie-down materials.

When advancing sheaths including the Evolution® Mechanical Dilator Sheath Set, use proper sheath technique and maintain adequate tension on the catheter/lead (via a Locking Stylet or directly) to prevent damage to vessel walls.

If excessive scar tissue or calcification prevents safe advancement of sheaths, consider an alternate approach.

Excessive force with sheaths used intravascularly may result in damage to the vascular system requiring surgical repair.

If the catheter/lead breaks, evaluate fragment; retrieve as indicated.

If hypotension develops, rapidly evaluate; treat as appropriate.

Due to rapidly evolving catheter/lead technology, this device may not be suitable for the removal of all types of catheters/leads. If there are questions or concerns regarding compatibility of this device with particular catheters/leads, contact the catheter/lead manufacturer.

POTENTIAL ADVERSE EVENTS

Potential adverse events related to the procedure of intravascular extraction of catheters/leads include (listed in order of increasing potential effect):

- dislodging or damaging nontargeted catheter/lead
- chest wall hematoma
- thrombosis
- arrhythmias
- acute bacteremia
- acute hypotension
- pneumothorax
- stroke
- migrating fragment from catheter/object
- pulmonary embolism
- laceration or tearing of vascular structures or the myocardium
- hemopericardium/pericardial effusion
- cardiac tamponade
- hemothorax
- cardiac arrest
- death

SUGGESTED INSTRUCTIONS FOR USE

Evolution® Shortie Mechanical Dilator Sheath Set

NOTE: In addition to closely following the "Suggested Instructions for Use" for this device, refer to the "Clinical Considerations for Removing Leads or Catheters: Superior Approach" shown in the next section of this insert. When using dilator sheaths or sheath sets, including the Evolution® and Evolution® Shortie Mechanical Dilator Sheath Set, the following precautions should be followed:

Prior to using dilator sheaths, it is essential to carefully inspect the extravascular catheter/lead tract to ensure removal of all suture sleeves, sutures, and tie-down materials.

When advancing the Evolution® Shortie Mechanical Dilator Sheath Set, use proper sheath technique and maintain adequate tension on the catheter/lead (via a Locking Stylet or directly) to prevent damage to vessel walls.

If excessive scar tissue or calcification prevents safe advancement of dilator sheaths, consider an alternate approach.

Excessive force with dilator sheaths, including the Evolution® Shortie Mechanical Dilator Sheath Set used intravascularly may result in damage to the vascular system requiring surgical repair.

WARNING: When using dilator sheaths or sheath sets, do not insert more than one sheath set into a vein at a time. Severe vessel damage, including venous wall laceration requiring surgical repair, may occur.

1. With the inner Evolution® Shortie Mechanical Dilator Sheath placed inside the appropriate outer sheath for telescopic action, advance the sheaths over the lead as described in Part III of **“Clinical Considerations for Removing Leads or Catheters: Superior Approach”** shown in the next section of this insert. To rotate the inner Evolution® Shortie Mechanical Dilator Sheath, slowly squeeze and release the trigger activation handle (one squeeze every one to two seconds max). Repeat as appropriate. Immediately end the use of the device when the first of either, minimal entry into the vessel has been achieved or fluoroscopic monitoring confirms that a mid-clavicular location along the cardiac lead has been reached. Exchange the Evolution® Shortie Mechanical Dilator Sheath Set for the appropriately sized Evolution® Mechanical Dilator Sheath Set.

WARNING: Do not attempt to negotiate the Evolution® Shortie Mechanical Dilator Sheath Set past a bend in the vessel as vessel wall or cardiac lead insulation damage may occur.

2. To free the catheter/lead from tissue encapsulation refer to Part III of **“Clinical Considerations for Removing Leads or Catheters: Superior Approach”** shown in the next section of this insert.

CLINICAL CONSIDERATIONS FOR REMOVING LEADS OR CATHETERS: SUPERIOR APPROACH

Clinical experience in the removal of leads has identified several considerations for lead removal techniques using superior or femoral approaches. Physicians highly experienced with lead removal techniques have suggested the following considerations for removal of leads using a superior approach.

Part I. General Considerations for Accessing the Catheter/Lead

1. Surgically expose the proximal end of the indwelling catheter/lead and remove the catheter/lead from its connections (if connected).
2. Remove all suture and tie-down materials.

NOTE: Because implant procedures vary dramatically, catheters/leads may be secured in a variety of ways at the time of implant. Carefully inspect the site to ensure all suture sleeves, sutures, and tie down materials, which may be hidden or blind have been removed.

3. Cut off all proximal fittings, if present, using clippers or other cutters. It is important to cut the catheter/lead very close to the connector (but past any crimp joints) leaving as long a portion of the indwelling catheter/lead to work with as possible. Avoid closing off the interior lumen (or coil) of the catheter/lead when cutting it.

NOTE: Do not pull on the catheter/lead because it may stretch, distort, or break, making subsequent removal more difficult. Damage to the lead may prevent passage of a Locking Stylet through the lumen and/or make dilation of scar tissue more difficult.

Part II. General Considerations for Using a Locking Stylet

4. Consider passing a LIBERATOR® Locking Stylet (Cook) through the inner lumen of the catheter/lead to stabilize the catheter/lead during dilation of surrounding tissues. Closely follow the “Suggested Instructions for Use” for the Locking Stylet to:
 - A. Expose the inner coil of the catheter/lead
 - B. Check patency of the coil lumen
 - C. Determine the appropriate size of Locking Stylet based on the inner diameter of the catheter/lead
 - D. Advance the Locking Stylet to the distal end of the catheter/lead
 - E. Lock the Locking Stylet in place

WARNING: Be aware that a lead that has a J-shape retention wire that occupies its inner lumen (rather than being outside the coil) may not be compatible with the Locking Stylet. Insertion of the Locking Stylet into such a lead may result in protrusion and possible migration of the J-shape retention wire.

5. Unless the catheter/lead insulation is damaged, degraded or too thin, tie a ligature or use a compression coil (Cook) at the proximal end of the catheter/lead, compressing the insulation against the coil and Locking Stylet to help prevent the coil and insulation from stretching. The ligature can be tied to the loop handle or to the suture tie loop.
6. For an active fixation catheter/lead, attempt to unscrew the catheter/lead by rotating the catheter/lead (and Locking Stylet if used) counterclockwise.
7. Gently pull back on the catheter/lead to see if it is still engaged in tissue. If the catheter/lead is sufficiently loose in the tissue, gently pull on the Locking Stylet (if used) and catheter/lead to remove it.

NOTE: If a Locking Stylet has not been used, be aware that damage to the catheter/lead caused by pulling on it may prevent subsequent passage of a Locking Stylet through the lumen and/or make dilation of scar tissue more difficult.

NOTE: If removing a chronic pacing lead, be aware that if it is freed spontaneously during the extraction procedure, the lead tip may become trapped in the upper vasculature. Dilator sheaths, advanced

at least to the innominate vein, are often necessary to extract the lead tip through the scar tissue at the site of venous entry, and to avoid a venotomy.

8. If the catheter/lead is not removed from the vessel with gentle pulling, dilator sheaths, including the Evolution® Shortie Mechanical Dilator Sheath Set, or other retrieval devices may help separate the catheter/lead from any tissue encapsulation as follows:

Collapse the proximal loop of the Locking Stylet (if used).

Advance dilator sheaths (or other retrieval device) over collapsed proximal loop or the proximal wire extender of the Locking Stylet (if applicable).

NOTE: If using dilator sheaths or dilator sheath sets, including the Evolution® and Evolution® Shortie Mechanical Dilator Sheath Set (Cook), consider the "General Considerations" shown below and, in addition, closely follow the "Suggested Instructions for Use" for the particular type of sheath or sheath set.

Part III. Considerations for Use of Dilator Sheaths

General Considerations

Evolution® Shortie Mechanical Dilator Sheath Sets and Byrd Telescoping Stainless Steel Dilator Sheath Sets are used to minimally enter the vessel. Upon vessel entry, these sheath sets may be exchanged for the Evolution Mechanical Dilator Sheath Set or Byrd Dilator Sheath Set Polypropylene which then may be advanced over the catheter/lead in a telescoping manner, while maintaining constant tension.

If scar tissue makes vessel entry difficult when using a Byrd Dilator Sheath Set Polypropylene or Evolution® Mechanical Dilator Sheath Set, consider using the Evolution® Shortie Mechanical Dilator Sheath Set or Byrd Telescoping Stainless Steel Dilator Sheath Set to obtain vessel access (refer to the "Suggested Instructions for Use" for the Byrd Telescoping Stainless Steel Dilator Sheath Set).

The size of the dilator sheaths used should be large enough that the sheaths can be advanced over the catheter/lead without causing the catheter/lead to buckle or its insulation to bunch up, but the sheaths should not be overly loose.

Always use fluoroscopic monitoring when advancing the dilator sheaths. Never use sheaths on more than one catheter or lead at a time. **Always maintain adequate tension on the catheter/lead to support the maneuvering of the dilator sheaths and properly guide them within the vessel.** With too little tension, the sheaths may rupture the vein. Too much tension may cause a myocardial avulsion.

NOTE: Follow closely the "Suggested Instructions for Use" for the particular type of sheath or sheath set being used.

9. With the inner Evolution® or Evolution® Shortie Mechanical Dilator Sheath placed inside the appropriate outer sheath for telescopic action, insert the proximal free end of the catheter/lead into the distal end of the inner Evolution® or Evolution® Shortie Mechanical Dilator Sheath. Advance the catheter/lead until it completely exits the opposite (proximal) end of the sheath set.

NOTE: Consider using a hooked stylet or small basket to aid in drawing the catheter/lead through the sheath set.

10. Apply adequate retracting pressure or tension on the catheter/lead and/or its wire guide or Locking Stylet (if used). **This is critical to safe passage of the sheath set over the catheter/lead. If tension is inadequate, the catheter/lead may buckle, precluding the sheath set from advancement along the appropriate path.**

NOTE: If the catheter/lead is entangled, other retrieval devices such as grasping devices, J guide wires, deflecting catheters or baskets may be helpful to grasp and straighten the catheter/lead before using dilator sheaths.

NOTE: If using Evolution® Shortie Mechanical Dilator Sheath Set or Stainless Steel Dilator Sheath Set, when minimal entry into the vessel has been verified, **exchange the Evolution® Shortie Mechanical Dilator Sheath Set or stainless steel dilator sheath(s) for the Evolution® Mechanical Dilator Sheath Set** as described in the "Suggested Instructions for Use".

NOTE: The Evolution® Shortie Mechanical Dilator Sheath Set should be used only to minimally enter the vessel.

11. With the catheter/lead in tension, and under fluoroscopic guidance, advance the inner sheath along the length of the catheter/lead and, if applicable, into the blood vessel. Rotate the inner Evolution® Mechanical Dilator Sheath, by squeezing and releasing the trigger activation handle. Repeat as appropriate. Advance the outer sheath over the inner sheath, keeping it within the vessel.
12. In a telescopic manner, alternately advance the inner, and then outer, sheaths while maintaining adequate tension on the catheter/lead and/or its wire guide or Locking Stylet, until the catheter/lead is free of tissue restriction. In general, the inner sheath should not be advanced more than 4 cm beyond the outer sheath.

NOTE: Always monitor fluoroscopically while advancing the sheaths to avoid shearing the catheter/lead or rupturing the vessel. Continue to maintain adequate tension on the catheter/lead (or Locking Stylet). Adjust the angle of entry to keep the sheaths aligned with the catheter/lead in the vessel and to minimize sheath curvature during advancement. Rotating the inner Evolution® Mechanical Dilator Sheath using the trigger activation handle and its outer sheath during advancement may facilitate progress through exuberant scar tissue.

NOTE: If the insulation separates and exposes the inner coil, watch the tips of the outer sheath to be sure that it does not split or tear, because this could lead to vessel damage. Advance a catheter or small sheath the same size as the insulation over the bare wire, if necessary.

NOTE: If the sheaths fail to progress after initial success, or if advancing the sheaths was difficult, remove the sheaths one at a time to inspect the tips; if the tips are distorted, exchange them for a new set.

NOTE: When tracking the lead around a bend, keep the point of the tip (if angled) on the inside of the bend if applicable.

NOTE: Be careful when trying to advance the inner sheath past the ring electrode of a bipolar lead: if it gets caught, the lead may break.

13. For cardiac leads, if the lead has not been freed by the time the sheaths near the myocardium, position the outer sheath so that the blunt end is directed toward the myocardium. Pull the inner sheath back several centimeters.

Apply countertraction with the outer sheath as follows:

Firmly hold the outer sheath one centimeter from the heart wall and steadily pull the Locking Stylet back; the pacing tip will be pulled into the sheath. Rotation of the sheath may help dislodge the tip.

14. When the catheter/lead has become unentangled and freed from tissue attachment, remove the catheter/lead, Locking Stylet if used, and Evolution® Mechanical Dilator Sheath Set (or other retrieval device) if used.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

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Cook Ireland LTD
Cook Vandergrift Inc.
Vandergrift, PA 15690 U.S.A.