EN Extraction Balloon

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FR Ballonnet d'extraction

DE Extraktionsballon

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ENGLISH

INTENDED USE

This device is used for endoscopic removal of stones in the biliary system and for contrast injection.

NOTES

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

If the package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Notify Cook for return authorization.

Do not use this device for any purpose other than the stated intended use.

Store in a dry location, away from temperature extremes.

Use of this device restricted to a trained healthcare professional.

CONTRAINDICATIONS

Contraindications include those specific to ERCP and any procedures to be performed in conjunction with balloon stone extraction.

Use of this natural latex rubber balloon is contraindicated in patients with a known hypersensitivity to latex.

POTENTIAL COMPLICATIONS

Potential complications associated with ERCP include, but are not limited to: pancreatitis, cholangitis, aspiration, perforation, hemorrhage, infection, sepsis, allergic reaction to contrast or medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

Additional complications that can occur during endoscopic balloon extraction include, but are not limited to: stone impaction, localized inflammation, pressure necrosis.

PRECAUTIONS

Refer to package label for minimum channel size required for this device.

The coordination of endoscope channel size with compatible devices is essential in obtaining optimal results during a procedure. Refer to product package label for the minimum channel size required for this device.

The wire guide diameter and the inner lumen of the wire-guided device must be compatible.

Assessment of the size of the stone and ampullary orifice must be made to determine the necessity of sphincterotomy. In the event sphincterotomy is required, all necessary precautions must be followed.

Injection of contrast during ERCP must be monitored fluoroscopically. Overfilling of the pancreatic duct may cause pancreatitis.

SYSTEM PREPARATION

Verify balloon integrity prior to use by attaching the enclosed syringe to stopcock and inflating balloon with air only. If any leakage is detected, **do not use**. Please notify Cook for return authorization.

Flush injection port with contrast to expel all air.

Note: The stopcock is in the open position, allowing access to the balloon, when the two arms are aligned with the catheter and syringe. To maintain balloon inflation, turn the stopcock arm 90°.

INSTRUCTIONS FOR USE

I. If using Fusion Quattro Balloon and a prepositioned short wire. (See fig. 1)

1. Unlock short wire from Wire Guide Locking Device and advance tip of balloon catheter onto prepositioned wire guide, ensuring wire guide exits catheter at Zip port.
2. Introduce deflated balloon into endoscope accessory channel and relock wire guide. Then advance device in short increments until it is endoscopically visualized exiting endoscope.
3. Using fluoroscopic monitoring, position deflated balloon above stone to be removed.
Note: If more than one stone is to be removed, extract one stone at a time.
4. After verifying desired position of balloon, inflate balloon with air only. **Note for multiple sizing balloons:** Inflate balloon using fluoroscopic monitoring until balloon is visualized occluding duct. If desired, adjust size of balloon by using reference marks on syringe. To achieve smallest balloon size, inflate balloon to next largest size and gently pull back on syringe to initial size. Lock stopcock.
5. Using fluoroscopic visualization and keeping endoscope elevator open, gently withdraw the inflated balloon toward the papilla. **Warning:** Do not exert excessive pressure on ampulla while extracting stones. If stone does not pass easily, reassess need for sphincterotomy.
6. Once balloon is endoscopically visualized in duodenum, turn stopcock to open position and deflate balloon.
7. Repeat extraction process, one stone at a time, until duct is clear.

Note: The previously placed wire guide may be left in position, in order to facilitate introduction of other wire-guided devices.

If wire guide is to remain in place while device is withdrawn, utilize following steps:

8. Remove device from endoscope accessory channel until resistance is achieved. Unlock wire guide from Wire Guide Locking Device and perform a short segment wire guide exchange. Relock wire guide into Wire Guide Locking Device and remove device from endoscope accessory channel.

II. If using Fusion Extraction Balloon and a prepositioned short wire. (See fig. 2)

1. Unlock short wire from Wire Guide Locking Device and advance tip of balloon catheter onto prepositioned wire guide ensuring wire guide exits catheter at IDE port.

REFER TO STEPS 2-7 IN "SECTION I", THEN RESUME WITH STEP 2 BELOW:

Note: Previously placed wire guide may be left in position, in order to facilitate introduction of other wire-guided devices.

2. Prior to removing device, utilize reference marks on catheter to ensure IDE port is within ductal system.
3. Fluoroscopically visualize radiopaque band at IDE port. Retract wire guide until radiopaque distal tip of wire guide passes band; a disengagement from wire guide lumen will occur.
4. Advance disengaged wire guide to maintain ductal access.
5. Lock wire guide into Wire Guide Locking Device and remove device from endoscope accessory channel.

III. If using Proximal Wire Port (PWP) and a prepositioned long wire guide. (See fig. 1)

Note: For best results, wire guide should be kept wet.

1. Remove stylet.
2. Advance tip of balloon catheter over prepositioned long wire guide, ensuring wire guide exits catheter at PWP.
3. Advance deflated balloon in short increments through accessory channel using standard long wire exchange technique until it is visualized exiting endoscope.

REFER TO STEPS 3-7 IN "SECTION I", THEN RESUME WITH STEP 4 BELOW:

4. Remove balloon using standard long wire exchange technique.

Upon completion of procedure, dispose of device per institutional guidelines for biohazardous medical waste.

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Wilson-Cook Medical, Inc.
Cook Ireland Ltd.