EN
3

Hemospray® Endoscopic Hemostat

FR
6

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






ENGLISH

INTENDED USE

This device is used for hemostasis of nonvariceal gastrointestinal bleeding.

NOTES

Do not use this device for any purpose other than stated intended use.

Store in a dry location, away from temperature extremes.

Use of this device is restricted to a trained healthcare professional.

DEVICE DESCRIPTION

Hemospray is an inert, bentonite powder developed for endoscopic hemostasis. The powder is delivered by use of a carbon dioxide powered delivery system and through a catheter inserted through the working channel of an endoscope which provides access to the site of the bleed. Each device contains approximately 20g of powder.

CONTRAINDICATIONS

Those specific to primary endoscopic procedure to be performed in gaining access to desired target site. Also contraindicated in patients who have gastrointestinal fistulas, are suspected of having a gastrointestinal perforation, or are at high risk of gastrointestinal perforation during endoscopic treatment.

WARNING

Keep catheter tip at least 1 cm away from bleeding site to minimize risk of embolization.

Hemospray has not been approved for use in pediatric populations, no safety or effectiveness data exists and would be considered off-label usage, to do so is at the professional risk of the surgeon/healthcare professional.

Patients with gastrointestinal bleeding that are on antithrombotic medication may be at an increased risk of rebleeding. Follow the

relevant clinical guidelines for management of antithrombotic agents for endoscopic procedures.

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working conditions, do not use. Please notify Cook for return authorization.

PRECAUTIONS

For best results, Hemospray must be delivered to the source of bleeding.

Ensure gastrointestinal lumen is not distended because Hemospray adds volume in excess of insufflation volumes during procedure. Closely monitor bowel distension and balance insufflation and Hemospray volumes as necessary.

Use of more than (3) Hemospray devices per patient may result in impaction in colon and is not recommended.

Product contains 16 g CO₂ cartridge. Contents under pressure. Do not puncture or heat above 120° F / 49° C. Do not inhale or discharge towards face/body. Keep CO₂ dispenser and cartridge out of reach of children.

Hemospray is inert and non-toxic.

As a granular material, unintentional exposure to the powder may cause potential irritation to the skin, eyes and lungs.

In the event of unintended exposure to the powder refer to the following First Aid measures:

Skin: Wash with soap and water until clean.

Eyes: Flush with water until irritation ceases.

Inhalation: Move to area free from powder. If symptoms of irritation persist, contact physician. Inhalation may aggravate existing respiratory illness.

Refer to package label for minimum channel size required for this device.

All endoscopic hemostatic therapies, including Hemospray, have an associated risk of rebleeding, particularly in situations where the cause of bleeding is an unresolved underlying disease. After hemostasis has been achieved, monitor patients for rebleeding per the relevant Clinical Guidelines.

Although not seen in clinical practice, there remains a theoretical risk of aspiration of Hemospray powder resulting in respiratory complications. It is prudent to restrict the use of Hemospray to 5 cm below the UES.

PROCEDURAL PRECAUTIONS

Like other modalities, Hemospray may not be effective for all types of bleeds. Gastrointestinal bleeding may exacerbate existing comorbidities, increasing the potential for adverse events including patient mortality.

POTENTIAL COMPLICATIONS

Those associated with gastrointestinal endoscopy include, but are not limited to: perforation, hemorrhage, aspiration, fever, infection, allergic reaction to medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

Use of Hemospray in the presence of bowel obstruction and/or an anastomosis may pose a risk of injury due to over-distention.

Hemospray may occlude ducts and orifices which communicate with the main bowel lumen. Use caution when using Hemospray in the vicinity of these orifices.

Others include, but are not limited to: powder impaction in colon or embolization.

When spraying in retroflexed position, Hemospray powder may adhere to the outside of the endoscope. This may result in difficulty repositioning/ removing the endoscope, particularly if passing through a strictured area.

MRI SAFETY INFORMATION

While Adverse Effects related to MR safety have not been observed in clinical use, Hemospray has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of Hemospray in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

SYSTEM PREPARATION

1. Remove device from package and attach catheter to handle, ensuring connection is secure (*See Fig 1*).
2. Activate CO₂ cartridge by turning red activation knob until it stops.
Note: Do not over rotate knob as this could damage the device (*See Fig. 2*). **Note:** Do not test device prior to insertion into endoscope accessory channel as this may increase risk of catheter occlusion.

INSTRUCTIONS FOR USE

1. Before inserting catheter into accessory channel, identify bleeding site, remove as much blood as possible, then flush accessory channel with air. **Caution:** Ensure gastrointestinal lumen is not distended because Hemospray adds volume during procedure.
2. Slowly advance catheter through accessory channel in short increments until catheter tip is visualized endoscopically (*See Fig. 3*).
Precaution: To avoid catheter occlusion, **do not** place catheter directly in contact with blood and/or mucosa, including any pooled blood and **do not** aspirate blood while catheter is in accessory channel.
3. To ensure proper visibility, catheter tip should be 1-2 cm away from bleeding site at all times (*See Fig. 4*).
4. To allow powder deployment, turn red valve to open position (*See Fig. 5*). **Note:** Device is now active and ready for use. **Do not** press trigger button until powder deployment is desired.
5. To deploy powder, hold handle upright and depress red trigger button for 1-2 seconds and release (*See Fig. 6*). Continue applying

in short bursts until site is completely covered with powder and no active bleeding is visualized. **Note:** Allow time for hemostasis before aspirating. Hemostasis may not be immediate because the powder takes time to absorb. Multiple applications may be required; however no more than (3) Hemospray devices should be applied per patient.

Note: Irrigation of treated site is not required after deployment. If irrigation is desired, allow time for hemostasis before irrigating with water. **Note:** To avoid potential endoscope occlusion, do not aspirate powder into endoscope channel.

6. Prior to removing Hemospray device from patient, turn red valve to closed position (*See Fig. 7*). To remove device from endoscope, grasp catheter and pull slowly. **Note:** If catheter becomes occluded, turn red valve to closed position, remove catheter from endoscope and replace with extra catheter provided in package. Then repeat steps 1-6 of "*Instructions for Use*." If more powder is needed and spray canister and/or CO₂ cartridge is empty, a new device is required.

Upon completion of procedure, depressurize device by rotating activation knob counterclockwise until CO₂ cartridge depressurizes completely. Dispose of device per institutional guidelines for biohazardous medical waste.

SUMMARY OF CLINICAL DATA

Clinical data summary information that was, in part, the basis for granting the de novo can be found on the Cook Medical website at CookMedical.com/HemosprayData

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A symbol glossary can be found at cookmedical.com/symbol-glossary

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If symbol appears on product label, X = maximum insertion portion width

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If symbol appears on product label, X = minimum accessory channel

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If symbol appears on product label, X = wire guide compatibility

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If symbol appears on product label, X = quantity per box

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This device is intended
for single use only.

Wilson-Cook Medical, Inc.
4900 Bethania Station Road
Winston-Salem, North Carolina 27105
USA

RxOnly STERILE R