## Retrograde Intubation Sets

Instructions for Use

EN  
5

## Retrograd intubationssæt

Brugsanvisning

DA  
6

## Retrograde Intubations-Sets

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DE  
7

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
**Fig. 5**

## RETROGRADE INTUBATION SETS

**CAUTION:** U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

### DEVICE DESCRIPTION

Retrograde Intubation Sets consist of the following components:

- Introducer needle
- Catheter introducer needle
- .038 inch wire guide with two positioning marks, 20 cm from each end
- Radiopaque catheter
- Syringe

### INTENDED USE

The Retrograde Intubation Set is intended to assist in the placement of an endotracheal tube during difficult or emergency airway access procedures where visualization of the vocal cords is not possible secondary to secretions, blood and/or anatomic anomalies.

Refer to the table below for sizing information.

Retrograde Intubation Set Size	For Use With Endotracheal Tube (ID)
6.0 Fr	2.5 mm or larger
11.0 Fr	4.0 mm or larger
14.0 Fr	5.0 mm or larger

### CONTRAINDICATIONS

- Ongoing coagulopathy
- Obscure cricothyroid anatomy
- Infection of cricothyroid membrane
- Mass (i.e., goiter)

### WARNINGS

None known

### PRECAUTIONS

- This product is intended for use by physicians trained and experienced in retrograde intubation techniques. Standard techniques for retrograde intubation should be employed.
- The 6.0 French set can only be used in placement of endotracheal tubes whose ID is 2.5 mm or larger.
- The 11.0 French set can only be used in placement of endotracheal tubes whose ID is 4 mm or larger.
- The 14.0 French set can only be used in placement of endotracheal tubes whose ID is 5 mm or larger.

### INSTRUCTIONS FOR USE

#### Illustrations

1. After prepping access site in standard fashion, advance either the introducer needle (attached to syringe) or the catheter introducer needle (attached to syringe) in a cephalad direction through cricothyroid membrane and into trachea. Confirm positioning via free-flow aspiration of air into syringe. **(Fig. 1)**
2. If using introducer needle, remove syringe, leaving needle in place. If using catheter introducer needle, remove needle and syringe, leaving catheter in place.
3. Advance "J" end of wire guide through catheter or needle and up into trachea in a cephalad direction until tip of wire guide can be retrieved through the patient's mouth or nose. **NOTE:** The black positioning marks on wire guide should be visible both at access site and where wire guide exits nares or mouth. **(Fig. 2)** This will ensure that enough wire is exposed orally or nasally for control of subsequent catheter introduction.
4. Remove catheter or needle, leaving wire guide in place.
5. Advance radiopaque catheter antegrade over wire guide via mouth or nose and into trachea until tenting is noted at cricothyroid access site. **(Fig. 3)**  
**NOTE:** Alternatively, after initial wire guide positioning, a fiberoptic bronchoscope may be used for direct visualization of endotracheal tube placement. Wire guide may be placed through suction port of fiberoptic scope.
6. With catheter in position, advance endotracheal tube over catheter and into position below the level of the vocal cords. **(Fig. 4)**  
**NOTE:** Always maintain control and position of wire guide during advancement of endotracheal tube.
7. Remove wire guide and catheter from endotracheal tube, advance endotracheal tube to optimal position, and inflate balloon cuff if provided.
8. Verify position and secure in standard fashion. **(Fig. 5)**

### HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or

undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

## REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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**Bloomington, IN 47404 U.S.A.**