
- EN**  
5 **EchoTip® Ultra Ultrasound Needle**  
Instructions for Use
- CS**  
8 **Ultrazvuková jehla EchoTip® Ultra**  
Návod k použití
- DA**  
11 **EchoTip® Ultra ultralydnål**  
Brugsanvisning
- NL**  
14 **EchoTip® Ultra ultrasone naald**  
Gebruiksaanwijzing
- FR**  
17 **Aiguille ultrasonore EchoTip® Ultra**  
Mode d'emploi
- DE**  
20 **EchoTip®-Ultra Ultraschallnadel**  
Gebrauchsanweisung
- EL**  
23 **Βελόνα υπερήχων EchoTip® Ultra**  
Οδηγίες χρήσης
- HU**  
26 **EchoTip® Ultra echogén tű**  
Használati utasítás
- IT**  
29 **Ago ecogenico EchoTip® Ultra**  
Istruzioni per l'uso
- NO**  
32 **EchoTip® Ultra ultralydnål**  
Bruksanvisning
- PL**  
35 **Igła ultrasonograficzna EchoTip® Ultra**  
Instrukcja użycia
- PT**  
38 **Agulha ecoguiada EchoTip® Ultra**  
Instruções de utilização
- RO**  
41 **Ac pentru ecografie Echotip® Ultra**  
Instrucțiuni de utilizare
- ES**  
44 **Aguja de ecografía EchoTip® Ultra**  
Instrucciones de uso
- SV**  
47 **EchoTip® Ultra ultraljudsnål**  
Bruksanvisning

1.


**Plunger**

Píst  
Stempel  
Plunjer  
Piston  
Kolben  
Έμβολο  
Dugattyú  
Stantuffo  
Stempel  
Tłoczek  
Émbolo  
Piston  
Émbolo  
Kolv

**Stopcock "open"**

Uzavírací kohout v otevřené poloze  
Hane "åben"  
Afsluitkraan 'open'  
Robinet ouvert  
Absperrhahn geöffnet  
Στρόφιγγα ανοικτή  
Elzárócsap „nyitott” állásban  
Rubinetto aperto  
Stoppekran åpen  
Kranik w pozycji otwartej  
Torneira de passagem "aberta"  
Robinet „deschis"  
Llave de paso «abierta»  
Kran "öppen"

2.

### Sliding Sheath Adjuster

Posuvné nastavovací zařizení sheathu  
 Glidende sheath-justering  
 Schuifafstelling van de sheath  
 Ajusteur de gaine coulissante  
 Verschiebbarer Schleuseneinsteller  
 Προσαρμογέας του συρόμενου θηκاريوύ  
 Elcsúszatható hüvelyhossz-beállító  
 Regolatore della guaina scorrevole  
 Glidende hylsejusteringsmekanisme  
 Przesuwany regulator koszulki  
 Regulator deslizante da bainha  
 Ajustor glisant al tecii  
 Ajustador deslizante de la vaina  
 Flyttbar justerare till höljet

### Sheath Reference Mark

Referenční značka sheathu  
 Referencemærke på sheath  
 Referentiemarkering van de sheath  
 Repère de référence sur la gaine  
 Schleusenbezugsmarke  
 Σήμανση αναφοράς θηκاريوύ  
 Hüvelyhossz referencialejele  
 Indicatore di riferimento della guaina  
 Referansemerke på hylsen  
 Znacznik referencyjny wysunięcia koszulki  
 Marca de referència da bainha  
 Marcaj de referință al tecii  
 Marca de referenția de la vaina  
 Höljets referensmarkering

### Needle Handle



Rukojeť jehly  
 Nålegreb  
 Naaldhandvat  
 Poignée de l'aiguille  
 Nadelgriff  
 Λαβή βελόνας  
 Tű fogantyúja  
 Impugnatura dell'ago  
 Nålehåndtak  
 Uchwyt igły  
 Punho da agulha  
 Mănerul acului  
 Mango de la aguja  
 Nålhandtag

### Safety Ring

Pojistný prstenec  
 Sikkerhedsring  
 Veiligheidsring  
 Bague de sécurité  
 Sicherheitsring  
 Δακτύλιος ασφάλειας  
 Biztonsági gyűrű  
 Anello di sicurezza  
 Sikkerhetsring  
 Pierścień bezpieczeństwa  
 Anel de siguranță  
 Anillo de seguridad  
 Säkerhetslåsring

### Stylet Hub

Nástavec stiletu  
 Stiletmuffe  
 Aanzetstuk van het stilet  
 Embase du stylet  
 Mandrinansatz  
 Ομφαλός σπειλεού  
 Mandrin végidoma  
 Pomello del mandrino  
 Stilettfeste  
 Glowica mandrynu  
 Conector do estilete  
 Acoperitoarea stiletului  
 Conector del estilete  
 Stilettfattning

## INTENDED USE

This device is used to sample targeted submucosal gastrointestinal lesions through the accessory channel of an ultrasound endoscope.

## NOTES

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

Do not use this device for any purpose other than stated intended use.

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Store in a dry location, away from temperature extremes.

Use of this device restricted to a trained healthcare professional.

## CONTRAINDICATIONS

Those specific to primary endoscopic procedure to be performed in gaining access to desired site.

Coagulopathy.

## POTENTIAL COMPLICATIONS

Those associated with gastrointestinal endoscopy include, but are not limited to: perforation, hemorrhage, aspiration, fever, infection, allergic reaction to medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

## PRECAUTIONS

Refer to package label for minimum channel size required for this device.

Needle must be retracted into sheath and thumbscrew on safety ring must be locked to hold needle in place prior to introduction, advancement or withdrawal of device. Failure to retract needle may result in damage to endoscope.

## SYSTEM PREPARATION

### Illustrations

1. Examine syringe. It has two plunger locks that must be depressed to advance plunger. Tip of syringe has a Luer slip fitting with a stopcock on the side port. Air can be exchanged when stopcock is in "open" position, aligned with syringe. (See **fig. 1**)
2. Prepare syringe as follows:
  - a. With stopcock in "open" position, depress plunger locks and fully advance plunger into syringe.
  - b. Turn stopcock 90° to "closed" position.

- c. Pull plunger back until it is locked into place at desired setting, creating a suction.
  - d. Set prepared syringe aside until aspiration biopsy is desired.
3. Advance device into ultrasound endoscope to determine preferred sheath length. To adjust length, loosen thumbscrew lock on sliding sheath adjuster and slide until preferred length is attained. **Note:** Reference mark for sheath length will appear in sliding sheath adjuster window. (**See fig. 2**) Tighten thumbscrew on sliding sheath adjuster to maintain preferred sheath length.

## INSTRUCTIONS FOR USE

## Illustrations

1. Identify desired biopsy site by endoscopic ultrasound.
2. Introduce ultrasound needle into accessory channel and advance device in small increments until Luer lock fitting at base of sliding sheath adjuster meets fitting on accessory channel. (**See fig. 3**)
3. Attach device to accessory channel port by rotating device handle until fittings are connected.
4. With ultrasound endoscope and device straight, adjust needle to desired length by loosening thumbscrew on safety ring, and advancing it until desired reference mark for needle advancement appears in the window of safety ring. (**See fig. 4**) Tighten thumbscrew to lock safety ring in place. **Note:** Number in safety lock ring window indicates extension of needle in centimeters.  
**Caution:** During needle adjustment or extension, ensure device has been attached to accessory channel. Failure to attach device prior to needle adjustment or extension may result in damage to endoscope.
5. While maintaining position of ultrasound endoscope, extend needle by advancing needle handle to pre-positioned safety ring.
6. Advance needle into lesion.
7. Remove stylet from needle by gently pulling back on plastic hub seated in metal fitting of needle handle. Preserve stylet for use if additional cell collection is desired.
8. Attach Luer lock fitting of previously prepared syringe securely onto fitting on handle.
9. Turn stopcock to "open" position aligned with syringe, allowing negative pressure in syringe to aspirate cells. Gently move handle in small increments back and forth within biopsy site. **Note:** Do not remove needle from biopsy site during cell collection.
10. Upon completion of cell collection, retract needle into sheath by pulling back on handle shaft. Lock thumbscrew on safety ring at 0 cm mark.
11. Disconnect Luer lock fitting from accessory channel and withdraw entire device from ultrasound endoscope.

12. Unlock thumbscrew on safety ring and extend needle. Exchange locking syringe for an empty syringe and attach new syringe to Luer lock fitting on handle.
13. Push in plunger to expel aspirate specimen, then prepare per institutional guidelines.
14. Aspirate remaining can be recovered for examination by flushing device.
15. For additional cell collection from same lesion, gently reinsert stylet into metal fitting on handle. **Note:** Prior to reinserting stylet, wipe with saline or sterile water. While supporting sheath at Luer lock fitting, advance stylet in small increments until stylet hub is engaged in fitting.
16. Additional samples may be obtained by repeating Step 2 of “System Preparation” through Step 14 of the “Instructions for Use.”

**Upon completion of procedure, dispose of device per institutional guidelines for biohazardous medical waste.**

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© COPYRIGHT COOK 2019 ©
**COOK IRELAND LTD.**