**EN**  
3

**Percutaneous Nephrostomy Set**

Instructions for Use

**DA**  
4

**Perkutant nefrostomisæt**

Brugsanvisning

**DE**  
5

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**EL**  
6

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Bruksanvisning

## PERCUTANEOUS NEPHROSTOMY SET

**CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).**

### DEVICE DESCRIPTION

The Percutaneous Nephrostomy Set includes the following:

- 22 gage, 15 cm long Chiba needle
- 18 gage, 15 cm long trocar needle
- .038 inch (0.97 mm) wire guide, 60 cm long with 3 mm Safe-T-J® tip
- 7, 8, and 9 French radiopaque dilators, 20 cm long
- 8.3 French radiopaque pigtail catheter, 30 cm long
- Fixation device with pull tie
- Connecting tube, 30 cm long with one-way stopcock and universal connector

### INTENDED USE

The Percutaneous Nephrostomy Set is intended for external urine drainage from the renal pelvis.

### CONTRAINDICATIONS

None known

### WARNINGS

Manipulation of this product requires fluoroscopic or ultrasound guidance.

### PRECAUTIONS

- The product is intended for use by physicians trained and experienced in diagnostic and interventional techniques. Standard techniques for placement of percutaneous drainage catheters should be employed.
- The potential effects of phthalates on pregnant/nursing women or children have not been fully characterized and there may be concern for reproductive and developmental effects.

### INSTRUCTIONS FOR USE

1. With patient prone, preliminary localization of the renal pelvis by plain film, delayed pyelogram, retrograde pyelogram, ultrasound or CT should be performed.

#### Skinny Needle Antegrade Pyelogram

2. After prepping skin, under fluoroscopic control, a 22 gage needle puncture is made two finger breadths lateral to the lateral margin of the paraspinous muscles and under the 12th rib.

**NOTE:** The puncture should be made as straight down as possible.

**NOTE:** Ultrasound could be used to guide the puncture, then fluoroscopy for contrast injection. Real-time ultrasound with the patient on a fluoroscopy table is most convenient.

3. An antegrade pyelogram is performed to provide visualization of the target collecting system, gauge depth, and distend the system for subsequent puncture.

**NOTE:** A urine specimen may be taken directly from the renal pelvis before contrast injection if desired.

4. The needle is withdrawn by grasping the needle at the skin so the depth of the collecting system can then be gauged. Alternatively, in the case of small, partially obstructed systems, the needle can be left in place and contrast infused during the next step. Alternatively, if needle position is suitable for nephrostomy, the 18 gage introducing needle can be placed next to it to the same depth.

#### Introducing Needle

5. With rubber-shod forceps placed on the 18 gage needle at the depth previously gauged, the renal pelvis is entered as before. Free flow of urine indicates proper positioning.
6. A .038 inch, 3 mm Safe-T-J® wire guide is introduced and advanced under fluoroscopic control until it is well inside the renal pelvis or down the ureter (3 inches).

#### Dilatation and Tube Placement

7. The introducing needle is withdrawn and the tract is dilated using the appropriate sequence of dilators passed over the wire guide, maintaining wire guide position utilizing fluoroscopic visualization.
8. Prior to introduction, the pigtail should be uncoiled. The nephrostomy catheter is introduced over the wire guide to a point where the pigtail portion is positioned in the renal pelvis, infundibulum or ureter as appropriate. To prevent accidental advancement or withdrawal of the wire guide, the location of its tip should be continuously monitored fluoroscopically.
9. Once catheter is properly positioned, the wire guide should be removed.

10. Secure the catheter to the skin with the disc and/or suture, then bandage well.

## HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

## REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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Bloomington, IN 47404 U.S.A.
Cook Medical Europe Ltd