Biodesign®

TISSUE GRAFT

FP0005-04K

COOK®

MEDICAL

Manufacturer

Temperature limit

Use-by-date

Do not re-use

Date of manufacture

Attention, see instructions for use

Keep dry

STERILE EO Sterilized using ethylene oxide

MR Magnetic resonance safe

MANUFACTURER

COOK BIOTECH
INCORPORATED
1425 Innovation Place
West Lafayette, IN 47906 U.S.A.
Phone: 812 339-2235
Toll Free: 800 457-4500
Toll Free Fax: 800 554-8335

EC REPRESENTATIVE

COOK IRELAND
O'Halloran Road
National Technological Park
Limerick, IRELAND
Phone: +353 61 334440

AUSTRALIAN SPONSOR

WILLIAM A. COOK
AUSTRALIA PTY. LTD.
Brisbane Technology Park
95 Brandl Street
Eight Mile Plains
Brisbane, QLD 4113 Australia
Phone: +61 7 38 41 11 88

www.cookmedical.com
© COOK BIOTECH, INC. 2019

February 2019

ENGLISH

BIODESIGN® TISSUE GRAFT

DEVICE DESCRIPTION:

The Cook® Biodesign® Tissue Graft is a dried single or multi-layered small intestinal submucosa (SIS) sheet. The graft is used to reinforce soft tissue and can be cut to size to accommodate the patient's anatomy. The Tissue Graft is provided sterile for single use only.

INTENDED USE

The Cook® Biodesign® Tissue Graft is intended for implantation to reinforce soft tissue. The graft is supplied sterile in peel-open packages and is intended for one-time use.

Rx ONLY This symbol means the following:

CAUTION: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

1-LAYER TISSUE GRAFT This symbol means the following: 1-Layer Tissue Graft

4-LAYER TISSUE GRAFT This symbol means the following: 4-Layer Tissue Graft

This product is intended for use by trained medical professionals.

CONTRAINDICATIONS

This graft is derived from a porcine source and should not be used for patients with known sensitivity to porcine materials.

PRECAUTIONS

- This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.
- Do not resterilize.** Discard all open and unused portions.
- Graft is sterile if the package is dry, unopened and undamaged. Do not use if the package seal is broken.
- Discard the graft if mishandling has caused possible damage or contamination, or if the graft is past its expiration date.
- Single-layer graft should not be used in applications requiring high strength.
- The graft may not have sufficient strength to support stresses encountered in some ventral hernias or large-area, body-wall repairs.
- Ensure that the graft is rehydrated prior to suturing or stapling.
- Patients undergoing radiation therapy may not experience normal wound healing.
- Graft performance has not been evaluated with suture spacing greater than 2 mm.**
- Ensure that all layers of the graft are secured when suturing or stapling.**

POTENTIAL COMPLICATIONS

Complications that can occur with the use of any surgical graft material may include, but are not limited to:

- Inflammation
- Allergic reaction
- Seroma formation
- Abscess
- Adhesion
- Fistula formation
- Induration
- Erosion
- Infection
- Tissue trauma
- Recurrence of tissue defect
- Delayed or failed incorporation
- Migration
- Extrusion
- Fever
- Pain
- Bleeding

If any of the following conditions occur and cannot be resolved, device removal should be considered:

- Infection
- Acute or chronic inflammation (initial application of surgical graft materials may be associated with transient, mild, localized inflammation)
- Allergic reaction
- Seroma formation

STORAGE

This graft should be stored in a clean, dry location at room temperature.

STERILIZATION

This graft has been sterilized with ethylene oxide, and should not be resterilized.

INSTRUCTIONS FOR USE

Required Materials

- A sterile dish (kidney dish or other bowl)
- Sterile forceps
- Rehydration fluid: at least 50 ml of room temperature sterile saline or sterile lactated Ringer's solution for each graft sheet

NOTE: Handle the graft using aseptic technique. Minimize contact with latex gloves.

- Using aseptic technique, remove the graft inner pouch from its outer pouch, and place the inner pouch in the sterile field.
- Open the inner pouch carefully, and aseptically remove the graft sheet with the sterile forceps.
- Place the graft sheet into the sterile dish in the sterile field. (Multiple graft sheets may be rehydrated simultaneously in the same dish.)
- Add to the dish at least 50 mL of the rehydration fluid for each graft sheet.
- Rehydrate the graft in the rehydration fluid until the desired handling characteristics are achieved. The graft should not be rehydrated for longer than five minutes.

- Prepare the graft site using standard surgical techniques.
- Using aseptic technique, trim the graft sheet to fit the site, providing a small allowance for overlap. **NOTE:** An alternative method is to cut the graft sheet to size prior to rehydration. If this method is selected, be sure to rehydrate the graft sheet prior to suturing or stapling it into place. See step 5.
- Using aseptic technique, transfer the graft sheet to the graft site and suture or staple into place, avoiding excess tension. **NOTE:** Surgical experience indicates that suturing or stapling graft sheets with close tissue approximation produces better outcomes.
- Complete the standard surgical procedure.
- Discard any unused portions of the graft sheet according to institutional guidelines for medical waste.

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