- EN**  
2 **Colon Decompression Set**  
Instructions for Use
- CS**  
4 **Souprava pro kolonoskopickou dekompresi**  
Návod k použití
- DA**  
6 **Colondekompressionssæt**  
Brugsanvisning
- NL**  
8 **Colondecompressieset**  
Gebruiksaanwijzing
- FR**  
10 **Set pour décompression du côlon**  
Mode d'emploi
- DE**  
12 **Kolon-Dekompressionsset**  
Gebrauchsanweisung
- EL**  
14 **Σετ αποσυμπίεσης παχέος εντέρου**  
Οδηγίες χρήσης
- IT**  
16 **Set di decompressione del colon**  
Istruzioni per l'uso
- PL**  
18 **Zestaw do odbarczania okrężnicy**  
Instrukcja użycia
- PT**  
20 **Conjunto de descompressão do cólon**  
Instruções de utilização
- ES**  
22 **Equipo de descompresión del colon**  
Instrucciones de uso
- SV**  
24 **Set för kolondekompression**  
Bruksanvisning

**INTENDED USE**

This device is used for treatment of acute non-toxic megacolon, pseudo-obstruction (Ogilvie's syndrome) and colonic strictures.

**NOTES**

Set includes a drainage catheter, guiding catheter and .035 inch wire guide.

Do not use this device for any purpose other than stated intended use.

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Store in a dry location, away from temperature extremes.

Use of this device restricted to a trained healthcare professional.

**CONTRAINDICATIONS**

Those for colonoscopy include, but are not limited to: fulminant colitis, possible existing perforation, acute severe diverticulitis, presence of barium or inadequate bowel preparation.

Those for colon decompression include, but are not limited to: ischemic bowel, recent rectal surgery, organic colon obstruction, recent surgical anastomosis, acute myocardial infarction.

**POTENTIAL COMPLICATIONS**

Those associated with colonoscopy include, but are not limited to: perforation, hemorrhage, fever, infection, allergic reaction to medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

Others associated with colon decompression include, but are not limited to: clogging of the tube, dislodgment of the tube as colonoscope is withdrawn, perforation.

**PRECAUTIONS**

Refer to package label for minimum channel size required for this device.

Careful history and physical examination must be performed prior to tube insertion in order to determine underlying cause of pseudo-obstruction.

Colon should not be prepped unless there is an inordinate amount of stool in rectum.

**INSTRUCTIONS FOR USE**

1. Advance colonoscope, insufflating minimal air, until it has reached cecum or as far in colon as possible. **Note:** Suction must be performed intermittently throughout procedure to help decompress colon and avoid further colonic distention.

2. Fluoroscopically observe location of colonoscope tip.
3. Insert .035 inch wire guide into biopsy channel of colonoscope until it exits tip.
4. Carefully withdraw colonoscope, making certain wire guide remains stationary within colon by periodic fluoroscopic monitoring.
5. Disassemble colon decompression tube and pre-loaded 6 FR guiding catheter.
6. Flush colon decompression tube and 6 FR guiding catheter with sterile water or lubricate with water-soluble lubricant, then reassemble and advance into colon over pre-positioned wire guide. Use fluoroscopic control to avoid migration of wire guide.
7. Once tube is advanced to end of wire guide, withdraw both guiding catheter and wire guide, leaving decompression tube in place.
8. Secure external portion of decompression tube to patient.

**Upon completion of procedure, dispose of device(s) per institutional guidelines for biohazardous medical waste.**

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© COPYRIGHT COOK 2017 ©
**COOK IRELAND LTD.**