**COOK®**

MEDICAL

**CE**

0123

**EN  
2****Disposable Snare****CS  
3****Očko na jedno použití****DA  
4****Engangsslynge****NL  
6****Disposable lus****FR  
8****Anse jetable****DE  
9****Einmalschlinge****EL  
11****Αναλώσιμος βρόχος****HU  
13****Egyszer használatos hurok****IT  
14****Ansa monouso****NO  
16****Slynge til engangsbruk****PL  
17****Pętla jednorazowego użycia****PT  
19****Laço descartável****ES  
20****Lazo desechable****SV  
22****Engångssnara**

\* 1 8 8 9 0 / 0 7 1 9 \*

**INTENDED USE**

This device is used with an electrosurgical unit for endoscopic polypectomy.

**NOTES**

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

If the package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Notify Cook for return authorization.

Do not use this device for any purpose other than the stated intended use.

Store in a dry location, away from temperature extremes.

Use of this device restricted to a trained healthcare professional.

**CONTRAINDICATIONS**

Contraindications include those specific to the primary endoscopic procedure to be performed in gaining access to the polypectomy site.

Contraindications to polypectomy include, but are not limited to: coagulopathy and insufficiently prepped bowel.

**POTENTIAL COMPLICATIONS**

Potential complications associated with gastrointestinal endoscopy include, but are not limited to: perforation, hemorrhage, aspiration, fever, infection, hypotension, allergic reaction to medication, respiratory depression or arrest, cardiac arrhythmia or arrest.

Complications which can occur with endoscopic polypectomy include, but are not limited to: transmural burns, thermal injury to the patient, explosion.

**PRECAUTIONS**

Refer to package label for minimum channel size required for this device.

Before using this device, follow recommendations provided by electrosurgical unit manufacturer to ensure patient safety through proper placement and utilization of patient return electrode. Ensure a proper path from patient return electrode to electrosurgical unit is maintained throughout the procedure.

Polypectomy should be performed under direct endoscopic visualization.

When using AcuSnare® Needle-Tip, needle tip of snare must be fully retracted into sheath prior to introduction, advancement, or withdrawal from endoscope. Failure to fully retract snare and needle may result in damage to endoscope.

**INSTRUCTIONS FOR USE**

1. Inspect active cord. Cord must be free of kinks, bends, breaks and exposed wires to allow for accurate transfer of current. If an abnormality is noted, do not use active cord.
2. With electrosurgical unit off, prepare equipment. Securely connect active cord to device handle and electrosurgical unit. Active cord fittings should fit snugly into both device handle and electrosurgical unit. Following instructions from electrosurgical unit manufacturer, position patient return electrode and connect it to electrosurgical unit.
3. Fully retract and extend snare to confirm smooth operation of device. **Note:** If using an AcuSnare, slide adjustable marker, located in handle, to establish a reference point indicating full retraction of snare into sheath and to set up reference points for establishing thickness of tissue being excised.
4. When polyp is in endoscopic view, introduce sheath and retracted snare into endoscope accessory channel. **Caution:** To ensure patient safety, power to electrosurgical unit should remain off until snare is properly positioned around the polyp.
5. Advance device, in small increments, until endoscopically viewed exiting endoscope.

6. Advance snare wire out of sheath and position it around polyp to be removed. **Warning:** When applying current, tissue must be isolated from surrounding mucosa. Failure to isolate tissue may cause fulguration of normal mucosa and/or perforation. Contact of snare wire with endoscope during electrosurgery may cause grounding, which could result in injury to patient and/or operator as well as damage to endoscope and/or snare wire. **Note:** When using AcuSnare Needle Tip, needle tip can be used to stabilize snare and aid in positioning.
7. Following electrosurgical unit manufacturer's instructions for settings, verify desired settings and activate electrosurgical unit. **Note:** Maximum rated input voltage for this device is 2 kVp-p for cut mode and 5 kVp-p for coagulation mode.
8. Proceed with polypectomy.
9. Upon completion of polypectomy, turn electrosurgical unit off. Retract snare into sheath and remove device from endoscope.
10. Retrieve polyp and prepare specimen per institutional guidelines.

**Upon completion of procedure, disconnect active cord from device handle and dispose of this device per institutional guidelines for biohazardous medical waste.** Disconnect active cord from electrosurgical unit. Wipe active cord with a damp cloth to remove all foreign matter. Store in a loose coil. **Note:** Wrapping active cord tightly may damage the device.

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**Wilson-Cook Medical, Inc.**
**Cook Medical Europe Ltd.**