COOK®

MEDICAL

EN
2

Achalasia Balloon

CS
3

Balónek pro pacienty s achalazií

DA
5

Achalasia ballon

NL
7

Achalasieballon

FR
10

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DE
12

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* 1 8 9 3 1 / 0 6 1 9 *

INTENDED USE

This device is used to dilate the lower esophageal sphincter in patients with achalasia. The device is indicated for adult use only.

NOTES

Do not use this device for any purpose other than the stated intended use.

Store in a dry location, away from temperature extremes.

Use of this device is restricted to a trained healthcare professional.

This device is supplied non-sterile and may not be sterilized by any method.

CONTRAINDICATIONS

Those specific to upper GI endoscopy.

Those specific to dilation include, but are not limited to: uncooperative patient, asymptomatic rings or strictures, inability to advance the balloon through the strictured area, coagulopathy, known or suspect perforation, severe inflammation or scarring near the dilatation site.

POTENTIAL COMPLICATIONS

Those associated with upper GI endoscopy include, but are not limited to: perforation, hemorrhage, aspiration, fever, infection, allergic reaction to medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

PRECAUTIONS

Inflate the balloon with air only. **Note:** Inflating with liquid can increase the risk of balloon damage, aspiration and/or perforation.

Do not pre-inflate the balloon.

WARNINGS

This device is supplied non-sterile, and therefore should not be used in a sterile field.

Do not exceed the recommended balloon pressure as specified on the molded junction and package label.

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

If the package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

EQUIPMENT REQUIRED

- Achalasia Balloon
- Three-way stopcock
- Compatible Wire Guide (not included)
- DS-60CC Inflation Syringe (not included)
- Inflation Device Handle (not included)

INSTRUCTIONS FOR USE

1. Introduce the endoscope into the patient and advance to the proximal margin of the stricture.
2. Introduce a wire guide, floppy tip first, into the accessory channel and advance in short increments until it is positioned well into the stomach.
3. Remove the endoscope, leaving the wire guide in place.
4. Attach a 60cc inflation syringe to the inflation device.
5. Attach the syringe extension line to the stopcock and then to the inflation hub of the balloon catheter.
6. With the plunger at the bottom of the syringe, turn the stopcock lever to the "open" position (pointing toward the side port) and then retract the plunger of the syringe back to aspirate all the residual air from the balloon.
7. While maintaining pressure, turn the stopcock lever to the "off" position (pointing toward the balloon).
8. Under fluoroscopic visualization, introduce the catheter over the wire guide, and advance in short increments until the dilator is in the proper position. **Caution:** Do not introduce the balloon without a wire guide as this may increase the risk for perforation.
9. Fluoroscopically confirm that the balloon remains in the proper position, then turn the stopcock lever 90° to the "open" position (pointing toward the side port).
10. Under fluoroscopic visualization, advance the plunger of the syringe to rapidly inflate the balloon with air up to the recommended pressure (as specified on the molded junction and the package label).
Note: Balloon inflation may be monitored endoscopically by placing the endoscope next to the balloon catheter.
11. When the balloon reaches the appropriate pressure, turn the stopcock lever to the "off" position (pointing toward the balloon) while maintaining balloon pressure. **Note:** In order to reach the recommended pressure, steps 10 & 11 may have to be repeated.
12. Proceed with dilation.
13. Turn the stopcock lever to the "open" position (pointing toward the side port) and deflate the balloon.
14. Carefully remove the deflated balloon, wire guide, and endoscope if applicable.

Upon completion of the procedure, dispose of the device per institutional guidelines for biohazardous medical waste.

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© 2019 Cook Medical
Wilson-Cook Medical, Inc.
Cook Medical Europe Ltd.