COOK®

MEDICAL

CE

0123

EN
3

Biliary Dilation Catheter

CS
4

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DA
5

Galdedilationskateter

NL
6

Galdilatatiekatheter

FR
7

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DE
9

Gallengang-Dilatationskatheter

EL
10

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

INTENDED USE

This device is used to dilate biliary strictures.

This device is indicated for adult use only.

NOTES

Store in a dry location, away from temperature extremes.

CONTRAINDICATIONS

Contraindications include those specific to ERCP.

WARNINGS

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Do not use this device for any purpose other than stated intended use.

Use of this device restricted to a trained healthcare professional.

POTENTIAL COMPLICATIONS

Potential complications associated with ERCP include, but are not limited to: pancreatitis, cholangitis, aspiration, perforation, hemorrhage, infection, sepsis, allergic reaction to contrast or medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

PRECAUTIONS

Refer to package label for minimum channel size required for this device.

If unexpected resistance is met during the procedure, do not continue advancement of catheter without first determining the cause of the unexpected resistance and taking remedial action.

If bleeding occurs during the procedure, stop the procedure immediately and take remedial action.

INSTRUCTIONS FOR USE (See fig. 1)**I. If using Zip port and a pre-positioned short wire.**

1. Advance tip of catheter over pre-positioned wire guide ensuring wire guide exits catheter at zip port, then unlock short wire from Wire Guide Locking Device.
2. Introduce catheter into endoscope accessory channel and relock wire guide. Then advance catheter until it is endoscopically visualized exiting endoscope.
3. Proceed with cannulation.
4. Under fluoroscopic monitoring, advance catheter to appropriate position and dilate stricture or papillary orifice. **Note:** Dilation is achieved once distal radiopaque band has passed through stricture or papillary orifice.

Note: Previously placed wire guide may be left in position, in order to facilitate introduction of other wire guided devices.

If wire guide is to remain in place while device is withdrawn, utilize following steps:

5. Upon completion of dilation, remove the catheter from accessory channel until resistance is achieved. Unlock wire guide from Wire Guide Locking Device and perform wire guide exchange. Relock wire guide and remove catheter from wire guide.

II. If using Proximal Wire Port (PWP) and a pre-positioned long wire guide. (See fig. 1)

Note: For best results, wire guide should be kept wet.

1. Remove stylet.
2. Advance tip of catheter over pre-positioned wire guide until endoscopically visualized exiting endoscope.

REFER TO STEPS 3-4 IN "SECTION I" THEN RESUME WITH STEP 3 BELOW.

3. Remove the catheter using standard long wire exchange technique.

Upon completion of procedure, dispose of device per institutional guidelines for biohazardous medical waste.

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© 2019 Cook Medical
Wilson-Cook Medical, Inc.
Cook Medical Europe Ltd.