EN

Biliary Cytology Brush

CS

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DA

Cytologibørste til galdevejene

NL

Biliaire cytologieborstel

FR

Brosse de cytologie biliaire

DE

Gallengang-Zytologiebürste

EL

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ENGLISH

INTENDED USE

This device is used for collection of cells in the biliary system.

NOTES

Inspect the device with particular attention to kinks, bends or breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Do not use this device for any purpose other than the stated intended use.

Store in a dry location, away from temperature extremes.

Use of this device restricted to a trained healthcare professional.

CONTRAINDICATIONS

Contraindications include those specific to ERCP and any procedure to be performed in conjunction with cell collection.

POTENTIAL COMPLICATIONS

Potential complications associated with ERCP include, but are not limited to: pancreatitis, cholangitis, aspiration, perforation, hemorrhage, infection, sepsis, allergic reaction to contrast or medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

PRECAUTIONS

Refer to package label for minimum channel size required for this device.

The wire guide diameter and the inner lumen of the wire-guided device must be compatible.

Injection of contrast during ERCP must be monitored fluoroscopically. Overfilling of the pancreatic duct may cause pancreatitis.

INSTRUCTIONS FOR USE

I. If using 8 FR brush with Intra Ductal Exchange (IDE) port and pre-positioned short wire. (See Fig. 1)

1. Retract brush completely into sheath.
2. Unlock short wire from Wire Guide Locking Device and advance tip of device over pre-positioned wire guide ensuring that wire guide exits catheter at IDE port.
3. Introduce device into endoscope accessory channel and relock wire guide. Then advance in short increments until device is endoscopically visualized exiting scope.
4. Advance device under fluoroscopic monitoring to desired position.
5. Advance cytology brush out of sheath, then brush area to be sampled to obtain adequate cellular material.
6. After obtaining specimen, retract brush into sheath.

Note: Previously placed wire guide may be left in position, in order to facilitate introduction of other wire guided devices.

If wire guide is to remain in place while device is withdrawn, utilize following steps:

7. Prior to removing device, utilize reference marks on catheter to ensure that IDE port is within ductal system.
8. Fluoroscopically visualize radiopaque band at IDE port. Retract wire guide until radiopaque distal tip of wire guide passes band; a disengagement from wire guide lumen will occur.
9. Advance disengaged wire guide to maintain ductal access.
10. Lock wire guide into Wire Guide Locking Device and remove device from endoscope accessory channel.
11. Prepare specimen per institutional guidelines.

II. If using Proximal Wire Port (PWP) on 8 FR brush and pre-positioned long wire.

Note: For best results, wire guide should be kept wet.

1. Retract brush completely into sheath.
2. Advance tip of device over pre-positioned wire guide ensuring that wire guide exits catheter at PWP.

REFER TO STEPS 3-6 IN "SECTION I", THEN RESUME WITH STEP 3 BELOW:

3. Remove device using standard long wire exchange technique.
4. Prepare specimen per institutional guidelines.

III. If using 6 FR brush with Intra Ductal Exchange (IDE) port and pre-positioned short wire.

1. Retract brush into sheath ensuring the distal tip is proximal to IDE port.
2. Unlock short wire from Wire Guide Locking Device and advance tip of device over pre-positioned wire guide ensuring that wire guide exits catheter at IDE port.
3. Introduce device into endoscope accessory channel and relock wire guide. Then advance in short increments until device is endoscopically visualized exiting scope.
4. Advance device under fluoroscopic monitoring to desired position.

Note: Previously placed wire guide may be left in position in order to facilitate introduction of other wire guided accessories.

5. Utilize reference marks on catheter to ensure that IDE port is within the ductal system.
6. Fluoroscopically visualize radiopaque band at IDE port. Retract wire guide until radiopaque distal tip of wire guide passes band; a disengagement from wire guide lumen will occur.
7. Advance disengaged wire guide and lock into Wire Guide Locking Device to maintain ductal access.
8. Advance cytology brush out of sheath and then brush area to be sampled to obtain adequate cellular material.
9. After obtaining specimen, retract brush into sheath and remove device from endoscope.
10. Prepare specimen per institutional guidelines.

Upon completion of procedure, dispose of device per institutional guidelines for biohazardous medical waste.

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Wilson-Cook Medical, Inc.
Cook Ireland Ltd.