

EN 3	EchoTip® Ultra Fiducial Needle
CS 5	Referenční jehla EchoTip® Ultra
DA 7	EchoTip® Ultra referencenål
NL 9	EchoTip® Ultra naald voor markerplaatsing
FR 12	Aiguille fiduciaire EchoTip® Ultra
DE 14	EchoTip® Ultra Markierungsnadel
EL 17	Βελόνα δεικτών αναφοράς EchoTip® Ultra
HU 19	EchoTip® Ultra referenciamarker-elhelyező tű
IT 22	Ago EchoTip® Ultra per il posizionamento di marker di riferimento
NO 24	EchoTip® Ultra rammemarkørnål
PL 26	Igła znacznika EchoTip® Ultra
PT 29	Agulha para marcadores fiduciais EchoTip® Ultra
ES 31	Aguja fiducial EchoTip® Ultra
SV 33	EchoTip® Ultra referensmarkörnål

1.

- | | | |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <p>1 Sliding Sheath Adjuster
Posuvné nastavovací zařízení sheathu
Glidende sheath-justering
Schuifregeling van de sheath
Ajusteur de gaine coulissant
Verschiebbarer Schleuseneinsteller
Συρόμενος προσαρμογέας θηκαριού
Elcsúsztható hüvelyhossz-beállító
Regolatore della guaina scorrevole
Glidende hylsejusteringsmekanisme
Przesuwany regulator koszulki
Regulador deslizante da bainha
Ajustador deslizante de la vaina
Glidande hylsjusterare</p> | <p>2 Needle Handle
Rukojet' jehly
Nålegræb
Naaldhandvat
Poignée de l'aiguille
Nadelgriff
Λαβή βελόνας
Tű fogantyúja
Impugnatura dell'ago
Nålehåndtak
Uchwyt igły
Punho da agulha
Mango de la aguja
Nålhandtag</p> | <p>3 Stylet
Stylet
Stilet
Stilet
Stylet
Mandrin
Στείλεός
Mandrin
Mandrino
Stilet
Mandryn
Estilete
Estilete
Mandräng</p> |
| <p>4 Sheath Reference Mark
Referenční značka sheathu
Referencemærke på sheath
Referentiemarkering van de sheath
Repère de référence sur la gaine
Schleusenbezugsmarke
Σήμανση αναφοράς θηκαριού
Hüvelyhossz referencijele
Indicatore di riferimento della guaina
Referansemerke på hylsen
Znacznik referencyjny wysunięcia koszulki
Marca de referència da bainha
Marca de referencia de la vaina
Hylsans referensmarkering</p> | <p>5 Safety Ring
Pojistný prsteneč
Sikkerhedsring
Veiligheidsring
Bague de sécurité
Sicherheitsring
Δακτύλιος ασφάλειας
Biztonsági gyűrű
Anello di sicurezza
Sikkerhetsring
Pierścień bezpieczeństwa
Anel de sigurança
Anillo de seguridad
Säkerhetsring</p> | |
| <p>6 Stylet Hub
Nástavec styletu
Stiletmuffe
Aanzetstuk van het stilet
Embase du stylet
Mandrinansatz
Ομφαλός στειλεού
Mandrin végidoma
Pomello del mandrino
Stilettfeste
Złącza mandrynu
Conector do estilete
Conector del estilete
Mandrängfattning</p> | | |

ENGLISH

INTENDED USE

This device is used to implant fiducials under endoscopic ultrasound to radiographically mark soft tissue for future therapeutic procedures.

NOTES

This device is designed for single use only. Attempts to reprocess, sterilize, and/or reuse may lead to device failure and/or transmission of disease.

Do not use this device for any purpose other than stated intended use.

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Store in a dry location, away from temperature extremes.

Use of this device restricted to a trained healthcare professional.

Device may not be used prior to training by manufacturer.

Fiducial migration may occur following placement. Patients who experience fiducial migration may require placement of additional fiducials prior to radiotherapy.

CONTRAINDICATIONS

Contraindications include those specific to the primary endoscopic procedure to be performed in gaining access to the desired target site.

WARNINGS

This needle is not intended for aspiration. Misuse may result in tissue damage.

Do not deploy fiducials under inadequate ultrasound imaging conditions as this may result in the fiducial being inadvertently placed in unintended locations such as the vasculature.

POTENTIAL COMPLICATIONS

Potential complications associated with gastrointestinal endoscopy include, but are not limited to: perforation, hemorrhage, aspiration, fever, infection, allergic reaction to medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

Potential complications from extra-luminal EUS guided procedures may include infection, hemorrhage, perforation and tumor seeding.

PRECAUTIONS

Refer to package label for minimum channel size required for this device.

Needle must be retracted into sheath and thumbscrew on safety ring must be locked to hold needle in place prior to introduction, advancement or withdrawal of device. Failure to retract needle may result in damage to endoscope.

Stylet must not be retracted at any point.

MRI INFORMATION

MR Conditional

Non-clinical testing has demonstrated that the fiducial marker is **MR Conditional** according to ASTM F2503. A patient with this marker can be scanned safely anytime after placement under the following conditions.

Static Magnetic Field

- Static magnetic field of 3 Tesla or less
- Highest spatial magnetic gradient field of 720 Gauss/cm or less

Non-clinical evaluation was conducted in a MR system (Excite, General Electric Healthcare) with a maximum spatial magnetic gradient field of 720 Gauss/cm as measured with a gaussmeter in the position of the static magnetic field pertinent to the patient (i.e., outside of scanner covering, accessible to a patient or individual).

MRI-Related Heating

- 1.5 and 3.0 Tesla Systems: It is recommended to scan in normal operation mode (whole body averaged specific absorption rate (SAR) = 2.0W/kg) ("Normal Operating Mode" is defined as the mode of operation of the MR system in which none of the outputs have a value that cause physiological stress to the patient) for 15 minutes of scanning (i.e., per scanning sequence).

In non-clinical testing, the fiducial marker produced maximum temperature rises of 1.5 °C during 15 minutes of MR imaging (i.e., for one scanning sequence) performed in a 3 Tesla MR System (Excite, HDx, GE Electric Healthcare, Milwaukee, WI, Software 14X.M5) at an MR system reported whole-body-averaged SAR of 2.9 W/kg (associated with a calorimetry measured whole body averaged value of 2.7 W/kg).

Image Artifacts

MR image quality may be compromised if the area of interest is within approximately 26 mm² of the position of the fiducial marker, as found during non-clinical testing using the sequence: T1-weighted, spin echo and gradient echo pulse sequence in a 3.0 Tesla MR system (Excite, HDx, General Electric Healthcare, Milwaukee, WI, Software 14X.M5). Therefore, it may be necessary to optimize MR imaging parameters for the presence of this metallic marker.

For US Patients Only

Cook recommends that the patient register the MR conditions disclosed in this IFU with the MedicAlert Foundation. The MedicAlert Foundation can be contacted in the following manners:

Mail: MedicAlert Foundation International, 2323 Colorado Avenue, Turlock, CA 95382

Phone: 888-633-4298 (toll free) / 209-668-3333 from outside the US

Fax: 209-669-2450

Web: www.medicalert.org

SYSTEM PREPARATION

1. To adjust sheath length, loosen thumbscrew lock on sliding sheath adjuster and slide handle until preferred length is attained (*See fig. 1*). Tighten thumbscrew.
2. Advance stylet until resistance is felt to ensure stylet is flush with proximal fiducial marker.
Note: Excessive pressure on stylet may result in premature deployment.

INSTRUCTIONS FOR USE

1. Identify desired site based on previous findings from endoscopy, radiography and/or CT scans.
2. Slowly introduce needle into accessory channel of endoscope and advance in short increments. Ensure needle is completely retracted and locked in place (*See fig. 2*).
Note: Bends or kinks in needle caused by improper introduction may result in the inability to deploy fiducials.
3. Attach device to endoscope accessory channel port.
4. With endoscope and device straight, advance needle to desired length by loosening thumbscrew on safety ring, and advancing it until desired reference mark appears in window of safety ring (*See fig. 3*). Tighten thumbscrew to lock safety ring in place.
Note: The number in safety lock ring window indicates extension of needle in centimeters.

Caution: Failure to attach device prior to needle adjustment or extension may result in damage to endoscope.

5. Advance needle into target site by advancing needle handle to pre-positioned safety ring.
6. To place a fiducial, depress thumb ring while stabilizing stylet. Continue applying pressure until tactile feedback and ultrasound visualization indicates fiducial has been deployed. Fluoroscopy may also be used to aid in visualizing placement of fiducials. **Note:** Depending on endoscope position and severity of angulation, more than one fiducial may deploy per site.
7. The needle tip can slowly be retracted and advanced into another site to place additional fiducials as needed. **Note:** To prevent premature deployment, retract needle slowly into sheath.
8. Upon completion of procedure, retract needle completely into sheath by pulling back on needle handle and lock thumbscrew on safety ring to hold needle in place.
9. Disconnect Luer lock fitting of device from accessory channel and withdraw entire device from the endoscope.

Upon completion of procedure, dispose of device per institutional guidelines for biohazardous medical waste.

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© 2019 Cook Medical
Wilson-Cook Medical, Inc.
Cook Medical Europe Ltd.