

EN 5	Evolution® Colonic Stent System - Uncovered Instructions for Use
BG 11	Стент система за колон Evolution® - Непокрита Инструкции за употреба
CS 18	Systém kolonického stentu Evolution® - nekrytý Návod k použití
DA 24	Evolution® kolonstentsystem - Udækket Brugsanvisning
NL 30	Evolution® colonstentsysteem - niet-bekleed Gebruiksaanwijzing
ET 37	Käärsoole stendisüsteem Evolution® - ilma katteta Kasutusjuhised
FR 43	Système d'endoprothèse colique Evolution® - non couverte Mode d'emploi
DE 50	Evolution® Kolonstentsystem - Nicht Gecovert Gebrauchsanweisung
EL 57	Σύστημα ενδοπρόσθεσης κόλου Evolution® - ακάλυπτο Οδηγίες χρήσης
HU 64	Evolution® colonsztentrendszer - Fedetlen Használati utasítás
IT 70	Sistema con stent colico Evolution® - Non rivestito Istruzioni per l'uso
LV 77	Evolution® resnās zarnas stenta sistēma - bez apvalka Lietošanas instrukcija
LT 83	„Evolution®“ gaubtinės žarnos stento sistema - nedengta Naudojimo nurodymai
NO 89	Evolution® stentsystem for kolon - utildekket Bruksanvisning
PL 95	System Evolution® do stentowania okrężnicy - Nieosłonięty Instrukcja użycia
PT 102	Sistema de stent de cólon Evolution® - não coberto Instruções de utilização
RO 109	Sistem pentru stentare colonică Evolution® - fără strat de acoperire Instrucțiunile de utilizare
SK 115	Systém kolonického stentu Evolution® - nekrytý Návod na použitie
ES 121	Sistema de stent colónico Evolution® - sin cubierta Instrucciones de uso
SV 128	Evolution® stentsystem för kolon - otäckat Bruksanvisning


1.


2.


3.

4.

5.

EVOLUTION® COLONIC STENT SYSTEM - UNCOVERED

IMPORTANT INFORMATION

Please review prior to use.

INTENDED USE

This device is used for palliative treatment of colonic obstruction or colonic strictures caused by malignant neoplasms, and to relieve large bowel obstruction prior to colectomy in patients with malignant strictures.

STENT DESCRIPTION

This flexible, self-expanding stent is constructed of a single, woven, nitinol wire. The stent foreshortens due to its design. The total length of the stent is indicated by radiopaque markers on the inner catheter, indicating the actual length of the stent at nominal stent diameter.

This Stent is supplied sterile and is intended for single use only.

DELIVERY SYSTEM DESCRIPTION

The stent is mounted on an inner catheter, which accepts a .035 inch wire guide, and is constrained by an outer catheter. An endoscopically and fluoroscopically visible yellow marker defines the proximal end of the stent when constrained in the catheter. A pistol-grip delivery handle allows stent deployment or recapture.

This Delivery System is supplied sterile and is intended for single use only.

NOTES

Do not use this device for any purpose other than stated intended use.

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Use of this device is restricted to a trained healthcare professional.

Store in a dry location, away from temperature extremes.

CONTRAINDICATIONS

Those specific to GI endoscopy and any procedure to be performed in conjunction with stent placement.

Additional contraindications include, but are not limited to: enteral ischemia, suspected or impending perforation, intra-abdominal abscess/perforation, inability to pass wire guide or stent through obstructed area, patients for whom endoscopic procedures are contraindicated, significant coagulopathy, benign disease.

POTENTIAL COMPLICATIONS

Those associated with GI endoscopy include, but are not limited to: perforation, hemorrhage, fever, infection, allergic reaction to medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

Additional complications include, but are not limited to: intestinal perforation, pain, inadequate Stent expansion, stent misplacement and/or migration, tumor ingrowth or overgrowth, stent occlusion, ulcerations, pressure necrosis, erosion of the luminal mucosa, septicemia, foreign body sensation, bowel impaction, diarrhea, constipation, peritonitis, symptoms of tenesmus or urgency/incontinence, death (other than due to normal disease progression).

PRECAUTIONS

Refer to product package label for the minimum channel size required for this device.

A complete diagnostic evaluation must be performed prior to use to determine proper stent size.

If wire guide or stent cannot advance through obstructed area, do not attempt to place stent.

Stent should be placed endoscopically with fluoroscopic monitoring.

The stent should only be placed with the Cook delivery system, which is provided with each stent.

This device is intended for palliative treatment only. Alternate methods of therapy should be investigated prior to placement.

After stent placement, alternative methods of treatment such as chemotherapy and radiation should not be administered as this may

increase risk of stent migration due to tumor shrinkage, stent erosion, and/or mucosal bleeding.

Long-term patency of this device has not been established. Periodic evaluation is advised.

WARNINGS

The stent is not intended to be removed or repositioned after stent placement and is considered a permanent implant. Attempts to remove or reposition stent after placement may cause damage to surrounding tissue or mucosa. Stent cannot be retrieved after the deployment threshold has been passed. Corresponding marks on outer catheter and delivery handle indicate when threshold has been passed.

This stent contains nickel, which may cause an allergic reaction in individuals with nickel sensitivity.

To minimize pain and tenesmus, the stent end nearest to the anal canal/ anus should be placed 2 cm above the anal canal or 6 cm from the anus.

The device should be used with caution and only after careful consideration in patients with:

- Patients with radiation colitis or proctitis.
- Patients with elevated bleeding times, coagulopathies.

MRI INFORMATION

This symbol means the device is MR Conditional.

Non-clinical testing has demonstrated that the Evolution® Colonic **Stent** is MR Conditional according to ASTM F2503. A patient with this **stent** can be scanned safely anytime after placement under the following conditions

- Static magnetic field of 3.0 Tesla or less.
- Spatial magnetic gradient field of 720 Gauss/cm. Non-clinical evaluation was conducted in a MR system (Excite, General Electric Healthcare) with a maximum spatial magnetic gradient field of 720 Gauss/cm as measured with a gaussmeter in the position of the static

magnetic field pertinent to the patient (i.e., outside of scanner covering, accessible to a patient or individual).

MRI-Related Heating

- 1.5 and 3.0 Tesla Systems: It is recommended to scan in normal operation mode (whole body averaged specific absorption rate (SAR) = 2.0W/kg) ("Normal Operating Mode" is defined as the mode of operation of the MR system in which none of the outputs have a value that cause physiological stress to the patient) for 15 minutes of scanning (i.e., per scanning sequence).

1.5 Tesla Temperature Rise

In non-clinical testing, the Evolution® Colonic **Stent** produced maximum temperature rises of 1.9 °C during 15 minutes of MR imaging (i.e., for one scanning sequence) performed in a 1.5 Tesla MR System (Magnetom, Siemens Medical Solutions, Malvern, PA, Software Numaris/4, Version Syngo MR 2002B DHHS Active-shielded) at an MR system reported whole-body-averaged SAR of 2.9 W/kg (associated with a calorimetry measured whole body averaged value of 2.1 W/kg).

3.0 Tesla Temperature Rise

In non-clinical testing, the Evolution® Colonic **Stent** produced maximum temperature rises of 2.7 °C during 15 minutes of MR imaging (i.e., for one scanning sequence) performed in a 3 Tesla MR System (Excite, HDx, General Electric Healthcare, Milwaukee, WI, Software 14X.M5) at an MR system reported whole-body-averaged SAR of 2.9 W/kg (associated with a calorimetry measured whole body averaged value of 2.7 W/kg).

Image Artifacts

MR image quality may be compromised if the area of interest is in within the lumen or within approximately 3 mm of the position of the Evolution® Colonic **Stent**, as found during non-clinical testing using the sequence: T1-weighted, spin echo and gradient echo pulse sequence in a 3.0 Tesla MR system (Excite, General Electric Healthcare, Milwaukee, WI, Software 14X.M5). Therefore, it may be necessary to optimize MR imaging parameters for the presence of this metallic **stent**.

Additional Information: The safety of performing an MRI procedure in a patient with overlapping colonic stents or other MRI-conditional device(s) in direct contact with this device has not been determined. Performing MRI in such situations is not recommended.

For US Patients Only

Cook recommends that the patient register the MR conditions disclosed in this IFU with the MedicAlert Foundation. The MedicAlert Foundation can be contacted in the following manners:

Mail: MedicAlert Foundation International
2323 Colorado Avenue
Turlock, CA 95382

Phone: 888-633-4298 (toll free)
209-668-3333 from outside the US

Fax: 209-669-2450

Web: www.medicalert.org

EQUIPMENT REQUIRED

Stent and delivery system

.035 inch wire guide

SELECTION OF STENT

A stent length of at least 4 cm longer than the stricture is recommended (e.g. approximately 2 cm longer at either end of the stricture)

PATIENT PREPARATION

1. Place wire guide through the stricture.

INSTRUCTIONS FOR USE

Illustrations

1. Remove protective tubing from stent.

2. Introduce delivery system in short increments over wire guide, into accessory channel, until it is endoscopically visualized exiting the scope.

3. To ensure that stent will bridge stricture after deployment, fluoroscopically position radiopaque markers on inner catheter beyond extremities of stricture to be crossed. **Note:** yellow marker on the delivery system located at proximal end of stent can also serve as an endoscopic/fluoroscopic reference for positioning of the proximal (nearest user) end of the stent relative to the stricture.

4. Confirm desired stent position fluoroscopically and deploy stent by removing red safety guard from the handle. (See Fig. 1)

5. Continue deploying stent by squeezing trigger. **Note:** Each trigger squeeze will deploy stent by an equal amount.

6. If stent repositioning is required during deployment, it is possible to recapture stent. **Note:** It is **not** possible to recapture stent after passing point-of-no-return, indicated when red marker on top of introducer has passed the point-of-no-return indicator on handle. (See Fig. 2) The stent can be recaptured a maximum of five times.

7. To recapture stent, push direction button on side of delivery system to opposite side. (See Fig. 3) **Note:** Hold thumb on button when squeezing trigger for first time to recapture.

8. Continue squeezing trigger as required to recapture stent by desired amount.

9. To resume deployment, push button to opposite side again and hold thumb on button for first stroke while squeezing trigger. (See Fig. 4)

10. When stent point-of-no-return has been passed, pull safety wire out of delivery handle near wire guide port. (See Fig. 5)

11. Continue deploying stent by squeezing trigger.

12. After deployment, fluoroscopically confirm full stent expansion. Once full expansion is confirmed, introduction system can be safely removed.

Upon completion of procedure, dispose of device per institutional guidelines for biohazardous medical waste.

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© COPYRIGHT COOK 2015 ©
COOK IRELAND LTD.