**COOK****MEDICAL****CE**  
0088

- EN**  
2-4  
Emergency Transtracheal Airway Catheter with Enk Oxygen Flow Modulator  
Instructions for Use
- DA**  
4-6  
Transtrakealt nødluftvejskateter med Enk modulator til iltflow  
Brugsanvisning
- DE**  
6-9  
Transtrachealer Katheter mit Enk Sauerstoff-Flowmodulator  
Gebrauchsanweisung
- EL**  
9-12  
Καθετήρας επείγουσας διατραχειακής αεροφόρου οδού με διαμορφωτή ροής οξυγόνου Enk  
Οδηγίες χρήσης
- ES**  
12-14  
Catéter transtraqueal de urgencia para la vía aérea con modulador de flujo de oxígeno Enk  
Instrucciones de uso
- FR**  
14-17  
Cathéter d'intubation transtrachéale d'urgence avec modulateur de débit d'oxygène d'Enk  
Mode d'emploi
- IT**  
17-20  
Catetere transtracheale per la ventilazione d'emergenza con modulatore del flusso di ossigeno Enk  
Istruzioni per l'uso
- NL**  
20-22  
Acute transtracheale luchtwegkatheter met Enk zuurstoffflowmodulator  
Gebruiksaanwijzing
- NO**  
22-25  
Transtrakealt akutt-luftveiskateter med Enk oksygenflowmodulator  
Bruksanvisning
- PT**  
25-27  
Cateter transtraqueal de emergência para vias aéreas com modulador do fluxo de oxigénio Enk  
Instruções de utilização
- SV**  
27-30  
Akut transtrakeal luftvägskateter med Enk syrgasflödesmodulator  
Bruksanvisning

C - T - D T J V E F M - R E V O

# EMERGENCY TRANSTRACHEAL AIRWAY CATHETER WITH ENK OXYGEN FLOW MODULATOR

**CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).**

## DEVICE DESCRIPTION

The Emergency Transtracheal Airway Catheter is a 6.0 French reinforced FEP catheter, introduced over a 15 gage needle.

The Enk Oxygen Flow Modulator has an oxygen supply connector at one end and a transtracheal catheter connector at the opposite end.

## INTENDED USE

The Emergency Transtracheal Airway Catheter is intended for emergency airway access when conventional endotracheal intubation cannot be performed.

The Enk Oxygen Flow Modulator allows manually controlled oxygen flow and is intended for emergency transtracheal catheter ventilation when conventional ventilation by mask or endotracheal tube cannot be performed.

## CONTRAINDICATIONS

None known

## WARNINGS

- **Firm control must be applied to the transtracheal catheter to maintain position fully within the trachea at all times during use. High oxygen pressure may cause an unsecured catheter to become dislodged from the airway, resulting in subcutaneous emphysema.**
- The Emergency Transtracheal Airway Catheter with Enk Oxygen Flow Modulator is not specifically designed for pediatric applications. Use with pediatric patients should be determined by attending physician.

## PRECAUTIONS

- Patients with airway injuries may have significant spinal injuries. Whenever possible, the cervical spine should be immobilized before beginning the procedure. Care should always be exercised to avoid additional spinal injuries.
- Whenever possible and appropriate, aseptic technique and local anesthetic should be used for the procedure.
- Do not reinsert the needle into the catheter once the needle is withdrawn. The needle may exit catheter wall.
- Physicians not experienced in the technique of transtracheal catheter ventilation should not use this device.

- The product is intended for use by physicians trained and experienced in emergency critical care procedures and transtracheal catheter ventilation techniques. Standard techniques for placement and monitoring of transtracheal catheters should be employed.
- The potential effects of phthalates on pregnant/nursing women or children have not been fully characterized and there may be concern for reproductive and developmental effects.

## **POTENTIAL ADVERSE EVENTS**

- Aspiration
- Barotrauma
- Bleeding
- Catheter dislodgement
- Hematoma
- Pneumomediastinum
- Subcutaneous emphysema
- Tracheoesophageal fistula

## **INSTRUCTIONS FOR USE**

### **Emergency Transtracheal Airway Catheter**

1. Prep the patient's anterior neck.
2. Identify the cricothyroid membrane between the cricoid and thyroid cartilages, and stabilize the trachea with the thumb and index finger.
3. Cannulate the trachea with the catheter tip facing caudad. A loss of resistance is felt when the trachea is entered.
4. Aspirate air into a water-filled 5 cc syringe to confirm catheter position within the tracheal lumen.
5. Remove the needle and advance the catheter caudad. Aspirate again to confirm position.
6. Connect the catheter to oxygen source, as directed in the following section.

### **Enk Oxygen Flow Modulator**

1. Connect the oxygen flow modulator to the tubing of an oxygen source.
2. Attach the Luer lock of the oxygen flow modulator to the transtracheal airway catheter.
 

**NOTE:** The pressure of the oxygen source must be high enough to create a pulsed gas flow through the flow modulator and the attached transtracheal catheter. Oxygen pressure of a common wall outlet (50 psi or 3.5 bar) is adequate.
3. Select an oxygen flow of at least 15, and up to 30, liters/minute.
4. Occlude and release all openings in the tubing of the oxygen flow modulator intermittently by using the thumb and index finger (e.g., 100 cycles per minute).

**CAUTION: Continuous closure of the openings may cause barotrauma.**

5. Observe the chest and the upper abdomen for slight outward and inward movements to confirm pulmonary inflation and exhalation.

**NOTE:** The pulsed oxygen flow should result in slight outward and inward movements of the chest and/or upper abdomen. In case of upper airway obstruction, gas discharge from the patient's lungs through the released openings of the oxygen flow modulator will require more time.

## HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

## REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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Bloomington, IN 47404 U.S.A.