
EN 9	EchoTip® Ultra Endobronchial HD Ultrasound Needle for Olympus Scopes Instructions for Use
BG 13	Ендобронхиална ултразвукова игла с висока разделителна способност EchoTip® Ultra за ендоскопи Olympus Инструкции за употреба
CS 18	Endobronchiální jehla EchoTip® Ultra pro ultrazvuk s vysokým rozlišením, pro endoskopy Olympus Návod k použití
DA 22	EchoTip® Ultra endobronkial HD ultralydsnål til Olympus skoper Brugsanvisning
NL 26	EchoTip® Ultra endobronchiale HD ultrasone naald voor Olympus endoscopen Gebruiksaanwijzing
ET 31	Endobronhiaalne kõrglahutusega ultrahelinõel EchoTip® Ultra endoskoopidele Olympus Kasutusjuhised
FR 35	Aiguille ultrasonore HD endobronchique EchoTip® Ultra pour les endoscopes Olympus Mode d'emploi
DE 40	EchoTip® Ultra endobronchiale HD-Ultraschallnadel für Olympus-Endoskope Gebrauchsanweisung
EL 45	Ενδοβρογχική βελόνα υπερήχων υψηλής ανάλυσης EchoTip® Ultra για ενδοσκόπια Olympus Οδηγίες χρήσης
HU 50	EchoTip® Ultra endobronchialis, nagy felbontású, echogén tű Olympus endoszkópokhoz Használati utasítás
IT 55	Ago ecografico endobronchiale ad alta definizione EchoTip® Ultra per endoscopi Olympus Istruzioni per l'uso
LV 60	EchoTip® Ultra endobronhiālā augstas precizitātes ultraskaņas adata Olympus endoskopiem Lietošanas instrukcija
LT 64	„EchoTip® Ultra“ endobronchinio didelės skiriamosios gebos ultragarso adata, skirta naudoti su „Olympus“ endoskopais Naudojimo nurodymai
NO 68	EchoTip® Ultra endobronkial HD-ultralydnål for Olympus-skop Bruksanvisning
PL 72	Igła EchoTip® Ultra do ultrasonografii wewnątrzozkrzelowej o wysokiej rozdzielczości do endoskopów Olympus Instrukcja użycia
PT 77	Agulha de ultrassons endobrônquica EchoTip® Ultra de alta definição para endoscópios Olympus Instruções de utilização
RO 82	Ac ecografic HD endobronşic EchoTip® Ultra pentru endoscoape Olympus Instrucţiuni de utilizare
SK 86	Endobronchiálna ultrazvuková ihla s vysokým rozlišením EchoTip® Ultra pre endoskopy Olympus Návod na použitie
ES 90	Aguja de ecografía endobronquial de alta definición EchoTip® Ultra para endoscopios Olympus Instrucciones de uso
SV 95	EchoTip® Ultra endobronkial HD ultraljudsnål för Olympus endoskop Bruksanvisning

Fig. 1

Biopsy Port Seal

Уплътнение на отвора за биопсия
Těsnění bioptického portu
Forsøgling på biopsiport
Afsluiting biopsiepoort
Biopsiapordi tihend
Capuchon étanche de l'orifice pour biopsie
Abdichtung des Biopsieanschlusses
Στεγανωτικό παρέμβυσμα θύρας βιοψίας
Biopszianyílás tömítése
Guarnizione del raccordo per biopsia
Biopsijas pieslēgvietais blīve
Biopsijos angos uždoris
Biopsiportforsøgling
Uszczelka portu biopsji
Vedante do orifício de biopsie
Sigiliul portului de biopsie
Uzáver bioptického portu
Junta del acceso de biopsia
Biopsiportstättning

Luer Connector

Конектор тип Luer
Konektor Luer
Luer-konnektor
Luer-connector
Luer-tüüpi ühendus
Connecteur Luer
Luer-Anschluss
Συνδετικό Luer
Luer-csatlakozó
Connettore Luer
Luer tipa savienotājs
Luerio jungtis
Luer-kobling
Łącznik Luer
Conector Luer
Conector Luer
Luerov konektor
Conector Luer
Luer-anslutning

Locking Clip

Блокиращ клипс
Aretovací svorka
Låseklemme
Borgclip
Kinnitusklamber
Pince de verrouillage
Arretierclip
Κλιπ ασφάλισης
Rögzítőkarocs
Clip di bloccaggio
Noslégskava
Fiksuojamasis spaustukas
Låseklips
Zacisk blokujący
Grampo de fixação
Clemă blocantă
Poistná objímka
Abrazadera de fijación
Låsklämma

Adapter Top

Горна част на адаптера
Horní strana adaptéru
Øverste del af adapter
Adaptertop
Adapteri ülemine osa
Haut de l'adaptateur
Oberteil des Adapters
Κορυφή του προσαρμογέα
Adapter teteje
Sommità dell'adattatore
Adaptera augšdaļa
Adapterio viršus
Adaptertopp
Część górna łącznika
Parte de cima do adaptador
Capacul adaptorului
Vrch adaptéra
Parte superior del adaptador
Adapters övre del

Adapter Base

Долна част (основа) на адаптера
Základna adaptéru
Nederste del af adapter
Adapterbasis
Adapteri alumine osa
Base de l'adaptateur
Unterteil des Adapters
Βάση του προσαρμογέα
Adapter alapja
Base dell'adattatore
Adaptera pamatne
Adapterio pagrindas
Adapterbase
Podstawa łącznika
Base do adaptador
Baza adaptorului
Podstavec adaptéra
Base del adaptador
Adapters bas

Fig. 2


Fig. 3



Fig. 4

Plunger

Бутало
Píst
Stempel
Plunjer
Kolb
Piston
Kolben
Εμβολο
Dugattyú
Stantuffo
Virzulis
Stúmoklis
Stempel
Tłoczek
Êmbolo
Piston
Piest
Êmbolo
Kolv

Stopcock "open"

Кранче "отворено"
„Otevřený“ uzavírací kohout
Stophane "åben"
Afsluitkraan 'open'
Sulgemisventiil avatud
Robinet ouvert
Absperrhahn geöffnert
Στρόφιγγα ανοικτή
Elzárócsap nyitott állásban
Rubinetto aperto
Slægkrans atvërts
Kranelis „atviras“
Stoppekran åpen
Kranik w pozycji otwartej
Torneira de passagem "aberta"
Robinet „deschis“
Otvorený kohútik
Llave de paso «abierta»
Kran "öppen"

Fig. 5



Fig. 6

Sliding Sheath Adjuster

Плъзгащ се регулатор на защитната обвивка
Posuvné nastavovací zařizení sheathu
Glidende sheath-justering
Schuifregeling van de sheath
Liughülsi regulaator
Ajusteur de gaine coulissant
Verschiebbarer Schleuseneinsteller
Συρόμενος προσαρμογέας θηκκαριού
Elcsúsztatható hüvelyhossz-beállító
Regolatore scorrevole della guaina
Pärbidämais ievadapvalka regulētājs
Slankusis vamzdelio reguliatorius
Glidende hylsejusteringsmekanisme
Przesuwany regulator koszulki
Regulador deslizante da bainha
Ajustor glisant al tecii
Posuvný nastavovač puzdra
Ajustador deslizante de la vaina
Glidande hylsjusterare

Sheath Reference Mark

Означение за дължината на защитната обвивка
Referenční značka sheathu
Referencemærke på sheath
Referentiemarkering van de sheath
Hülsi tähis
Repère de référence sur la gaine
Schleusenbezugsmarke
Σήμανση αναφοράς θηκκαριού
Hüvelyhossz referencijele
Indicatore di riferimento della guaina
Ievadapvalka atskaites žymē
Vamzdelio atskaitos žymė
Referansemerke på hylsen
Znacznik referencyjny wysunięcia koszulki
Marca de referenția da bainha
Marcaj de referință al tecii
Referenčná značka puzdra
Marca de referenția de la vaina
Hylsans referensmarkering

Needle Handle

Дръжка на иглата
Rukojeť jehly
Nålegreb
Naaldhandvat
Nöela käepide
Poignée de l'aiguille
Nadelgriff
Λαβή βελόνας
Tű fogantyúja
Impugnatura dell'ago
Ađatas rokturis
Adatos kotelis
Nålehåndtak
Uchwył igły
Punho da agulha
Månerul acului
Rukoväť ihly
Mango de la aguja
Nålhandtag



Safety Ring

Предпазен пръстен
Pojistný prstenec
Sikkerhedsring
Veiligheidsring
Kaitserõngas
Bague de sécurité
Sicherheitsring
Δακτύλιος ασφαλείας
Biztonsági gyűrű
Anello di sicurezza
Drošības gredzens
Apsauginis žiedas
Sikkerhetsring
Pierścień bezpieczeństwa
Anel de siguranța
Inel de siguranță
Bezpečnostný krúžok
Anillo de seguridad
Säkerhetslåsring

Stylet Hub

Втулка на стилето
Nástavec stiletu
Stiletmuffe
Aanzetstuk van het stilet
Stiletkinnitus
Embase du stylet
Mandrinansatz
Ομφαλός του στείλεου
Mandrin végidoma
Pomello del mandrino
Stilettes galvina
Zondo kaištis
Stilletfeste
Głowica mandrynu
Conector do estilete
Acoperitoarea stiletului
Hrdlo stiletu
Conector del estilete
Mandrångfatning

Fig. 7

ECHOTIP® ULTRA ENDOBRONCHIAL HD ULTRASOUND NEEDLE FOR OLYMPUS SCOPES

CAUTION: U.S. Federal Law restricts this device to sale by or on the order of a physician.

INTENDED USE

This device is used to sample targeted submucosal and extramural lesions within or adjacent to the tracheobronchial tree or gastrointestinal tract through the accessory channel of an ultrasound endoscope for Fine Needle Aspiration (FNA).

DEVICE DESCRIPTION

The device consists of a handle with adjustable components to allow the user to adjust the extension of the needle and sheath. The device has a 4.1 Fr outer sheath for protection of the needle and the scope, and it is attached to the accessory channel of an Olympus endobronchial ultrasound (EBUS) scope using an adapter. The device has a maximum needle extension length of 5 cm and sheath adjustment length of 3 cm. The device is supplied with a 10 ml syringe to utilize the standard vacuum technique for sample collection, if desired, and an accessory channel adapter. The device is available in two different needle sizes: 22 and 25 ga.

CONTRAINDICATIONS

Contraindications include those specific to the primary endoscopic ultrasound procedure to be performed in gaining access to the desired position to visualize the targeted site.

Relative contraindications include, but are not limited to: coagulopathy.

POTENTIAL ADVERSE EVENTS

Potential adverse events associated with bronchial endoscopy and gastrointestinal endoscopy include, but are not limited to: allergic reaction to medication, allergic reaction to nickel, aspiration, cardiac arrhythmia or arrest, damage to blood vessels, death, discomfort, fever, hemorrhage, hypotension, infection, inflammation, nerve damage, pain, perforation, pneumoperitoneum, respiratory depression or arrest, sepsis, septicemia/bacteremia and tumor seeding.

WARNINGS

Not for use in the heart or vascular system.

The tip of the needle and stylet are sharp and could cause injury to the patient or user if not used with caution.

This device is designed for single use only. Attempts to reprocess, re sterilize, and/or reuse may lead to device failure and/or transmission of disease.

If the package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook Medical for return authorization.

PRECAUTIONS

Refer to package label for minimum channel size required for this device.

The needle must be retracted into the sheath and the thumbscrew on the safety ring must be locked at the 0 cm mark to hold the needle in place prior to introduction, advancement or withdrawal of the device. Failure to retract the needle may result in damage to the scope.

Ensure the stylet is fully inserted when advancing the needle into the target site.

Do not use this device for any purpose other than stated intended use.

Use of this device is restricted to a trained healthcare professional (Rx only).

When targeting multiple sites, adequate care should be taken to minimize the risk of tumor seeding.

This device is intended for use with an Olympus EBUS scope.

PLACEMENT OF ACCESSORY CHANNEL ADAPTER

Illustrations

1. Examine the adapter and note features (**Fig. 1**).
2. Place the adapter, in the open position, in line with the scope accessory channel metal hub.
3. Slide the adapter on to the accessory channel ensuring that the base of the adapter slides under the accessory channel metal hub (**Fig. 2**).
4. Close the adapter ensuring it is locked (**Fig. 3**).
5. The adapter will now function as a Luer lock to which the needle can be connected. **Note:** If irrigation of accessory channel is desired/necessary, irrigate through closed adapter prior to insertion or after removal of needle.

SYSTEM PREPARATION

1. Examine syringe. It has two plunger locks that must be depressed to advance plunger. The tip of the syringe has a Luer lock fitting with a stopcock on the side port. Air can be exchanged when the stopcock is in "open" position, aligned with syringe (**Fig. 4**).
2. Prepare syringe as follows:
 - a) With stopcock in "open" position, depress plunger locks and fully advance plunger into syringe.
 - b) Turn stopcock to "closed" position.
 - c) Pull plunger back until it is locked into place at desired setting, creating a vacuum.
 - d) Set prepared syringe aside until aspiration is required.

INSTRUCTIONS FOR USE

1. Identify the desired target site by endobronchial ultrasound.
2. With the needle retracted into the sheath and the thumbscrew on the safety ring locked at the 0 cm mark to hold the needle in place, introduce the ultrasound needle into the accessory channel of the scope (**Fig. 5**).
Caution: If resistance is encountered on needle introduction, reduce angulation of the scope until smooth passage is allowed.
3. Advance the device, in small increments, until the Luer lock fitting at the base of the sliding sheath adjuster meets the Luer fitting of the accessory channel port.
4. Attach the device to the accessory channel port by rotating the device handle clockwise until the fittings are connected.
5. Adjust the sheath to the desired position, ensuring that it is visible on the endoscopic image, emerging from the working channel of the scope. To adjust length, loosen the thumbscrew lock on the sliding sheath adjuster and slide until the preferred length is attained. **Note:** The reference mark for the sheath length will appear in the sliding sheath adjuster window (**Fig. 6**). Tighten the thumbscrew on the sliding sheath adjuster to maintain the preferred sheath length.
6. While maintaining position of ultrasound scope, set the needle extension to the desired length by loosening the thumbscrew on the safety ring, and advancing it until the desired reference mark for needle advancement appears in the window of the safety ring (**Fig. 7**). Tighten the thumbscrew to lock the safety ring in place. **Note:** The number in the safety lock ring window indicates the extension of the needle in centimeters. **Caution:** During needle adjustment or extension, ensure the device has been attached to the accessory channel of the scope. Failure to attach the device prior to needle adjustment or extension may result in damage to the scope.
7. Extend the needle, by advancing the needle handle of the device to the pre-positioned safety ring, into the target site. **Caution:** If excessive resistance is encountered on needle advancement, retract the needle into the sheath with the thumbscrew locked at the 0 cm mark, reposition the scope and attempt needle advancement from another angle. Failure to do so may result in needle breakage, device damage or malfunction.
8. Standard vacuum syringe techniques may be applied for sample collection (Steps 9-11) or, if desired, other techniques that may or may not incorporate the use of this stylet may be used.
9. Remove the stylet from the ultrasound needle by gently pulling back on the plastic hub seated in the Luer fitting of the needle handle. Preserve the stylet for use if additional sample collection is required later.
10. Attach the Luer lock fitting of the previously prepared syringe securely to the Luer fitting on the needle handle.

11. Turn the stopcock to the “open” position aligned with the syringe, allowing the negative pressure in the syringe to aspirate the sample. Move the needle handle back and forth within the target site. **Note:** Do not remove the needle from the target site during sample collection.

12. Upon completion of sample collection, if using the syringe, turn the stopcock to the “closed” position. Retract the needle completely into the sheath by pulling back on the needle handle and lock the thumbscrew on the safety ring at the 0 cm mark to hold the needle in place.

13. Disconnect the Luer lock fitting of the device from the accessory channel port, by rotating the device handle counterclockwise, and withdraw the entire device from the scope.

14. Unlock the thumbscrew on the safety ring and extend the needle.

15. Use a new syringe or the stylet to expel the specimen, then prepare per institutional guidelines.

16. Remaining material may be recovered for examination by reinserting the stylet or flushing the device.

17. For additional sample collection from the same target site, gently reinsert the stylet into the Luer fitting on the needle handle. **Note:** Prior to reinserting the stylet, wipe with saline or sterile water. While supporting the stylet at the Luer fitting, advance the stylet in small increments, until the hub is engaged in the fitting.

18. Additional samples may be obtained by repeating Step 2 through 16 of the “Instructions for Use”.

19. Upon completion of the procedure, disconnect the Luer lock fitting of the device from the accessory channel port, by rotating the device handle counterclockwise, and withdraw the entire device from the scope.

20. To remove the adapter from the scope, pull back on the locking clip, lift the top and slide off the scope.

Upon completion of the procedure, dispose of the device, including the adapter, per institutional guidelines for biohazardous medical waste.

HOW SUPPLIED

Store in a dry location, away from temperature extremes.

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© COPYRIGHT COOK 2019 ©
COOK IRELAND LTD.