EN 5 **EchoTip® Ultra Celiac Plexus Neurolysis Needle**

Instructions for Use

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
Fig. 3

INTENDED USE

This device is used to deliver neurolysing agents to the celiac plexus under guidance by endoscopic ultrasound.

NOTES

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

Do not use this device for any purpose other than the stated intended use.

If the package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Store in a dry location, away from temperature extremes.

(Use of this device is restricted to or on the order of a licensed healthcare practitioner.) U.S. Federal Law restricts this device to sale by or on the order of a physician.

CONTRAINDICATIONS

Contraindications include those specific to the primary endoscopic procedure to be performed in gaining access to the desired injection site. Contraindications to Celiac Plexus Neurolysis include, but are not limited to: coagulopathy, use of anticoagulants, aspirin or NSAIDS.

POTENTIAL COMPLICATIONS

Potential complications associated with gastrointestinal endoscopy include, but are not limited to: perforation, hemorrhage, aspiration, fever, infection, allergic reaction to medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest. Potential complications associated with celiac plexus neurolysis include but are not limited to: transient diarrhea, hypotension, transient increasing epigastric pain, lower extremity weakness and paresthesias, paraplegia, renal puncture and prolonged gastroparesis.

PRECAUTIONS

Refer to package label for minimum channel size required for this device. A thorough understanding of the technical principles, clinical applications and risks associated with endoscopic ultrasound and celiac plexus neurolysis is necessary before using this device. Use of prophylactic antibiotics is recommended for endoscopic ultrasound-guided celiac plexus neurolysis. Use of total anesthesia may be needed in some cases. The Celiac Plexus Neurolysis Needle should only be used by, or under the supervision of, physicians thoroughly trained in endoscopic ultrasound and celiac plexus neurolysis. The needle must be retracted into the sheath and the thumbscrew on the safety ring locked to hold the needle in place prior to introduction, advancement or withdrawal of the device. Failure to retract the needle may result in damage to the endoscope.

SYSTEM PREPARATION

Advance the device into the endoscope to determine preferred sheath length. To adjust length, loosen the thumbscrew lock on the sliding sheath adjuster and slide until the preferred length is obtained. **Note: The reference mark for the sheath length will appear**

in the sliding sheath adjuster window. (See fig. 1) Tighten thumbscrew on sliding sheath adjuster to maintain the preferred sheath length.

INSTRUCTIONS FOR USE

Illustrations

1. Introduce the ultrasound endoscope and advance to the appropriate position for neurolysis injection into the celiac plexus.
 2. Identify the desired injection site based on endoscopic ultrasound.
 3. With the needle retracted into the sheath and the thumbscrew on the safety ring locked to hold the needle in place, introduce the ultrasound needle into the accessory channel of the endoscope. (See fig. 2)
 4. Advance the device, in small increments, until the Luer lock fitting at the base of the sliding sheath adjuster meets the fitting on the endoscope accessory channel.
 5. Attach the device to the endoscope accessory channel port by rotating the device handle until the fittings are connected.
 6. With the endoscope and device straight, advance the needle to the desired length by loosening the thumbscrew on the safety ring, and advancing it until the desired reference mark for needle advancement appears in the window of the safety ring. (See fig. 3) Tighten the thumbscrew to lock the safety ring in place. **Note: The number in the safety lock ring window indicates the extension of the needle in centimetres. Caution: During needle adjustment or extension, ensure the device has been attached to the accessory channel of the endoscope. Failure to attach the device prior to needle adjustment or extension may result in damage to the endoscope.**
 7. While maintaining the position of the endoscope, extend the needle by advancing the needle handle of the device to the pre-positioned safety ring.
 8. Advance the needle into the desired injection site. **Note: The EchoTip needle allows visualization during ultrasound.**
 9. Remove the stylet from the ultrasound needle by gently pulling back on the plastic hub seated in the metal fitting of the needle handle. Preserve the stylet for use if additional neurolysing injections are required later.
 10. Inadvertent puncture of blood vessels may be determined by attempting to aspirate blood. A luer lock syringe may be attached to the luer lock fitting on the handle of the device for this purpose.
 11. Injection of a local anaesthetic and neurolysing agents should be performed with a luer lock syringe.
 12. For additional neurolysis injections in alternative sites on the celiac plexus, gently reinsert the stylet wire into the metal fitting on the handle. **Note: Prior to reinserting stylet, wipe with saline or sterile water. While supporting the sheath at the Luer lock fitting, advance the stylet in small increments, until the stylet hub is engaged in the fitting.**
 13. Upon completion of the procedure, retract the needle completely into the sheath by pulling back on the needle handle and lock the thumbscrew on the safety ring to hold the needle in place.
 14. Disconnect the Luer lock fitting of the device from the accessory channel and withdraw the entire device from the endoscope.
 15. Additional injections may be performed by repeating Steps 2-14.
- Upon completion of the procedure, dispose of the device per institutional guidelines for biohazardous medical waste.**

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© COPYRIGHT COOK 2017 ©
COOK IRELAND LTD.