
EN 3	Pushing Catheter Instructions for Use
CS 5	Tlačný katetr Návod k použití
DA 7	Skubbekateter Brugsanvisning
NL 9	Pushing-katheter Gebruiksaanwijzing
FR 11	Cathéter pousse-prothèse Mode d'emploi
DE 13	Platzierungskatheter Gebrauchsanweisung
EL 15	Καθετήρας ώθησης Οδηγίες χρήσης
HU 18	Tolókatéter Használati utasítás
IT 20	Catetere di spinta Istruzioni per l'uso
NO 23	Skyvekateter Bruksanvisning
PL 25	Cewnik popychający Instrukcja użycia
PT 28	Cateter posicionador Instruções de utilização
ES 31	Catéter empujador Instrucciones de uso
SV 34	Pushing-kateter Bruksanvisning

Removable Stylet

Snímatelný stylet
 Aftagelig stilet
 Verwijderbaar stilet
 Stylet amovible
 Entfernbarer Mandrin
 Αφαιρούμενος σπειλεός
 Eltávolítható mandrin
 Mandrino rimovibile
 Stilet som kan fjernes
 Usuwalny mandryn
 Estilete amovível
 Estilete extraíble
 Borttagbar mandräng

Proximal Wire Port (PWP)

Proximální vodičí port (PWP)
 Proksimal kateterlederport (PWP)
 Proximale voerdraadpoort (PWP)
 Orifice pour guide proximal (PWP)
 Proximaler Drahtport (PWP)
 Εγγύς θύρα σύρματος (PWP)
 Proximális drótnyílás (PWP)
 Ingresso prossimale della guida (PWP)
 Proksimal vaierport (PWP)
 Proksymalny port dla przewodnika (PWP)
 Orificio proximal do fio (PWP)
 Acceso proximal de la guía (PWP)
 Proximal ledarport (PWP)

IDE Port

Port IDE
 IDE-port
 IDE-poort
 Orifice IDE
 IDE-Port
 Θύρα IDE
 IDE nyílás
 Accesso IDE
 IDE-port
 Port IDE
 Orificio de IDE
 Acceso de IDE
 IDE-port

INTENDED USE

Used for endoscopic biliary stent placement to drain obstructed bile ducts.

NOTES

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

If the package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Notify Cook for return authorization.

Do not use this device for any purpose other than the stated intended use.

Store in a dry location, away from temperature extremes.

Use of this device restricted to a trained healthcare professional.

CONTRAINDICATIONS

Contraindications include those specific to ERCP.

POTENTIAL COMPLICATIONS

Potential complications associated with ERCP include, but are not limited to: pancreatitis, cholangitis, aspiration, perforation, hemorrhage, infection, sepsis, allergic reaction to contrast or medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

Additional complications that may occur with biliary stent placement include, but are not limited to: trauma to the biliary tract or duodenum, obstruction of the pancreatic duct, stent migration.

PRECAUTIONS

Refer to package label for minimum channel size required for this device.

The wire guide diameter and the inner lumen of the wire-guided device must be compatible.

Only Cook biliary and pancreatic stents are recommended for use with this device.

Do not use this device if stent cannot advance through strictured area.

INSTRUCTIONS FOR USE (SEE FIG. 1)

Illustrations

I. If using Intra Ductal Exchange (IDE) port and a pre-positioned short wire:

1. Unlock short wire from Wire Guide Locking Device and load stent and positioning sleeve over pre-positioned wire guide. **Note:** Use positioning sleeve included with this device.

2. Advance pushing catheter over pre-positioned wire guide, ensuring that wire guide exits pushing catheter at IDE port.
3. As the front flap of stent enters accessory channel, keep positioning sleeve over the back flap until completely in accessory channel. Slide positioning sleeve over pushing catheter and wire guide, keeping it clear of accessory channel.
4. Relock wire guide.
5. With elevator open, advance device in short increments until visualized exiting endoscope.
6. Advance device into appropriate duct.
7. Fluoroscopically and endoscopically confirm desired stent placement.
8. Retract pushing catheter until IDE port is visible at Wire Guide Locking Device and resistance is felt.
9. Unlock wire guide and perform standard exchange over remaining 6 cm of pushing catheter. Once distal tip exits Wire Guide Locking Device, relock wire guide.

II. If using Proximal Wire Port (PWP) and a pre-positioned long wire guide:

Note: For best results, wire guide should be kept wet.

1. Remove stylet.
2. Load stent and positioning sleeve onto pre-positioned wire guide.
3. Advance pushing catheter over pre-positioned wire guide, ensuring that wire guide exits pushing catheter at PWP.
4. As front flap of stent enters accessory channel, keep positioning sleeve over back flap until completely in accessory channel. Slide positioning sleeve over pushing catheter, keeping it clear of accessory channel.

REFER TO STEPS 5-7 IN "SECTION I", THEN RESUME WITH STEP 5 BELOW:

5. Remove pushing catheter using standard long wire exchange technique.

Upon completion of procedure, dispose of device per institutional guidelines for biohazardous medical waste.

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© COPYRIGHT COOK 2019 ©
COOK IRELAND LTD.