EN 6	EchoTip Ultra with High Definition Ultrasound Access Needle Instructions for Use
BG 9	Игла EchoTip Ultra за ултразвуков достъп с висока разделителна способност Инструкции за употреба
CS 13	Přístupová jehla EchoTip Ultra s vysokým ultrazvukovým rozlišením Návod k použití
DA 17	EchoTip Ultra med ultralyds-indføringsnål med høj opløsning Brugsanvisning
NL 21	EchoTip Ultra met high-definition ultrasone introductienaald Gebruiksaanwijzing
ET 25	Kõrglahutusega ultraheli-juurdepääsunõelaga EchoTip Ultra Kasutusjuhised
FR 29	Aiguille de ponction pour écho-endoscopie avec haute définition EchoTip Ultra Mode d'emploi
DE 33	EchoTip Ultra mit High-Definition-Ultraschall-Zugangsnadel Gebrauchsanweisung
EL 37	EchoTip Ultra με βελόνα εισαγωγέα υπερήχων υψηλής ανάλυσης Οδηγίες χρήσης
HU 41	EchoTip Ultra nagy felbontású ultrahangos hozzáférési tűvel Használati utasítás
IT 45	EchoTip Ultra con ago per accesso ecoguidato ad alta definizione Istruzioni per l'uso
LV 49	EchoTip Ultra ar augstas precizitātes ultraskaņas piekļuves adatu Lietošanas instrukcija
LT 53	„EchoTip Ultra“ su didelės skiriamosios gebos ultragarso prieigos adata Naudojimo nurodymai
NO 57	EchoTip Ultra med ultralydtilgangsnål med høy oppløsning Bruksanvisning
PL 61	EchoTip Ultra z igłą dostępową do endoskopów ultrasonograficznych wysokiej rozdzielczości Instrukcja użycia
PT 65	Agulha introdutora de ultra-sons EchoTip Ultra de alta definição Instruções de utilização
RO 69	Ac EchoTip Ultra pentru acces cu ultrasunete de înaltă definiție Instrucțiuni de utilizare
SK 73	Přístupová ihla EchoTip Ultra s vysokým ultrazvukovým rozlišením Návod na použitie

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


3.

4.

ECHOTIP ULTRA WITH HIGH DEFINITION ULTRASOUND ACCESS NEEDLE

CAUTION: U.S. Federal Law restricts this device to sale by or on the order of a physician.

INTENDED USE

This device is used to access or sample submucosal and extramural lesions of the gastrointestinal tract, access of the following: the intra- or extra-hepatic bile ducts, pancreatic ducts, cystic duct, gallbladder or for delivery of injectable materials into tissues through the accessory channel of an ultrasound endoscope.

NOTES

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

Do not use this device for any purpose other than stated intended use.

This needle is intended to be used with legally marketed drugs and devices.

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify COOK Medical for return authorization.

This device can be used together with all ECHO compatible products.

Aspiration syringe not intended to be used for injection.

Store in a dry location, away from temperature extremes.

Use of this device restricted to a trained healthcare professional.

CONTRAINDICATIONS

Those specific to primary endoscopic procedure to be performed in gaining access to desired site.

Coagulopathy

POTENTIAL COMPLICATIONS

Those associated with gastrointestinal endoscopy include, but are not limited to: perforation, hemorrhage, aspiration, fever, infection, allergic reaction to medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest, damage to blood vessels, nerve damage and acute pancreatitis.

WARNING

Not for use in the heart or vascular system.

PRECAUTIONS

When performing fine needle injection (FNI), ensure that injectable material

does not contact untargeted site.

When targeting multiple sites, replace device for each site.

Refer to package label for minimum channel size required for this device.

Needle must be retracted into sheath and thumbscrew on safety ring must be locked to hold needle in place prior to introduction, advancement or withdrawal of device. Failure to retract needle may result in damage to endoscope.

SYSTEM PREPARATIONS

1. For FNI and/or Access, skip to step 4.
2. For fine needle aspiration (FNA), examine syringe. It has two plunger locks that must be depressed to advance plunger. Tip of syringe has a Luer lock fitting with a stopcock on the side port. Air can be exchanged when stopcock is in "open" position, aligned with syringe. (See fig. 1)
3. Prepare syringe as follows:
 - a. With stopcock in "open" position, depress plunger locks and fully advance plunger into syringe.
 - b. Turn stopcock 90° to "closed" position.
 - c. Pull plunger back until it is locked into place at desired setting, creating suction.
 - d. Set prepared syringe aside if aspiration biopsy is desired.
4. Advance device into ultrasound endoscope to determine preferred sheath length. To adjust length, loosen thumbscrew lock on sliding sheath adjuster and slide until preferred length is attained.
Note: Reference mark for sheath length will appear in sliding sheath adjuster window. (See fig. 2) Tighten thumbscrew on sliding sheath adjuster to maintain preferred sheath length.

INSTRUCTIONS FOR USE

1. Identify desired biopsy/injection site by endoscopic ultrasound.
2. Introduce ultrasound needle into accessory channel and advance device in small increments until Luer lock fitting at base of sliding sheath adjuster meets fitting on accessory channel. (See fig. 3)
3. Attach device to endoscope accessory channel port by rotating device handle until fittings are connected.
4. With ultrasound endoscope and device straight, adjust needle to desired length by loosening thumbscrew on safety ring, and advancing it until desired reference mark for needle advancement appears in window of safety ring. (See fig. 4) Tighten thumbscrew to lock safety ring in place. **Note:** Number in safety lock ring window indicates extension of needle in centimeters. **Caution:** During needle adjustment or extension, ensure device is attached to accessory channel. Failure to attach device prior to needle adjustment or extension may result in damage to endoscope.
5. While maintaining position of ultrasound endoscope, extend needle by advancing needle handle to pre-positioned safety ring.
6. Advance needle into lesion, intra- or extra-hepatic bile ducts, pancreatic duct, cystic duct or gallbladder.
7. Remove stylet from needle by rotating stylet cap anti clockwise and gently

pulling it back from metal fitting of handle. Preserve stylet for use if additional punctures are desired later.

8. For maintaining access a 0.035" wire guide (Cook Medical recommended) can be advanced through needle. Skip to step 12. **Note:** a 0.035 wire guide can also be advanced post aspiration. For FNA, attach Luer lock fitting of previously prepared syringe securely onto fitting on handle. For FNI, follow manufacturer's guidelines for use of injectable materials and skip to step 12.

9. Turn stopcock to "open" position aligned with syringe, allowing negative pressure in syringe to aspirate cells. Gently move handle in small increments, back and forth within biopsy site. **Note:** Do not remove needle from biopsy site during cell collection

10. Upon completion of FNA, turn stopcock to "close" position.

11. Retract needle into sheath by pulling back on handle shaft. Lock thumbscrew on safety ring at 0cm mark.

12. Disconnect Luer lock fitting from accessory channel and withdraw entire device from ultrasound endoscope. **Note:** If a wireguide has been advanced, withdraw needle over wireguide.

13. The remaining steps are for FNA only.

14. Unlock thumbscrew on safety ring and extend needle. Exchange locking syringe for an empty syringe; attach to Luer lock fitting on handle.

15. Push in plunger to expel aspirate specimen, then prepare per institutional guidelines.

16. Remaining aspirate may be recovered for examination by reinserting stylet or flushing device.

17. For additional aspiration from same targeted site, gently reinsert stylet into metal fitting on the needle handle. **Note:** Prior to reinserting the stylet, wipe with saline or sterile water. While supporting the stylet at the Luer lock fitting, advance the stylet in small increments, until the hub is engaged in the fitting. Rotate stylet cap clockwise until the fittings are connected.

18. Additional samples may be obtained by repeating Step 2 through 16 of the "Instructions for Use".

Upon completion of procedure, dispose of device per institutional guidelines for biohazardous medical waste.

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© COPYRIGHT COOK 2016 ©
COOK IRELAND LTD.