- EN 4** Melker Cuffed and Uncuffed Emergency Cricothyrotomy Catheter Sets
Instructions for Use
- CS 5** Soupravy katetrů Melker s manžetou a bez manžety pro naléhavou koniotomii
Návod k použití
- DA 6** Melker nødcricothyreotomi katetersæt med og uden cuff
Brugsanvisning
- DE 7** Melker-Kanülen-Sets für Notfall-Krikothyreotomie, mit und ohne Manschette
Gebrauchsanweisung
- EL 9** ΣΕΤ καθετήρα επείγουσας κρικοθυρεοειδοτομής Melker με και χωρίς δακτύλιο (cuff)
Οδηγίες χρήσης
- ES 10** Equipos de catéter de cricotirotomía de urgencia Melker con y sin balón
Instrucciones de uso
- FR 12** Sets de cathéter de Melker pour cricothyrotomie en urgence avec et sans ballonnet
Mode d'emploi
- IT 13** Set con catetere cuffiato e non cuffiato per cricotirotomia d'emergenza Melker
Istruzioni per l'uso
- NL 15** Melker acute cricothyreotomiekathetersets met en zonder cuff
Gebruiksaanwijzing
- NO 16** Melker katetersett til akutt krikotyreotomi med og uten mansjett
Bruksanvisning
- PT 17** Conjuntos de cateteres para cricotirotomia de emergência com e sem cuff Melker
Instruções de utilização
- SV 19** Melker kateterset för akut krikotyreotomi med och utan cuff
Bruksanvisning


Fig. 1

- | | |
|-------------------------|--------------------------|
| A. Thyroid Cartilage | A. Cartilage thyroïde |
| B. Access Site | B. Site d'accès |
| C. Cricoid Cartilage | C. Cartilage cricoïde |
| A. Štítná chrupavka | A. Cartilagine tiroidea |
| B. Přístupové místo | B. Sito di accesso |
| C. Prstencová chrupavka | C. Cartilagine cricoidea |
| A. Thyreoideabrusk | A. Schilddrüsenknorpel |
| B. Indføringssted | B. Introductieplaats |
| C. Cricoideabrusk | C. Cricoïdknorpel |
| A. Cartilago thyreoidea | A. Skjoldbrusk |
| B. Zugangsstelle | B. Tilgangssted |
| C. Cartilago cricoidea | C. Ringbrusk |
| A. Θυρεοειδής χόνδρος | A. Cartilagem tiróide |
| B. Σημείο προσπέλασης | B. Local de acceso |
| C. Κρικοειδής χόνδρος | C. Cartilagem cricóide |
| A. Cartílago tiroides | A. Sköldbrusk |
| B. Lugar de acceso | B. Åtkomstområde |
| C. Cartílago cricoides | C. Ringbrusk |


Fig. 2




Fig. 3


Uncuffed Airway Catheter
 Katetr pro dýchací cesty bez manžety
 Luftvejskateter uden cuff
 Trachealkanüle ohne Manschette
 Καθετήρας αεραγωγών χωρίς δακτύλιο
 Catéter para la vía aérea sin balón
 Cathéter d'intubation sans ballonnet
 Catetere di ventilazione non cuffiato
 Luchtwegkatheter zonder cuff
 Luftveiskateter uten mansjett
 Cateter para vias aéreas sem cuff
 Luftvägskateter utan kuff


Cuffed Airway Catheter
 Katetr pro dýchací cesty s manžetou
 Luftvejskateter med cuff
 Trachealkanüle mit Manschette
 Καθετήρας αεραγωγών με δακτύλιο
 Catéter para la vía aérea con balón
 Cathéter d'intubation à ballonnet
 Catetere di ventilazione cuffiato
 Luchtwegkatheter met cuff
 Luftveiskateter med mansjett
 Cateter para vias aéreas com cuff
 Luftvägskateter med kuff

Fig. 4

Uncuffed Airway Catheter
 Katetr pro dýchací cesty bez manžety
 Luftvejskateter uden cuff
 Trachealkanüle ohne Manschette
 Καθετήρας αεραγωγών χωρίς δακτύλιο
 Catéter para la vía aérea sin balón
 Cathéter d'intubation sans ballonnet
 Catetere di ventilazione non cuffiato
 Luchtwegkatheter zonder cuff
 Luftveiskateter uten mansjett
 Cateter para vias aéreas sem cuff
 Luftvägskateter utan kuff

Cuffed Airway Catheter
 Katetr pro dýchací cesty s manžetou
 Luftvejskateter med cuff
 Trachealkanüle mit Manschette
 Καθετήρας αεραγωγών με δακτύλιο
 Catéter para la vía aérea con balón
 Cathéter d'intubation à ballonnet
 Catetere di ventilazione cuffiato
 Luchtwegkatheter met cuff
 Luftveiskateter med mansjett
 Cateter para vias aéreas com cuff
 Luftvägskateter med kuff

Fig. 5

MELKER CUFFED AND UNCUFFED EMERGENCY CRICOTHYROTOMY CATHETER SETS

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

Melker Cuffed and Uncuffed Emergency Cricothyrotomy Catheter Sets consist of an airway catheter, scalpel, syringe, introducer needle, catheter introducer needle, wire guide and curved dilator.

NOTE: Emergency airway catheter with 6 mm inner diameter allows use of standard positive pressure ventilation techniques and will also permit spontaneous patient breathing in adults.

INTENDED USE

Melker Cuffed and Uncuffed Emergency Cricothyrotomy Catheter Sets are intended to establish emergency airway access when endotracheal intubation cannot be performed. Airway access is achieved using percutaneous entry (Seldinger) technique via the cricothyroid membrane. Subsequent dilation of the tract and tracheal entrance site permits passage of the emergency airway catheter.

CONTRAINDICATIONS

No absolute contraindications known.

WARNINGS

- The Melker Cuffed and Uncuffed Emergency Cricothyrotomy Catheter Sets are not specifically designed for pediatric applications. Use with pediatric patients should be determined by attending physician.
- Consideration should be given to the following medical and anatomic conditions:
 - Distorted anatomy
 - Subcutaneous abscess
 - Hematoma
 - Post-operative scarring/radiation
 - Coagulopathies or systemic thrombolytic therapy

PRECAUTIONS

- This product is intended for use by physicians trained and experienced in emergency airway management techniques. Standard techniques for percutaneous placement of airway catheters should be employed.
- Patients in need of cricothyrotomy may have significant spinal injury. In patients who have sustained significant trauma, the cervical spine should be immobilized throughout the procedure, if possible.
- Whenever possible and appropriate, utilize aseptic technique and local anesthetic for the procedure.

INSTRUCTIONS FOR USE

1. Identify the cricothyroid membrane between the cricoid and thyroid cartilages. **(Fig. 1)**
2. Carefully palpate the cricothyroid membrane and, while stabilizing the cartilage, make a vertical incision in the midline. **NOTE:** Ensure that the incision is sufficiently large to allow passage of the dilator and airway catheter.
3. Attach the supplied syringe to either the introducer needle or the catheter introducer needle and advance the needle through the incision into the airway at a 45 degree angle to the frontal plane in a caudad direction, in the midline. **(Fig. 2)** **NOTE:** Entrance into the airway can be confirmed by aspiration on the syringe, resulting in free air return.
4. If using the catheter introducer needle, remove the syringe and needle, leaving the catheter in place. If using the introducer needle, remove only the syringe, leaving the needle in place.
5. Advance the soft, flexible end of the wire guide through the catheter or needle and into the airway several centimeters.
6. Remove the catheter or needle, leaving the wire guide in place. **(Fig. 3)**
7. Advance the handled dilator, tapered end first, into the connector end of the airway catheter until the handle stops against the connector. **NOTE:** This step may be performed prior to beginning the procedure. Use of lubrication on the surface of the dilator may enhance fit and placement of the airway catheter.
8. Advance the airway catheter/dilator assembly over the wire guide until the proximal stiff end of the wire guide is completely through and visible at the handle end of the dilator. **It is important to continually visualize the proximal end of the wire guide during the airway insertion procedure to prevent its inadvertent loss into the trachea.**
9. Maintaining wire guide position, continue to advance the airway catheter/dilator assembly over the wire guide with a reciprocating motion completely into the trachea. **(Fig. 4)** **Take care not to advance the tip of the dilator beyond the tip of the wire guide within the trachea.**
10. Remove the wire guide and dilator simultaneously.

11. If placing a standard (uncuffed) airway catheter, proceed to step #12. If placing a cuffed airway catheter, inflate the cuff with air using a syringe. **(Fig. 5) NOTE:** An 8-10 cc volume in the cuff will yield a cuff diameter of 22-29 mm. The inflation and deflation procedure is at the discretion of the physician.
WARNING: Inflation of the cuff with more than 20 cc is not recommended.
12. Fix the airway catheter in place with tracheostomy tape strip in standard fashion.
13. Connect the airway catheter, using its standard 15 mm connector, to an appropriate ventilatory device.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2016
Bloomington, IN 47404 U.S.A.