EN
2

Nasal Jejunal Feeding Tube

CS
4

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DA
6

Nasal-jejunal ernæringssonde

NL
7

Nasale jejunum-voedingssonde

FR
9

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
Nasal jejunal matnings slang

ENGLISH

These recommendations are designed to serve only as a general guideline. They are not intended to supersede institutional protocols or professional clinical judgement concerning patient care.

RxOnly

INTENDED USE

This device is intended to provide short-term enteral access for delivery of nutrition and/or medications to the small bowel.

A thorough understanding of the technical principles, clinical applications and risks associated with nasal jejunal feeding tube placement is necessary before using the device. The Nasal Jejunal Feeding Tube should only be used by or under the supervision of personnel trained in standard nasal gastric tube placement procedure.

CONTRAINDICATIONS

Contraindications to placement and use of a nasal jejunal feeding tube include, but are not limited to:

- Esophageal varices
- Gastric varices
- INR (international normalized ratio) > 1.3 (at time of insertion and/or expected at time of removal)
- Anticoagulated patients (anticoagulated at time of insertion and/or expected to be anticoagulated at time of removal)
- Pathologic coagulopathies
- History of bleeding disorder(s)
- Small or large bowel obstruction
- Ischemic bowel
- Peritonitis
- Esophageal stricture or obstruction
- Gastric obstruction

- Recent nasal, oral, esophageal or gastric surgery, or trauma
- Deviated septum
- Inability to pass the feeding tube through the nares
- Uncooperative patient

POTENTIAL ADVERSE EVENTS

Potential adverse events associated with placement and use of a nasal jejunal feeding tube include, but are not limited to:

- Bleeding
- Clogged or leaking feeding tube
- Sinusitis
- Premature displacement of the tube
- Aspiration
- Nasal irritation
- Sore throat

PRECAUTIONS

The coordination of the endoscope accessory channel size and compatible devices is essential in obtaining optimal results during a procedure. Wilson-Cook Nasal Jejunal Feeding Tubes are available in 8 French, requiring a minimum channel size of 3.2 mm, and 10 French, requiring a minimum channel size of 3.7 mm.

This device is not intended for use beyond 30 days.

WARNINGS

None known.

INSTRUCTIONS FOR USE

1. Pre-load wire guide, floppy tip first, through nasal jejunal feeding tube.
2. Advance pre-loaded feeding tube and wire guide unit through accessory channel of previously placed endoscope until it reaches third or fourth portion of duodenum.
3. Once feeding tube is in proper position, slowly withdraw endoscope making sure to keep feeding tube in place. **Note:** Distance markings are observed both endoscopically and externally to assure proper positioning of feeding tube.
4. Once endoscope is withdrawn completely, gently remove wire guide from feeding tube.
5. Dispose of wire guide per institutional guidelines for biohazardous medical waste.
6. Generously lubricate blunt end of nasal transfer tube with water soluble lubricant.
7. Advance nasal transfer tube through nose until it exits in back of throat. (*See fig. 1*) **Note:** If difficulty is encountered while attempting to advance nasal transfer tube, reaccess nares.
8. Grasp tip of nasal transfer tube and pull it out through mouth, leaving other end of nasal transfer tube in place through nose.
9. If placing 8 French feeding tube, thread tip of feeding tube through oral end of nasal transfer tube and advance until it exits through nasal end. If placing 10 French feeding tube, securely connect tip of feeding tube to oral end of nasal transfer tube. (*See fig. 2*)
10. Gently pull nasal transfer tube through nose until feeding tube emerges through nostril. (*See fig. 3*)
11. Dispose of nasal transfer tube per institutional guidelines for biohazardous medical waste.
12. Cut feeding tube to desired external length, then place barbed fitting onto feeding tube.
13. Secure feeding tube in place using hypoallergenic tape.
14. Confirm proper tube placement using x-ray prior to initiating feedings.

HOW SUPPLIED

Supplied sterile by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use product if there is doubt as to whether product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and/or their published literature. Refer to your local Cook sales representative for information on available literature.

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© 2019 Cook Medical
Wilson-Cook Medical, Inc.
Cook Medical Europe Ltd.