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## CANTATA® 2.5 AND 2.8 MICROCATHETERS AND CANTATA® 2.5 AND 2.8 DUO MICROCATHETERS

**CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).**

### DEVICE DESCRIPTION

The Cantata 2.5 and Cantata 2.5 Duo Microcatheters are shapeable-tip, braided, kink-resistant 2.5 French microcatheters with hydrophilic coating. The microcatheters have an inner diameter of 0.021 inch, which allows them to accept a wide variety of embolization coils and appropriately sized embolization particles or spheres.

The Cantata 2.8 and Cantata 2.8 Duo Microcatheters are shapeable-tip, braided, kink-resistant 2.8 French microcatheters with hydrophilic coating. The microcatheters have an inner diameter of 0.025 inch, which allows them to accept a wide variety of embolization coils and appropriately sized embolization particles or spheres.

The Cantata 2.5 Duo and Cantata 2.8 Duo Microcatheters offer an additional radiopaque marker band placed 3 cm proximal to the distal marker band.

### Cantata 2.5 and Cantata 2.5 Duo Microcatheter Flow Rates (ml/sec)\*

The infusion dynamic pressure should not exceed 1200 psi (82.7 bar). The static pressure should not exceed 300 psi (20.6 bar).

| PSI (bar)   | Infusion Medium                                  | Catheter Length |        |        |        |
|-------------|--|-----------------|--------|--------|--------|
|             |  | 100 cm          | 110 cm | 135 cm | 150 cm |
| 100 (6.8)   | 100% Saline                                      | 1.2             | 1.2    | 1.1    | 1.0    |
|             | 100% Omnipaque <sup>***</sup> (Iohexol) 300mg/mL | 0.2             | 0.2    | 0.1    | 0.1    |
|             | 100% ISOVUE <sup>****</sup> (iopamidol) 300mg/mL | 0.3             | 0.3    | 0.2    | 0.2    |
| 300 (20.6)  | 100% Saline                                      | 1.7             | 1.7    | 1.5    | 1.4    |
|             | 100% Omnipaque <sup>***</sup> (Iohexol) 300mg/mL | 0.6             | 0.6    | 0.5    | 0.4    |
|             | 100% ISOVUE <sup>****</sup> (iopamidol) 300mg/mL | 0.8             | 0.7    | 0.6    | 0.5    |
| 600 (41.3)  | 100% Saline                                      | 2.5             | 2.4    | 2.1    | 2.0    |
|             | 100% Omnipaque <sup>***</sup> (Iohexol) 300mg/mL | 1.3             | 1.2    | 1.0    | 0.8    |
|             | 100% ISOVUE <sup>****</sup> (iopamidol) 300mg/mL | 1.6             | 1.4    | 1.2    | 1.0    |
| 1200 (82.7) | 100% Saline                                      | 3.9             | 3.8    | 3.4    | 3.2    |
|             | 100% Omnipaque <sup>***</sup> (Iohexol) 300mg/mL | 2.6             | 2.4    | 2.0    | 1.7    |
|             | 100% ISOVUE <sup>****</sup> (iopamidol) 300mg/mL | 3.1             | 2.9    | 2.3    | 2.0    |

### Cantata 2.8 and Cantata 2.8 Duo Microcatheter Flow Rates (ml/sec)\*

The infusion dynamic pressure should not exceed 1000 psi (68.9 bar). The static pressure should not exceed 300 psi (20.6 bar).

| PSI (bar)  | Infusion Medium                                   | Catheter Length |        |        |        |
|------------|---|-----------------|--------|--------|--------|
|            |   | 100 cm          | 110 cm | 135 cm | 150 cm |
| 100 (6.8)  | 100% Saline                                       | 1.9             | 1.9    | 1.7    | 1.6    |
|            | 100% Omnipaque <sup>***</sup> (Iohexol) 300 mg/mL | 0.4             | 0.3    | 0.2    | 0.2    |
|            | 100% ISOVUE <sup>****</sup> (iopamidol) 300mg/mL  | 0.5             | 0.5    | 0.4    | 0.3    |
| 300 (20.6) | 100% Saline                                       | 2.8             | 2.7    | 2.4    | 2.3    |
|            | 100% Omnipaque <sup>***</sup> (Iohexol) 300 mg/mL | 1.0             | 1.0    | 0.8    | 0.7    |
|            | 100% ISOVUE <sup>****</sup> (iopamidol) 300mg/mL  | 1.3             | 1.2    | 1.0    | 0.9    |

| PSI<br>(bar)   | Infusion Medium                                      | Catheter Length |        |        |        |
|----------------|--|-----------------|--------|--------|--------|
|                |  | 100 cm          | 110 cm | 135 cm | 150 cm |
| 600<br>(41.3)  | 100% Saline  | 4.1             | 3.9    | 3.5    | 3.3    |
|                | 100% Omnipaque <sup>***</sup><br>(Iohexol) 300 mg/mL | 2.1             | 1.9    | 1.6    | 1.3    |
|                | 100% ISOVUE <sup>****</sup><br>(Iopamidol) 300mg/mL  | 2.5             | 2.4    | 2.0    | 1.7    |
| 1000<br>(68.9) | 100% Saline  | 5.8             | 5.5    | 5.0    | 4.6    |
|                | 100% Omnipaque <sup>***</sup><br>(Iohexol) 300mg/mL  | 3.4             | 3.2    | 2.6    | 2.2    |
|                | 100% ISOVUE <sup>****</sup><br>(Iopamidol) 300mg/mL  | 4.2             | 3.9    | 3.3    | 2.9    |

\* Microcatheter flow rates are approximate values. Omnipaque and ISOVUE warmed to 37° C. Saline at room temperature (21° C).

\*\*Omnipaque is a registered trademark of Amersham Health, Princeton, New Jersey.

\*\*\* ISOVUE is a registered trademark of Bracco Diagnostics Inc.

## INTENDED USE

The Cantata 2.5, Cantata 2.8, Cantata 2.5 Duo, and Cantata 2.8 Duo Microcatheters are intended for use in small vessel or superselective anatomy for diagnostic and interventional procedures including peripheral and coronary use.

## CONTRAINDICATIONS

None known

## WARNINGS

- Damage to the catheter may result from exposure to heat or steam.
- Do not exceed maximum dynamic or static pressures. Exceeding indicated pressures may result in microcatheter damage.
- If resistance is felt between the vessel and microcatheter, or between the microcatheter and wire guide, verify the cause of resistance by fluoroscopy or remove the microcatheter and wire guide as a unit. Damage may occur to the microcatheter and/or wire guide.
- If microcatheter occlusion is suspected, do not flush the microcatheter. Over-pressurization may occur, resulting in microcatheter rupture.

## PRECAUTIONS

- This product is intended for use by physicians trained and experienced in small vessel access and interventional procedures. Standard techniques for placement of percutaneous microcatheters should be employed.
- The Cantata 2.8 and Cantata 2.8 Duo Microcatheters fit through any angiographic catheter that accepts a 0.038 inch wire guide.
- The Cantata 2.5 and Cantata 2.5 Duo Microcatheters fit through any angiographic catheter that accepts a 0.035 inch wire guide.
- To avoid vascular trauma, maintain a wire guide in the microcatheter lumen and beyond the distal tip of the microcatheter during manipulation.
- When using the microcatheter for infusion of embolic agents, thrombolytic therapies, or other agents, refer to the appropriate manufacturer's suggested Instructions for Use for proper use guidelines.
- The potential effects of phthalates on pregnant/nursing women or children have not been fully characterized and there may be concern for reproductive and developmental effects.

## INSTRUCTIONS FOR USE

1. Using standard technique, place appropriate angiographic catheter or guiding catheter into the vascular system. A rotating hemostatic valve with side-arm adapter may be attached to the introducer device to prevent back bleeding and provide a continuous flush of heparinized saline or sterile water.
2. Remove the microcatheter spiral holder from its sealed packaging.
3. Attach a syringe filled with heparinized saline solution or sterile water to the Luer lock fitting of the microcatheter holder.
4. Inject enough solution to wet the microcatheter surface entirely. This will activate the hydrophilic coating on the microcatheter surface.  
**NOTE:** The surface of the microcatheter may become dry after removal from the holder. Additional wetting with heparinized saline or sterile water will renew the hydrophilic effect.
5. Attach a second hemostatic valve with side-arm adapter to the microcatheter, purge any air and flush with heparinized saline or sterile water.  
**NOTE:** The tip of the microcatheter can be reshaped with fingers or by using a 19 gage cannula. Do not use heat or steam.
6. Carefully insert an appropriately sized wire guide into the microcatheter and completely close the valve around the wire guide.

7. Introduce the microcatheter and wire guide assembly through the hemostatic valve. Gently close the valve to prevent backflow while still allowing back-and-forth movement of the microcatheter.
8. Using fluoroscopy, introduce the microcatheter and wire guide assembly into the vascular system, making sure the wire guide is always ahead of the microcatheter. Final positioning is accomplished by short advances of the wire guide and microcatheter until the desired position is achieved and then confirmed by fluoroscopic visualization.
9. Monitor microcatheter placement and position during use.

## HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

## REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2019
Bloomington, IN 47404 U.S.A.
Cook Medical Europe Ltd