
- EN**
5 **EchoTip ProCore® HD Ultrasound Biopsy Needle**
Instructions for Use
- BG**
8 **Игла за биопсия под ехографски контрол с висока разграничителна способност EchoTip ProCore®**
Инструкции за употреба
- CS**
12 **Bioptická jehla EchoTip ProCore® používající ultrazvuk s vysokým rozlišením**
Návod k použití
- DA**
16 **EchoTip ProCore® ultralydbiopsinål med høj opløsning**
Brugsanvisning
- NL**
20 **EchoTip ProCore® HD ultrasone biopsienaald**
Gebruiksaanwijzing
- ET**
24 **HD-ultraheli-biopsianõel EchoTip ProCore®**
Kasutusjuhised
- FR**
27 **Aiguille de biopsie ultrasonore à haute définition EchoTip ProCore®**
Mode d'emploi
- DE**
31 **EchoTip ProCore® HD Ultraschall-Biopsienadel**
Gebrauchsanweisung
- EL**
35 **Βελόνα βιοψίας υπερήχων υψηλής ευκρίνειας EchoTip ProCore®**
Οδηγίες χρήσης
- HU**
39 **EchoTip ProCore® nagy felbontású ultrahangos biopsziás tű**
Használati utasítás
- IT**
43 **Ago di biopsia per ecografia ad alta definizione EchoTip ProCore®**
Istruzioni per l'uso
- LV**
47 **EchoTip ProCore® augstas izšķirtspējas ultraskaņas biopsijas adata**
Lietošanas instrukcija
- LT**
51 **EchoTip ProCore® didelės skiriamosios gebos ultragarsinė biopsijos adata**
Naudojimo nurodymai
- NO**
55 **EchoTip ProCore® HD ultralydbiopsinål**
Bruksanvisning
- PL**
59 **Igła biopsyjna do USG o wysokiej rozdzielczości EchoTip ProCore®**
Instrukcja użycia
- PT**
63 **Agulha de biopsia de ultra-sons HD EchoTip ProCore®**
Instruções de utilização
- RO**
67 **Ac ecografic HD pentru biopsie EchoTip ProCore®**
Instrucțiuni de utilizare
- SK**
71 **Ultrazvuková bioptická ihla s vysokým rozlišením EchoTip ProCore®**
Návod na použitie
- ES**
75 **Aguja de biopsia por ecografía de alta definición EchoTip ProCore®**
Instrucciones de uso
- SV**
79 **EchoTip ProCore® ultraljudsbiopsinål med hög upplösning**
Bruksanvisning


Plunger

Бутало
Pist
Stempel
Plunjer
Kolb
Piston
Kolben
Έμβολο
Dugattyú
Stantuffo
Virzulis
Stüimoklis
Stempel
Tloczek
Èmbolo
Piston
Piest
Èmbolo
Kolv

Stopcock "open"

Кранче "отворено"
Otevřený uzavírací kohout
Hane "åben"
Afsluitkraan 'open'
Sulgemisventiil avatud
Robinet ouvert
Absperrhahn geöffnet
Στρόφιγγα ανοικτή
Elzárócsap „nyitott” állásban
Rubinetto aperto
Slægkrans atvërts
Kranelis „atviras"
Stoppekran åpen
Kranik w pozycji otwartej
Torneira de passagem "aberta"
Robinet „deschis"
Otvorený kohútik
Llave de paso «abierta»
Kran "öppen"

Sliding Sheath Adjuster

Плъзгащ се регулатор на защитната обвивка
 Posuvné nastavovací zařízení sheathu
 Glidende sheath-justering
 Schuifregeling van de sheath
 Liughülsi regulaator
 Ajusteur de gaine coulissant
 Verschiebbarer Schleuseneinsteller
 Συρόμενος προσαρμογέας θηκαριού
 Elcsúszatható hüvelyhossz-beállító
 Regolatore della guaina scorrevole
 Pärbidāmais ievadapvalka regulētājs
 Slankusis vamzdelio reguliatorius
 Glidende hylsejusteringsmekanisme
 Przesuwany regulator koszulki
 Regulador deslizante da bainha
 Ajustor glisant al tecii
 Posuvný nastavovač puzdra
 Ajustador deslizante de la vaina
 Glidande hylsjusterare

Needle Handle

Дръжка на иглата
 Rukojeť jehly
 Nålegreb
 Naaldhandvat
 Nõela käepide
 Poignée de l'aiguille
 Nadelgriff
 Λαβή βελόνας
 Tű fogantyúja
 Impugnatura dell'ago
 Adatas rokturis
 Adatos kotelis
 Nālehāndtak
 Uchwyty igły
 Punho da agulha
 Mănerul acului
 Růčka ihly
 Mango de la aguja
 Nålhandtag

Safety Ring

Предпазен пръстен
 Pojistný prsteneč
 Sikkerhedsring
 Veiligheidsring
 Kaitserõngas
 Bague de sécurité
 Sicherheitsring
 Δακτύλιος ασφάλειας
 Biztonsági gyűrű
 Anello di sicurezza
 Drošības gredzens
 Apsauginis žiedas
 Sikkerhetsring
 Pierścień bezpieczeństwa
 Anel de sigurança
 Inel de siguranță
 Bezpečnostný krúžok
 Anillo de seguridad
 Säkerhetslåsring

Sheath Reference Mark

Означение за дължината на защитната обвивка
 Referenční značka sheathu
 Referencemærke på sheath
 Referentiemarkering van de sheath
 Hülsi tähis
 Repère de référence sur la gaine
 Schleusenbezugsmarke
 Σήμανση αναφοράς θηκαριού
 Hüvelyhossz referenciája
 Indicatore di riferimento della guaina
 Ievadapvalka atskaites atzīme
 Vamzdelio atskaitos žymė
 Referansemerke på hylsen
 Znacznik referencyjny wysunięcia koszulki
 Marca de referència da bainha
 Marcaj de referință al tecii
 Referenčná značka puzdra
 Marca de referencia de la vaina
 Hylsans referensmarkering

Stylet Hub

Втулка на стилета
 Nástavec stiletu
 Stiletmuffe
 Aanzetstuk van het stilet
 Stileti kinnitus
 Embase du stylet
 Mandrinansatz
 Ομφαλός στείλειου
 Mandrin végidoma
 Connettore del mandrino
 Stiletes galviņa
 Zondo įvorė
 Stilettfeste
 Złączka mandrynu
 Conector do estilete
 Amboul stiletului
 Hrdlo sondy
 Conector del estilete
 Mandrängfattning

EchoTip ProCore® HD Ultrasound Biopsy Needle

IMPORTANT INFORMATION

Please review prior to use.

INTENDED USE

This device is used with an ultrasound endoscope for fine needle biopsy, (FNB), of submucosal lesions, mediastinal masses, lymph nodes and intraperitoneal masses within or adjacent to the gastrointestinal tract.

NOTES

Do not use this device for any purpose other than stated intended use.

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Store in a dry location, away from temperature extremes.

Use of this device restricted to a trained healthcare professional.

CONTRAINDICATIONS

Those specific to primary endoscopic procedure to be performed in gaining access to desired site.

Coagulopathy.

POTENTIAL COMPLICATIONS

Those associated with gastrointestinal endoscopy include, but are not limited to: perforation, hemorrhage, aspiration, fever, infection, allergic reaction to medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest, damage to blood vessels, nerve damage, and acute pancreatitis. Those associated with EUS needle biopsy include but are not limited to: pain, death, peritonitis, portal vein gas and thrombosis, pneumoperitoneum and tumor seeding of the needle tract.

WARNINGS

Not for use in the heart or vascular system.

The tip of the needle and stylet are sharp and could cause injury to the patient or user if not used with caution.

PRECAUTIONS

Refer to package label for minimum channel size required for this device.

Ensure the stylet is fully inserted when advancing the needle into the biopsy site.

When targeting multiple sites, replace device for each site.

Needle must be retracted into sheath and thumbscrew on safety ring must be locked to hold needle in place prior to introduction, advancement or withdrawal of device. Failure to retract needle may result in damage to endoscope.

SYSTEM PREPARATION

Illustrations

1. Examine syringe. It has two plunger locks that must be depressed to advance plunger. Tip of syringe has a Luer lock fitting with a stopcock on the side port. Air can be exchanged when stopcock is in “open” position, aligned with syringe. (See fig. 1)

2. Prepare syringe as follows:

- a. With stopcock in “open” position, depress plunger locks and fully advance plunger into syringe.
- b. Turn stopcock 90° to “closed” position.
- c. Pull plunger back until it is locked into place at desired setting, creating suction.
- d. Set prepared syringe aside until biopsy is desired.

3. Advance device into ultrasound endoscope to determine preferred sheath length. To adjust length, loosen thumbscrew lock on sliding sheath adjuster and slide until preferred length is attained. **Note:** Reference mark for sheath length will appear in sliding sheath adjuster window. (See fig. 2) Tighten thumbscrew on sliding sheath adjuster to maintain preferred sheath length.

INSTRUCTIONS FOR USE

1. Identify desired biopsy site by endoscopic ultrasound.

2. Introduce ultrasound needle into accessory channel and advance device in small increments until Luer lock fitting at base of sliding sheath adjuster meets fitting on accessory channel. (See fig. 3)

3. Attach device to accessory channel port by rotating device handle until fittings are connected.

4. With ultrasound endoscope and device straight, adjust needle to desired length by loosening thumbscrew on safety ring, and advancing

it until desired reference mark for needle advancement appears in the window of safety ring. (See fig. 4) Tighten thumbscrew to lock safety ring in place. **Note:** Number in safety lock ring window indicates extension of needle in centimeters. **Caution:** During needle adjustment or extension, ensure device is attached to accessory channel. Failure to attach device prior to needle adjustment or extension may result in damage to endoscope.

5. While maintaining position of ultrasound endoscope, extend needle by advancing needle handle to pre-positioned safety ring.

6. Advance needle into lesion.

7. Remove stylet from needle by gently pulling back on plastic hub seated in metal fitting of needle handle.

8. Attach Luer lock fitting of previously prepared syringe securely onto fitting on handle.

9. Turn stopcock to “open” position aligned with syringe, allowing negative pressure in syringe to facilitate biopsy. Gently move handle in small increments back and forth within biopsy site. **Note:** Do not remove needle from biopsy site during biopsy.

10. Upon completion of biopsy, turn stopcock to “close” position.

11. Retract needle into sheath by pulling back on handle shaft. Lock thumbscrew on safety ring at 0 cm mark.

12. Disconnect Luer lock fitting from accessory channel and withdraw entire device from ultrasound endoscope.

13. Unlock thumbscrew on safety ring and extend needle. Exchange locking syringe for an empty syringe; attach to Luer lock fitting on handle.

14. Push in plunger to expel specimen, then prepare per institutional guidelines.

15. Remaining material may be recovered for examination by reinserting stylet or flushing device.

Upon completion of procedure, dispose of device per institutional guidelines for biohazardous medical waste.

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© COPYRIGHT COOK 2016 ©
COOK IRELAND LTD.