COOK

MEDICAL

CE
0123

EN
5
BG
11
CS
18
DA
24
NL
30
ET
37
FR
43
DE
50
EL
57
HU
65
IT
72
LV
79
LT
85
NO
91
PL
97
PT
104
RO
111
SK
117
ES
123
SV
130

Zilver 635® Biliary Stent

Instructions for Use

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ZILVER 635® BILIARY STENT

EXPANDABLE METAL BILIARY STENT SYSTEM

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

WARNING: The safety and effectiveness of this device for use in the vascular system has not been established.

DEVICE DESCRIPTION

Zilver 635® Biliary Self-Expanding Stent

The Zilver Biliary stent is a self-expanding, open-ended cylinder of nitinol (**Fig. 1**). The design of the stent resists radial compression and easily conforms to bends in the bile duct anatomy without collapsing. There are radiopaque markers on each end of the stent to assist in fluoroscopic visualization of the stent position. It can be used in the Biliary tree from the intrahepatic ducts through the papilla. The device is recommended for use in the adult population.

Delivery System

The stent is compressed into a flexible delivery system (**Fig. 1**). There is a radiopaque marker at the distal end of the delivery system to aid in deployment of the stent under fluoroscopic guidance. The delivery system accepts a .035 inch wire guide.

INTENDED USE

This device is used in palliation of malignant neoplasms in the biliary tree.

CONTRAINDICATIONS

- Contraindications include those specific to ERCP and any procedure to be performed in conjunction with stent placement.
- Additional contraindications include, but are not limited to:
 - Inability to pass the wire guide or stent delivery system through the obstructed area
 - Very small intrahepatic ducts
 - Biliary duct strictures of benign etiology
 - Biliary obstruction preventing cholangiography
 - Concurrent perforated bile duct
 - Concurrent bile duct stones
 - Patients for whom endoscopic procedures are contraindicated
 - Patients with coagulopathy
 - Any use other than that specifically outlined in the Intended Use section of this booklet

WARNINGS

- The safety and efficacy of combined side-by-side with overlapping stents has not been established.
- This stent is **not intended to be removed** and is considered to be a permanent implant. Attempts to remove the stent after placement may cause damage to the surrounding mucosa.
- This device is **not** intended to be deployed through the wall of a previously placed or existing metal stent. Doing so could make it difficult or impossible to remove the delivery system.
- The stent contains nickel, titanium, and gold, which may cause allergic reaction in individuals with nickel, titanium, or gold sensitivity.
- These metal biliary stents are not intended to be repositioned or removed after deployment in the bile duct. In case of accidental deployment or improper placement (immediately following deployment), leave the stent in place.

PRECAUTIONS

- This device is intended for use by physicians trained and experienced in biliary stenting. Standard ERCP techniques for placement of biliary metal stents should be employed.
- This device is not compatible with the Cook Medical THSF wire guide.
- Refer to the package label for the minimum channel size required for this device.
- Prior to stent placement, a complete diagnostic evaluation should be performed to determine the appropriate stent length and diameter. The stent length chosen should allow for additional length on both ends of the stricture. **NOTE:** In the event a single stent will not adequately cover the stricture, a second stent of the same diameter should be placed. In relation to the lesion site, the area of narrowing furthest away from the papilla should be stented first and the area nearest to the papilla should be stented second. This second stent should provide adequate overlap (at least 1 cm) with the initially placed stent to ensure a bridging of the stricture between the stents.
- This stent must be placed under fluoroscopic monitoring.
- This stent must only be placed using the delivery system provided.
- This stent is intended for palliative treatment only. Alternate methods of therapy should be investigated prior to placement.
- After stent placement, alternate methods of treatment such as chemotherapy and irradiation may increase the risk of: a) stent migration due to tumor shrinkage, b) stent erosion of the tissue, and/or c) mucosal bleeding.
- Long-term patency with this stent has not been established. Periodic evaluation of the stent is advised.
- If the wire guide or delivery system cannot be advanced through the obstructed area, do not attempt to place the stent.
- Assessment must be made to determine the necessity of sphincterotomy or balloon dilation prior to stent placement. In the event sphincterotomy or balloon dilation is required, all appropriate cautions, warnings and contraindications must be observed.
- Take care not to kink the device during use.
- Do not expose the delivery system to organic solvents (e.g., alcohol).
- Do not excessively torque the device.
- Use of balloon dilation after stent placement has not been evaluated.
- Following stent placement if resistance is met during the withdrawal of the delivery system, carefully remove the delivery system and wire guide as a unit.

POTENTIAL ADVERSE EVENTS

Potential adverse events associated with ERCP include, but are not limited to: pancreatitis, cholangitis, cholecystitis, cholestasis, aspiration, perforation, hemorrhage, infection, liver abscess, sepsis, allergic reaction to contrast or medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

Additional adverse events that can occur in conjunction with biliary stent placement include, but are not limited to: trauma to the biliary tract or duodenum, perforation, obstruction of the pancreatic duct, stent migration, stent occlusion, tumor ingrowth or excessive hyperplastic tissue ingrowth, tumor overgrowth, stent misplacement, pain, fever, nausea, vomiting, inflammation, recurrent obstructive jaundice, bile duct ulceration, death (other than due to normal disease progression).

MRI SAFETY INFORMATION

This symbol means the device is MR conditional

Nonclinical testing has demonstrated the Zilver 635® Biliary Stent is MR Conditional (for either a single stent, a pair of overlapping stents or a pair of side-by-side stents). A patient with this device can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 3 T or 1.5 T
- Maximum spatial field gradient of 1900 gauss/cm (19 T/m)
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2 W/kg (Normal Operating Mode)

Under the scan conditions defined above, nonclinical testing results indicate the Zilver 635® Biliary Stent is expected to produce a maximum temperature rise of less than 5.0°C after 15 minutes of continuous scanning.

In non-clinical testing, the image artifact caused by the device extends approximately 14 mm from the Zilver 635® Biliary Stent when imaged with a spin echo pulse sequence and a 3 T MR system.

For US Patients only

Cook recommends that the patient register the MR conditions disclosed in this IFU with the MedicAlert Foundation. The MedicAlert Foundation can be contacted in the following manners:

Mail:	MediAlert Foundation International 2323 Colorado Avenue Turlock, CA 95382, USA
Phone:	888-633-4298 (toll free) 209-668-3333 from outside the US
Fax:	209-669-2450
Web:	www.medicalert.org

INSTRUCTIONS FOR USE

Equipment Required

- Duodenoscope
- Zilver Biliary Stent delivery system with appropriate stent length and diameter
- .035 inch wire guide of appropriate length
- Fusion wire lock (optional)

IMPORTANT: This stent system is intended for single use only. **Do not** attempt to reload a deployed or partially deployed stent.

Procedure Preparation

Single Stent

1. With the duodenoscope in place, introduce a wire guide, floppy tip first.
2. Advance the wire guide until it is fluoroscopically visualized in position through the stricture.
3. Using the included syringe, prepare the previously selected stent delivery system and stent by irrigating the wire guide inner lumen and stent lumen with sterile water. Flushing the device with contrast media is not recommended.
4. After irrigation of the stent, replace the cap.

NOTE: Ensure that the red safety lock is not inadvertently removed prior to steps 1-4 above.

Side-by-Side Stents

NOTE: For side-by-side stenting, a duodenoscope with a minimum channel size of 4.2 mm is required.

1. With the duodenoscope in place, introduce a wire guide, floppy tip first.
2. Advance the wire guide until it is fluoroscopically visualized in position through the stricture.
3. With the duodenoscope in place, introduce the second wire guide, floppy tip first.
4. Advance the second wire guide until it is fluoroscopically visualized in position through stricture.
5. Using the included syringes, prepare the two previously selected stent delivery systems and stents by irrigating the wire guide inner lumens and stent lumens with sterile water. Flushing the device with contrast media is not recommended.
6. After irrigation of the stent, replace the cap.

NOTE: Ensure that the red safety lock is not inadvertently removed prior to steps 1-6 above.

Stenting Procedure

Single Stent

NOTE: For best results, keep the wire guide wet.

1. Unlock the wire guide, if using the fusion wire guide locking device.
2. Backload the device over a pre-positioned wire guide, ensuring the wire guide exits the wire guide hub.
3. Introduce the device in short increments into the accessory channel of the duodenoscope. With the elevator open, continue advancing the device in short increments until it is endoscopically visualized exiting the duodenoscope.
4. Under fluoroscopic guidance, with the elevator open, pass the stent delivery system through the papilla and into the common bile duct. Advance the device in short increments until the stent is fluoroscopically visualized through the stricture. Confirm desired stent position fluoroscopically. **NOTE:** For transpapillary placement, reference the yellow endoscopic marker for stent positioning. A stent bridging the papilla should extend beyond the papilla and into the duodenum approximately 0.5 cm after deployment.
5. Prior to deploying stent, remove the red safety lock. **NOTE:** Ensure that the red safety lock is not inadvertently removed prior to this step.
6. Under fluoroscopic guidance, with the elevator open, begin deployment of the stent by holding the wire guide hub stationary and slowly pulling back on the handle. Simultaneously maintain stent position during deployment by withdrawing the delivery system from the duodenoscope. **WARNING:** Do not advance the delivery system up into the bile duct after stent deployment has been initiated. **CAUTION:** The delivery system cannot be readvanced over the stent to recapture the stent during deployment.
7. Continue simultaneously pulling back on the handle and withdrawing the delivery system until you fluoroscopically confirm that the stent is completely deployed. For transpapillary placement, endoscopic confirmation of complete stent deployment may also be used.
8. While maintaining wire guide position, with the elevator open, carefully remove the delivery system from within the expanded stent. Remove the delivery system from the duodenoscope while monitoring the position of the stent fluoroscopically. **NOTE:** Fluoroscopically ensure the stent has expanded sufficiently to allow for safe removal of the tip.

Side-by-Side Stents

NOTE: For best results, keep wire guides wet.

1. Unlock the wire guide, if using the Fusion Wire Guide Locking Device.
2. Backload the device over a pre-positioned wire guide, ensuring the wire guide exits the wire guide hub
3. Introduce the device in short increments into the accessory channel of the duodenoscope. With the elevator open, continue advancing the device in short increments until it is endoscopically visualized exiting the duodenoscope.
4. Under fluoroscopic guidance, with the elevator open, pass the stent delivery system through the papilla and into the common bile duct. Advance the device in short

increments until the stent is fluoroscopically visualized through stricture. Confirm desired stent position fluoroscopically. **NOTE:** For transpapillary placement, reference the yellow endoscopic marker for stent positioning. A stent bridging the papilla should extend beyond the papilla and into the duodenum approximately 0.5 cm after deployment.

NOTE: Repeat steps 1-4 for second device. Both devices must be in position prior to the deployment of either stent.

5. Prior to deploying stent(s), remove red safety lock(s). **NOTE:** Ensure that the red safety lock(s) is/are not inadvertently removed prior to this step.

6. Under fluoroscopic guidance, with the elevator open begin deployment of stent(s) by holding the wire guide hub(s) stationary and slowly pulling back on the handle(s). Simultaneously maintain stent(s) position during deployment by withdrawing the delivery system(s) from the duodenoscope. **WARNING:** Do not advance the delivery system up into the bile duct after stent deployment has been initiated. **CAUTION:** The delivery system cannot be readvanced over stent to recapture the stent during deployment.

7. Continue simultaneously pulling back on the handle(s) and withdrawing the delivery system(s) until it is fluoroscopically confirmed that the stent(s) is/are completely deployed. For transpapillary placement, endoscopic confirmation of complete stent deployment may also be used.

8. While maintaining wire guide position, with the elevator open, carefully remove the delivery system from within the expanded stent. Remove the delivery system from the duodenoscope while monitoring the position of the stent fluoroscopically. **NOTE:** Fluoroscopically ensure the stent has expanded sufficiently to allow for safe removal of the tip.

NOTE: Repeat step 8 for removal of second delivery system.

Upon completion of procedure, dispose of device(s) per institutional guidelines for biohazardous medical waste.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred. If the package is opened or damaged when received, do not use the device. Visually inspect the device with particular attention to kinks, bends, or breaks. If an abnormality is detected that would prohibit proper working condition, do not use the device. Notify Cook for return authorization.

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© COPYRIGHT COOK 2017 ©
COOK IRELAND LTD.