- EN**
3 **Flexor®** Introducer / **Flexor®** Guiding Sheath
Instructions for Use
- DA**
4 **Flexor®** indfører / **Flexor®** styresheath
Brugsanvisning
- DE**
5 **Flexor®** Einführschleuse / **Flexor®**
Führungsschleuse
Gebrauchsanweisung
- EL**
7 Εισαγωγέας **Flexor®** / Οδηγό θηκάρι **Flexor®**
Οδηγίες χρήσης
- ES**
8 **Introductor Flexor®** / **Vaina guía Flexor®**
Instrucciones de uso
- FR**
10 **Introduceur Flexor®** / **Gaine d'introduction Flexor®**
Mode d'emploi
- IT**
11 **Introduttore Flexor®** / **Guaina di introduzione Flexor®**
Istruzioni per l'uso
- NL**
13 **Flexor®** introducer / **Flexor®** geleidesheath
Gebruiksaanwijzing
- PT**
14 **Introdutor Flexor®** / **Bainha guia Flexor®**
Instruções de utilização
- SV**
16 **Flexor®** införrarhylsa / **Flexor®** styrhylsa
Bruksanvisning

FLEXOR® INTRODUCER / FLEXOR® GUIDING SHEATH

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

Flexor Introducers and Guiding Sheaths are designed to perform as a guiding sheath and/or introducer sheath. Flexor Introducers and Guiding Sheaths incorporate a Flexor shaft with a hemostasis valve and are provided with one or more dilators. These devices are available in various sizes, lengths and configurations. Configurations include differences in shaft construction, such as varying shaft stiffness and distal tip material and shape, dilator material, hydrophilic coating and/or distal radiopaque markers. Refer to the product label for product specifications (i.e., sheath inner diameter, dilator inner diameter, distal curve configuration, hemostasis valve design).

INTENDED USE

Flexor Introducers and Guiding Sheaths are intended to introduce therapeutic or diagnostic devices into the vasculature, excluding coronary and neuro vasculature.

CONTRAINDICATIONS

None known

WARNINGS

- If resistance is encountered during advancement of Flexor sheath, assess cause of resistance and consider dilation of any restriction identified or consider alternate treatment strategy. If Flexor sheath is advanced through resistance, force to remove the sheath will be higher, increasing the risk of sheath material or hub separation upon withdrawal.
- Reinsertion of dilator prior to removal of Flexor sheath increases the strength of the sheath and lessens the risk of device separation. If resistance is anticipated or encountered during withdrawal of Flexor sheath, consider carefully reinserting the dilator prior to continuing removal.

PRECAUTIONS

- This product is intended for use by physicians trained and experienced in diagnostic and interventional techniques. Standard techniques for placement of vascular access sheaths should be employed.
- In order to ensure device compatibility, choose a sheath size large enough to accommodate the maximum outer diameter of any devices that will be placed through the sheath.
- All interventional or diagnostic instruments used with this product should move freely through the valve and sheath to avoid damage.
- Do not attempt to heat or reshape the device.
- When inserting, manipulating or withdrawing a device through the sheath, always maintain sheath position.
- When puncturing, suturing or incising the tissue near the sheath, use caution to avoid damaging the sheath.
- The potential effects of phthalates on children or pregnant/nursing women have not been fully characterized and there may be concern for reproductive and developmental effects.

POTENTIAL ADVERSE EVENTS

Adverse events that may be associated with the use of an introducer set include, but are not limited to:

- Bleeding
- Extravasation
- Hematoma
- Vessel laceration
- Vessel perforation
- Local inflammation
- Local pain
- Access site infection
- Distal embolization

INSTRUCTIONS FOR USE

Sheath Introduction

1. Ensure that the inner diameter (ID) of the sheath is appropriate for the maximum diameter of the instruments to be introduced.
2. Using the side-arm of the valve, flush the sheath by filling the sheath assembly completely with heparinized saline.
3. Flush the dilator with heparinized saline.
4. Insert the dilator completely into the sheath. If the sheath has a Tuohy-Borst valve, tighten the valve around the dilator.
5. If the device has hydrophilic coating, activate the coating by wetting the outer surface of the device with heparinized saline. **NOTE:** For best results, maintain wetted condition of device during placement.

6. Using standard Seldinger technique, access the target vessel with the appropriate needle.
7. Insert a wire guide into the vessel through the needle.
8. Leave the wire guide in place and remove the needle.
9. Insert the dilator/sheath assembly over the wire guide.
10. Remove the wire guide and dilator, then aspirate and flush through the sheath side-arm.
11. Insert appropriately sized device(s) as needed.

Sheath Removal

1. Insert a wire guide until its tip extends at least 10 cm past the tip of the sheath.
2. Remove the sheath. Avoid applying traction to hub during removal. If resistance is anticipated or encountered during withdrawal of the Flexor sheath, consider reinserting the dilator and removing the sheath and dilator as a unit.
3. Remove the wire guide.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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Bloomington, IN 47404 U.S.A.
Cook Medical Europe Ltd