EN	Needle Knife
CS	Jehlový papilotom
DA	Papillotom kanylekniv
NL	Naaldmespapillotoom
FR	Papillotome-couteau-aiguille
DE	Nadelmesser-Papillotom
EL	Βελονοτόμος
HU	Tűkés
IT	Papillotomo ad ago
NO	Nålekniv
PL	Nóż igłowy
PT	Papilótomo de agulha cortante
ES	Papilótomo de aguja/precorte
SV	Nålknivspapillotom

ENGLISH

INTENDED USE

This device is used for accessing the common bile duct when standard methods of cannulation have been exhausted.

NOTES

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

If the package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Notify Cook for return authorization.

Do not use this device for any purpose other than the stated intended use.

Store in a dry location, away from temperature extremes.

Use of this device restricted to a trained healthcare professional.

CONTRAINDICATIONS

Contraindications include those specific to ERCP and any procedures to be performed in conjunction with papillotomy.

Contraindications to papillotomy include, but are not limited to: coagulopathy.

POTENTIAL COMPLICATIONS

Potential complications associated with ERCP include, but are not limited to: pancreatitis, cholangitis, aspiration, perforation, hemorrhage, infection, sepsis, allergic reaction to contrast or medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

PRECAUTIONS

Refer to package label for minimum channel size required for this device.

The wire guide diameter and the inner lumen of the wire-guided device must be compatible.

Injection of contrast during ERCP should be monitored fluoroscopically. Overfilling of the pancreatic duct may cause pancreatitis.

The maximum rated input voltage for this device is 1.67 kVp-p.

Safe and effective electrosurgery is dependent not only on equipment design but, to a large extent, on factors under the control of the operator. Any electrosurgical accessory constitutes a potential electrical hazard to the patient and the operator. Possible adverse effects include, but are not limited to: fulguration, burns, nerve and/or muscle stimulation and cardiac arrhythmia.

In order to enhance safety and effectiveness, it is important that these cautions as well as the instructions supplied with the electrosurgical unit be read, understood and followed when using this device.

Before using this papillotome, follow the recommendations provided by the electrosurgical unit manufacturer to ensure patient safety through the proper placement and utilization of the patient return electrode. Ensure a proper path from the patient return electrode to the electrosurgical unit is maintained throughout the procedure.

Examine all components to be used during the procedure. Do not use a device that has been cut, burned or damaged. Damaged device insulation may cause unsafe currents in either the patient or operator. Do not use a device with a handle that has been cracked, broken or has a damaged metal active post.

Switch the electrosurgical unit to the "off" position when it is not in use.

If a non-insulated wire guide is used in the papillotome, it must be removed prior to applying electrosurgical current.

SYSTEM PREPARATION

1. Uncoil papillotome and carefully **remove precurved stylet wire** from distal tip of device. **Note:** Do not extend or retract cutting wire while device is coiled or stylet wire is in place, as this may cause damage to papillotome and render it inoperable. Extend and retract cutting wire only when papillotome is straight and stylet wire is removed.
2. Inspect active cord. Cord must be free of kinks, bends or breaks to allow for accurate transfer of current. If any abnormalities are noted, **do not use** active cord.
3. With electrosurgical unit off, prepare equipment. Securely connect active cord to device handle and electrosurgical unit. Active cord fittings should fit snugly into both device handle and electrosurgical unit. Following instructions from electrosurgical unit manufacturer, position patient return electrode and connect it to electrosurgical unit.

INSTRUCTIONS FOR USE (SEE FIG. 1)

I. If using Intra Ductal Exchange (IDE) port and short wire.

1. Attach Wire Guide Locking Device to endoscope accessory port. **Note:** Remove Wire Stop handle from its retaining clip.
2. Ensure that Wire Stop is disengaged from Proximal Wire Port (PWP) hub and ensure that distal tip of Wire Stop is proximal to IDE port.
3. Insert distal tip of wire guide into IDE port and advance until flush with distal tip of papillotome.
4. Advance Wire Stop handle until it reaches PWP hub and connect Luer lock securely to PWP hub.
5. With cutting wire fully retracted into catheter, advance tip of papillotome through cap of Wire Guide Locking Device and continue advancing until endoscopically visible.
6. Disengage Wire Stop by releasing handle from PWP hub and remove Wire Stop from catheter.
7. Position catheter and extend cutting wire to desired length. **Note:** Cutting wire extension may be adjusted by loosening thumbscrew and advancing or retracting adjustable ring on handle. **Caution:** It is essential to tighten thumbscrew after adjusting cutting wire length and before proceeding with papillotomy. **Caution:** Papillotome must fully exit endoscope. When applying current, contact of cutting wire with endoscope may cause grounding, which can result in patient injury, operator injury, a broken cutting wire, and/or damage to endoscope.
8. Following electrosurgical unit manufacturer's instructions, verify desired settings and activate electrosurgical unit.
9. Insert cutting wire into papillary orifice and, using a continuous motion, activate electrosurgical unit. **Caution:** It is essential to move cutting wire while applying current. Maintaining cutting wire in one position may cause excessive focal coagulation, charring of tissue and/or damage to cutting wire.
10. Once papillotomy is complete, turn electrosurgical unit off and retract cutting wire into sheath. Disconnect active cord from device handle and from electrosurgical unit. Wipe active cord with a damp cloth to remove all foreign matter. Store in a loose coil. **Note:** Wrapping active cord tightly may damage device.
11. Papillotome may then be used to cannulate ductal system.

Note: Previously placed wire guide may be left in position, in order to facilitate introduction of other wire-guided devices. If wire guide is to remain in place while device is withdrawn, utilize following steps:

12. Prior to removing device, utilize reference marks on catheter to ensure that IDE port is within ductal system.
13. Fluoroscopically visualize radiopaque band at IDE port. Retract wire guide until radiopaque distal tip of wire guide passes band; a disengagement from wire guide lumen will occur.
14. Advance disengaged wire guide to maintain ductal access.
15. Lock wire guide into Wire Guide Locking Device and remove papillotome from endoscope accessory channel.

II. If using Proximal Wire Port (PWP) and a pre-positioned long wire guide.

Note: For best results, wire guide should be kept wet.

1. Remove Wire Stop.
2. With cutting wire fully retracted into catheter, advance papillotome over pre-positioned wire guide ensuring that wire exits catheter at PWP until endoscopically visible.
3. Continue advancing papillotome until it is endoscopically visible.

REFER TO STEPS 7-11 IN "SECTION I"; THEN RESUME WITH STEP 4 BELOW:

4. Remove papillotome using standard long wire exchange technique.

Upon completion of procedure, dispose of device per institutional guidelines for biohazardous medical waste.

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© 2013 Cook Medical
Wilson-Cook Medical, Inc.
Cook Ireland Ltd.