


Deliver a uniform,
radial neurolysing
spray injection.



EchoTip[®] Ultra

CELIAC PLEXUS NEUROLYSIS NEEDLE

When your patients are experiencing severe chronic abdominal pain, the more efficient your CPN needle, the better. The highly echogenic EchoTip Celiac Plexus Neurolysis Needle, with its unique cone-shaped tip and pattern of side holes, allows you to smoothly deliver a neurolysing agent to the celiac plexus.



EchoTip® Ultra

CELIAC PLEXUS NEUROLYSIS NEEDLE

Used to deliver neurolysing agents to the Celiac Plexus under endoscopic ultrasound guidance.

Supplied sterile and is disposable - intended for single use only.

- A uniform, radial spray of neurolysing agent may be introduced into the Celiac Plexus through the multiple side holes of the cone-tipped needle
- Naturally-contoured handle delivers enhanced control and stability
- Safety lock ring easily slides and locks at the desired needle extension
- Integrated sheath adjuster offers compatibility with a full range of ultrasound endoscopes

Order Number	Reference Part Number	Needle Gage	French Size Fr	French Size Fr	Tip w/ Side Holes	Minimum Accessory Channel mm
G49257	ECHO-20-CPN	20	6	0-8	Cone Tip	2.8

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

Customer Service

EU Website: cookmedical.eu
 EDI: cookmedical.eu/edi
 Distributors: +353 61239240, ssc.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001702, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 Finland: +358 972519996, fi.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 6950072804, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Iceland: +354 800 7615, is.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Norway: +47 23162968, no.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland - French: +41 448009609, fr.orders@cookmedical.com
 Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland - German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com

USA Website: cookmedical.com

EDI: cookmedical.com/edi

Americas:

Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335

E-mail: orders@cookmedical.com

Australia:

Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283

E-mail: cau.custserv@cookmedical.com



AI-ESC-IR-OHNS-PI-RH-SUR-A4