

Enhance your diagnostic capabilities.



EchoTip® Ultra

ENDOBONCHIAL HD ULTRASOUND NEEDLE

Cook's EchoTip® Ultra Endobronchial HD Ultrasound Needle is designed to enhance your diagnostic and staging capabilities. With its high echogenicity, the EchoTip® Ultra aids needle visualization, allowing you to more adequately assure the needle tip is within the target, which may improve cytological yields. The ergonomic handle and scope adapter can also optimise procedural performance, stability and security.



cookmedical.com



Olympus



Pentax

EchoTip® Ultra

ENDOBONCHIAL HD ULTRASOUND NEEDLE

Olympus - Used to sample targeted submucosal and extramural lesions within or adjacent to the tracheobronchial tree or gastrointestinal tract through the accessory channel of an ultrasound endoscope for Fine Needle Aspiration.

Pentax - Used to sample targeted submucosal and extramural lesions within or adjacent to the tracheobronchial tree through the accessory channel of an ultrasound endoscope for Fine Needle Aspiration.

This device is supplied sterile and is disposable - intended for single use only.

- High Definition dimpling pattern for needle visibility
- Naturally contoured handle offers needle control & stability
- On scopes with metal luer hubs, the removable sheath adjuster allows direct connection to the scope
- On scopes with metal non-luer hubs, the sturdy scope adaptor design allows:
 - secure connection and procedural stability
 - suction and injection capability due to the airtight valve
- Integrated sheath adjuster offers compatibility with a full range of ultrasound endoscopes

Order Number	Reference Part Number	Needle Size gage	Sheath Size Fr	Adjustable Needle Extension cm	Stylet Tip	Minimum Accessory Channel mm
G52011	ECHO-HD-22-EBUS-O*	22	4.1	0-5	beveled	2
G53408	ECHO-HD-22-EBUS-P**	22	4.1	0-5	beveled	2
G24062	ECHO-HD-25-EBUS-O*	25	4.1	0-5	beveled	2
G24061	ECHO-HD-25-EBUS-P**	25	4.1	0-5	beveled	2

*Needle attaches to scope adapter on Olympus Endobronchial Ultrasound scope with non-luer metal hub.

**Needle attaches to luer metal hub on Pentax Endobronchial Ultrasound scope.

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

Customer Service

EMEA: EDI - www.cookmedical.com/edi.do
 Distributors: +353 61239240, ssc.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001633, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 Finland: +358 972519996, fi.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 6950072804, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Norway: +47 23162968, no.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland - French: +41 448009609, fr.orders@cookmedical.com
 Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland - German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com

cookmedical.com

Americas: EDI - www.cookmedical.com/edi.do
 Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335
 E-mail: orders@cookmedical.com

Australia:
 Phone: +61 738411188, 1800777222, Fax: +61 738411288, 1800077283
 E-mail: cau.custserv@cookmedical.com

