

LEGENDARY PERFORMANCE COMES FROM WITHIN

SOFTWARE SOLUTIONS

In 1984, LABORIE introduced the world to computerized urodynamics. Thirty years later we maintain our unwavering passion to bring the most capable software tools to the market. As the most comprehensive Urodynamics software suite, UDS120™ and i-LIST™ offer a world-class experience including:

POWERFUL: Our intelligent i-LIST™ platform offers complex, retrospective patient and study data analysis at a whole new level.



STREAMLINED: Manages appointments, previous visits, patient information and test data providing security and efficiency.

FEATURE RICH: Preloaded with 45 software features, UDS120™ offers facilities of all sizes a full range of testing abilities.



AQUARIUS™ CT

CAPABLE - PROGRESSIVE - COMPREHENSIVE



USA:
Tel.: +1 802 857 1300
800 522 6743
Email: usmarketing@laborie.com

EUROPE/INTERNATIONAL:
Tel.: +31 53 4803700
Email: info@laborie.com

HEADQUARTERS:
Tel.: +1 905 612 1170
Email: marketing@laborie.com



ANNIVERSARY
www.laborie.com



Innovation for Health

AQUARIUS PLATFORM TECHNOLOGY

INTRODUCING AQUARIUS™CT

The Aquarius CT pushes the boundaries of innovation by delivering industry leading, state-of-the-art technology. It boasts advanced Urodynamics capabilities while offering efficiency through user-centric design.

Whether being used in a small clinic or large hospital, the Aquarius CT exceeds expectation.

CAPABLE: From basic Uroflowmetry to advanced Video Urodynamics, the Aquarius CT adapts to fit your needs.

PROGRESSIVE: Touch-screen interfaces deliver the most natural way of performing Urodynamics for the physician and patient.

COMPREHENSIVE: Advanced specifications together with a custom designed industrial computer delivers performance.

THE ANATOMY OF EXCELLENCE



QUALITY, INNOVATION, SERVICE, SATISFACTION

The AQUARIUS is the ultimate expression of our core company values. The development of this revolutionary platform is predicated on exhaustive attention to user feedback, clinical needs and technological advancement.

The industry's first universally modular product suite, the Aquarius platform celebrates LABORIE's commitment to Quality, Innovation, Service and Satisfaction like never before.

The height adjustable Workstations of the Aquarius platform deliver enhanced ergonomics at the push of a button.



The **Central Hub** provides unparalleled ease-of-use by enabling device management at the touch of a finger.

The fifth generation of **Urocap** has been thoughtfully redesigned for perfect integration within the Aquarius system platform.

The **Pump Module** features state-of-the-art silencing technology, a wide range of flow rates and multiple calibration line storage.

The gold-standard **Spinning Disk** is a highly accurate and reliable artefact-resistant flowmeter with continuous flow-through technology.

The core of the measurement system, the **Roam DX** delivers high fidelity performance in a state-of-the-art package.

Bubble Detection and Infusion Measurement modules heighten precision and procedural efficiency.

The New Aquarius **Uroflow Stand** marks dramatic upgrades in size, quality and usability over any predecessors.

The **UPX Wireless Puller** boasts the industry's smallest footprint and completely wireless battery operation. Available on a stand-alone rolling floor stand or a bracket mounted to the Aquarius system.