



Right hepatic duct: Fusion Titan 6 mm biliary dilation balloon





Left hepatic duct: Fusion Titan 8 mm biliary dilation balloon

Images courtesy of Prof. Andreas Püspök Medical University Vienna, Austria



COOK[®] Fusion[®] Titan[™] BILIARY DILATION BALLOON



Designed to improve dilation outcomes

The Fusion Titan is the only biliary dilation balloon on the market with a multiport balloon catheter within the balloon body. It is a design that may help to anchor the balloon within the stricture and potentially minimize balloon movement during dilation. The low-profile balloon delivers the expected radial force to achieve successful dilation outcomes. And, depending on the clinical preference, the Fusion Titan biliary dilation balloon can be monorailed or delivered over a traditional long wire guide.