



# EXPANSIVE OPTIONS for biliary stone management.




A surgeon saw a need. We listened. No matter your approach, our combined sets and individual devices can help you manage biliary stones.

DRAIN

Some lodged stones can't be removed at the time of a cholecystectomy. Our biliary **stent** is an alternative to external draining catheters and T-tubes. Placed within the common bile duct and the duodenum, it allows bile drainage around the lodged stones.




The double-flange stent is positioned across the ampulla to allow for bile drainage.




DIAGNOSE


Access, view, and diagnose bile duct stones. We have **cholangiography** devices that can be used in percutaneous approaches and laparoscopic approaches.




REMOVE

Our **stone management** devices are sold as kits and as individual pieces. You can have everything you may need, gathered together, to fit any patient's need.





CHOLANGIOGRAPHY




**Fanelli Cholangiography Catheter Set**  
Used for diagnostic evaluation of the bile ducts during laparoscopic cholecystectomy procedures.

Order Number	Reference Part Number	Catheter Fr	Catheter ID Inch	Wire Guide Diameter Inch	Wire Guide Length cm
G12950	C-FCC-100	5.0	.044	.021	90




**Franklin Endoscopic Cholangiography Set**  
Used for diagnostic evaluation of the bile ducts during percutaneous cholecystectomy procedures.

Order Number	Reference Part Number	Catheter Fr	Catheter Length cm	Wire Guide Diameter Inch	Wire Guide Length cm
<b>Radiopaque</b>					
G10456	C-ECS-400-60-FECS	4.0	60	.025	145




**Mixer Endoscopic Cholangiography Set**  
Used for diagnostic evaluation of the bile ducts during percutaneous cholecystectomy procedures.

Order Number	Reference Part Number	Catheter Fr	Catheter Length cm	Catheter End Hole Diameter Inch	Needle Stylet Length cm
G08045	C-OECS-400-MIXTER	4.0	50	.025	52
G09331	C-OECS-500-MIXTER	5.0	50	.035	52



**Olsen Endoscopic Cholangiography Set**  
Used for access and diagnostic evaluation of the bile ducts during percutaneous cholecystectomy procedures.


Order Number	Reference Part Number	Catheter Fr	Catheter Length cm	Introducer Sheath Needle gage	Introducer Sheath Needle Length cm
G05926	C-OECS-300	3.0	43	16	12
G05927	C-OECS-400	4.0	43	14	12



**Petelin Endoscopic Cholangiography Set**  
Used for access and diagnostic evaluation of the bile ducts during percutaneous cholecystectomy procedures.


Order Number	Reference Part Number	Catheter Fr	Catheter Length cm	Catheter Needle gage	Catheter Needle Length cm
G10664	C-OECS-400-PECS	4.0	60	14	7.5

STONE MANAGEMENT



**Common Bile Duct Exploration Set**

Order Number	Reference Part Number
G26908	C-CDES-100



**Appel-Berci Cystic Duct Introducer Set**  
Used for placement of choledochoscopes, cholangiography catheters and other instruments.


Order Number	Reference Part Number	Sheath Lumen Diameter mm	Sheath Length cm	Introducer Fr	Accepts Instrumentation Diameter Fr
G08297	C-CDIS-4.0-15-BERCI	4.0	15	12.0	12.0

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.




**ATB Advance® Balloon Catheter**

Order Number	Reference Part Number	Catheter Fr	Catheter Length cm	Accepts Wire Guide Diameter inch	Inflated Balloon Diameter mm	Balloon Length cm
G36347	C-ATB5-35-40-8-4.0-CDES	5.0	40	.035	8	4




**NCompass® Nitinol Stone Extractor**  
Used for extraction of stones or debris during biliary surgical procedures.

Order Number	Reference Part Number	Fr	Length cm	Basket Configuration	Basket Diameter cm
G36250	C-NTSE-2.4-115-NC3	2.4	115	3 wires, tipped	1.5
G36251	C-NTSE-2.4-115-NCT4	2.4	115	4 wires, tipless	1.5



**Nathanson Transcystic Bile Duct Stone Exploration Pack**  
Used for extraction of bile duct stones via a transcystic approach.

Order Number	Reference Part Number	Catheter Length cm
G03298	CDES-550-NATHANSON	70



**Strange Cholangiography Set**  
Used for cholangiography and bile duct stone retrieval using fluoroscopy. The device is indicated for adults only.


Order Number	Reference Part Number	Catheter Fr	Catheter Length cm
G09108	C-CSRC-7.5D-40-STRANGE	7.5	40



**Three-Way Multipurpose Adapter**

Order Number	Reference Part Number
G13070	C-MPA-200


STENTING



**Fanelli Laparoscopic Endobiliary Stent Set**  
Used for internal biliary drainage.

Order Number	Reference Part Number	Stent Fr	Stent Length Between Flaps cm	Stent Assembly Fr	Stent Assembly Length cm	Wire Guide Diameter inch	Wire Guide Length cm
G13699	C-FBSS-100	7.0	6	8.5	50	.035	150

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.



No matter your approach, we have the tools to help.

## Customer Service

EU Website: [cookmedical.eu](http://cookmedical.eu)  
 EDI: [cookmedical.eu/edi](http://cookmedical.eu/edi)  
 Distributors: +353 61239240, [ssc.distributors@cookmedical.com](mailto:ssc.distributors@cookmedical.com)  
 Austria: +43 179567121, [oe.orders@cookmedical.com](mailto:oe.orders@cookmedical.com)  
 Belgium: +32 27001702, [be.orders@cookmedical.com](mailto:be.orders@cookmedical.com)  
 Denmark: +45 38487607, [da.orders@cookmedical.com](mailto:da.orders@cookmedical.com)  
 Finland: +358 972519996, [fi.orders@cookmedical.com](mailto:fi.orders@cookmedical.com)  
 France: +33 171230269, [fr.orders@cookmedical.com](mailto:fr.orders@cookmedical.com)  
 Germany: +49 6950072804, [de.orders@cookmedical.com](mailto:de.orders@cookmedical.com)  
 Hungary: +36 17779199, [hu.orders@cookmedical.com](mailto:hu.orders@cookmedical.com)  
 Iceland: +354 800 7615, [IS.orders@cookmedical.com](mailto:IS.orders@cookmedical.com)  
 Ireland: +353 61239252, [ie.orders@cookmedical.com](mailto:ie.orders@cookmedical.com)  
 Italy: +39 0269682853, [it.orders@cookmedical.com](mailto:it.orders@cookmedical.com)  
 Netherlands: +31 202013367, [nl.orders@cookmedical.com](mailto:nl.orders@cookmedical.com)  
 Norway: +47 23162968, [no.orders@cookmedical.com](mailto:no.orders@cookmedical.com)  
 Spain: +34 912702691, [es.orders@cookmedical.com](mailto:es.orders@cookmedical.com)  
 Sweden: +46 858769468, [se.orders@cookmedical.com](mailto:se.orders@cookmedical.com)  
 Switzerland - French: +41 448009609, [fr.orders@cookmedical.com](mailto:fr.orders@cookmedical.com)  
 Switzerland - Italian: +41 448009609, [it.orders@cookmedical.com](mailto:it.orders@cookmedical.com)  
 Switzerland - German: +41 448009609, [de.orders@cookmedical.com](mailto:de.orders@cookmedical.com)  
 United Kingdom: +44 2073654183, [uk.orders@cookmedical.com](mailto:uk.orders@cookmedical.com)

USA Website: [cookmedical.com](http://cookmedical.com)  
 EDI: [cookmedical.com/edi.do](http://cookmedical.com/edi.do)  
**Americas:**  
 Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335  
 E-mail: [orders@cookmedical.com](mailto:orders@cookmedical.com)  
**Australia:**  
 Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283  
 E-mail: [cau.custserv@cookmedical.com](mailto:cau.custserv@cookmedical.com)

AI, ESC, IR, OHNS, PI, RH, SUR-A4