

Gain access when standard methods fail.



Fusion[®]
NEEDLE KNIFE

DUAL PLATFORM ERCP

ACCESS

STONE MANAGEMENT

STRICTURE MANAGEMENT



www.cookmedical.com



Fusion®

NEEDLE KNIFE

Used for accessing the common bile duct when standard cannulation methods have been exhausted.

- Aids in accessing difficult-to-cannulate ducts
- Adjustable 4 mm needle accommodates precise control and 15 cm from the distal tip of the catheter

Order Number	Reference Part Number	Catheter Size Fr	Catheter Tip Fr	Adjustable Cutting Wire Length mm	Wire Guide Diameter inch	Minimum Accessory Channel mm
G31539	FS-PRECURT	7	6	4	.035	4.2

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

Customer Service

EMEA: EDI - www.cookmedical.com/edi.do
 Distributors: +353 61239240, ssc.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001633, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 Finland: +358 972519996, fi.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 6950072804, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Norway: +47 23162968, no.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland - French: +41 448009609, fr.orders@cookmedical.com
 Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland - German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com

www.cookmedical.com

Americas: EDI - www.cookmedical.com/edi.do
 Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335
 E-mail: orders@cookmedical.com

Australia:
 Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283
 E-mail: cau.custserv@cookmedical.com

