


# Maximize cell collection and preservation.



## Cytology Brushes

To better meet your clinical needs, we offer one of the most complete selections of cytology brushes on the market. Uniquely designed to maximize cell collection and preservation, our colonic, esophageal and bronchial brushes are available in a wide array of lengths and configurations, allowing for varying procedural needs and clinician preference.



## Cytology Brushes

Used to collect cells from the bronchi and upper and lower gastrointestinal tracts.

This device is supplied sterile and is disposable - intended for single use only.

- Brush bristles designed for maximum tissue acquisition
- Multiple brush diameters and catheter lengths to meet all your clinical needs

Order Number	Reference Part Number	Brush Diameter mm	Sheath Size Fr	Length cm	Minimum Accessory Channel mm	Color Code	Box Quantity
G22124	BCB-5-120-1-S	1	5	120	2	Blue	10
G22178	BCB-5-120-2-S	2	5	120	2	Blue	10
G21386	BCB-5-120-3-S	3	5	120	2	Blue	10
G22409	ECB-5-180-2-S	2	5	180	2	Green	10
G22174	ECB-5-180-3-S	3	5	180	2	Green	10
G22108	CCB-7-240-3-S	3	7	240	2.8	Red	10

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

### Customer Service

EU Website: [cookmedical.eu](http://cookmedical.eu)  
 EDI: [cookmedical.eu/edi](http://cookmedical.eu/edi)  
 Distributors: +353 61239240, [ssc.distributors@cookmedical.com](mailto:ssc.distributors@cookmedical.com)  
 Austria: +43 179567121, [oe.orders@cookmedical.com](mailto:oe.orders@cookmedical.com)  
 Belgium: +32 27001633, [be.orders@cookmedical.com](mailto:be.orders@cookmedical.com)  
 Denmark: +45 38487607, [da.orders@cookmedical.com](mailto:da.orders@cookmedical.com)  
 Finland: +358 972519996, [fi.orders@cookmedical.com](mailto:fi.orders@cookmedical.com)  
 France: +33 171230269, [fr.orders@cookmedical.com](mailto:fr.orders@cookmedical.com)  
 Germany: +49 6950072804, [de.orders@cookmedical.com](mailto:de.orders@cookmedical.com)  
 Hungary: +36 17779199, [hu.orders@cookmedical.com](mailto:hu.orders@cookmedical.com)  
 Ireland: +353 61239252, [ie.orders@cookmedical.com](mailto:ie.orders@cookmedical.com)  
 Italy: +39 0269682853, [it.orders@cookmedical.com](mailto:it.orders@cookmedical.com)  
 Netherlands: +31 202013367, [nl.orders@cookmedical.com](mailto:nl.orders@cookmedical.com)  
 Norway: +47 23162968, [no.orders@cookmedical.com](mailto:no.orders@cookmedical.com)  
 Spain: +34 912702691, [es.orders@cookmedical.com](mailto:es.orders@cookmedical.com)  
 Sweden: +46 858769468, [se.orders@cookmedical.com](mailto:se.orders@cookmedical.com)  
 Switzerland - French: +41 448009609, [fr.orders@cookmedical.com](mailto:fr.orders@cookmedical.com)  
 Switzerland - Italian: +41 448009609, [it.orders@cookmedical.com](mailto:it.orders@cookmedical.com)  
 Switzerland - German: +41 448009609, [de.orders@cookmedical.com](mailto:de.orders@cookmedical.com)  
 United Kingdom: +44 2073654183, [uk.orders@cookmedical.com](mailto:uk.orders@cookmedical.com)

USA Website: [cookmedical.com](http://cookmedical.com)

EDI: [cookmedical.com/edi.do](http://cookmedical.com/edi.do)

#### Americas:

Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335

E-mail: [orders@cookmedical.com](mailto:orders@cookmedical.com)

#### Australia:

Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283

E-mail: [cau.custserv@cookmedical.com-A4](mailto:cau.custserv@cookmedical.com-A4)



AI, ESC, IR, OHNS, PI, RH, SUR-A4