#### From cytology to histology: Continuing to increase your predictive values.




### EchoTip ProCore®




# Get to the core issue throughout the entire GI tract.

EchoTip ProCore allows for collection of histology with cytology while providing access to the entire GI tract, where you can target smaller lesions and increase your yields, all while decreasing your needle passes.<sup>1</sup>











"Cytology is not enough in the diagnosis and treatment of all lesions. You need sufficient quality of material to correctly characterise tumours and conditions."

Marc Giovannini, MD, PhD, Endoscopic Unit, Paoli-Calmettes Institute Marseilles, France



## Leading innovation in EUS since 1998.

Over the years, you've told us what matters to you and we've listened. Incorporating extensive feedback provided by clinicians, we offer you our latest developments in needle flexibility and sampling capability, tailored to you and your pathologist's needs and heavily supported by published clinical data.





#### **Product Specifications**

Used with an ultrasound endoscope for fine needle biopsy (FNB) of submucosal lesions, mediastinal masses, lymph nodes and intraperitoneal masses within or adjacent to the gastrointestinal tract.

Supplied sterile and is disposable - intended for single use only.

Order Number	Reference Part Number	Needle Size gauge	Sheath Size Fr	Adjustable Needle Extension cm	Stylet Tip	Minimum Accessory Channel mm
G53585	ECHO-HD-19-C	19	4.8	0-8	recessed ball tip	2.8
G34785	ECHO-HD-3-20-C*	20	7.95	0-8	recessed ball tip	3.7
G55736	ECHO-HD-22-C	22	5.2	0-8	recessed ball tip	2.8
G55738	ECHO-HD-25-C	25	5.2	0-8	beveled	2.8

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

#### **Customer Service**

EU Website: cookmedical.eu EDI: cookmedical.eu/edi

Distributors: +353 61239240, ssc.distributors@cookmedical.com Austria: +43 179567121, oe.orders@cookmedical.com Belgium: +32 27001633, be.orders@cookmedical.com Denmark: +45 38487607, da.orders@cookmedical.com Finland: +358 972519996, fi.orders@cookmedical.com France: +33 171230269, fr.orders@cookmedical.com Germany: +49 6950072804, de.orders@cookmedical.com Hungary: +36 17779199, hu.orders@cookmedical.com Ireland: +353 61239252, ie.orders@cookmedical.com Italy: +39 0269682853, it.orders@cookmedical.com Netherlands: +31 202013367, nl.orders@cookmedical.com

Switzerland - French: +41 448009609, fr.orders@cookmedical.com Switzerland - Italian: +41 448009609, it.orders@cookmedical.com Switzerland - German: +41 448009609, de.orders@cookmedical.com USA Website: cookmedical.com EDI: cookmedical.com/edi.do

Americas:

Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335 E-mail: orders@cookmedical.com

Australia:

Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283 E-mail: cau.custserv@cookmedical.com





<sup>\*</sup>ReCoil stylet is only available on EchoTip ProCore 20ga.