

ERCP is all about smooth,
easy access.



Fusion[®] OMNI[™]
ERCP CATHETER WITH DOMETIP[®]

DUAL PLATFORM ERCP

ACCESS

STONE MANAGEMENT

STRICTURE MANAGEMENT



www.cookmedical.com



Fusion[®] OMNI[™]

ERCP CATHETER WITH DOMETIP[®]

Used for endoscopic cannulation of the ductal system.

- DomeTip designed with anatomy in mind for smooth, easy access
- Breakthrough channel technology- the option for physician wire guide control or the ability to use the device in an over-the-wire long wire technique
- The unique injection lumen configuration allows for ease of injection

Order Number	Reference Part Number	Catheter Size Fr	Catheter Length cm	Tip Shape	Wire Guide Diameter inch	Minimum Accessory Channel mm
G31590	FS-GT-OMNI	7	200	DomeTip	.035	4.2

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

Customer Service

EMEA: EDI - www.cookmedical.com/edi.do
 Distributors: +353 61239240, ssc.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001633, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 Finland: +358 972519996, fi.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 6950072804, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Norway: +47 23162968, no.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland - French: +41 448009609, fr.orders@cookmedical.com
 Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland - German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com

Americas: EDI - www.cookmedical.com/edi.do
 Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335
 E-mail: orders@cookmedical.com

Australia:
 Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283
 E-mail: cau.custserv@cookmedical.com



www.cookmedical.com